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Instructions For Use

Controller Full HD USB Image/Video Capture

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	Content	Page
0	Important Information about this Document.....	67
0.1	Scope of Validity, Identification, Purpose.....	67
0.2	Target Group.....	67
0.3	Using and Storing this Document.....	67
0.4	Supplementary Documents.....	67
1	General Information about the Device.....	68
1.1	Scope of Delivery.....	68
1.2	Product Description.....	68
1.2.1	Performance Characteristics and Function.....	68
1.2.2	Front View of the Controller.....	69
1.2.3	Rear View of the Controller.....	70
1.2.4	Menu - Functions, Access and Navigation.....	71
1.3	Usage.....	72
1.3.1	Intended Use.....	72
1.3.2	Indications.....	72
1.3.3	Contraindications.....	72
1.4	Conformity.....	73
1.4.1	Standards and Directives.....	73
1.4.2	Medical Device Classification.....	73
1.5	Marking.....	74
1.5.1	Pictograms and Information on the Device and Packaging.....	74
1.5.2	Pictograms in this Document.....	75
1.6	Service Department Contact Details.....	75
1.7	Obligation to Report Serious Incidents.....	75
2	General Safety Information.....	76
2.1	Safety Messages in this Document.....	76
2.1.1	Safety Messages at the Start of a Chapter.....	76
2.1.2	Safety Messages in the Body of the Text.....	77
2.2	Product Safety.....	77
2.2.1	Basic Safety Notices.....	77
2.2.2	Staff Qualifications.....	78
2.2.3	Electromagnetic Compatibility.....	79
2.2.4	Combination with Medical Electrical Equipment.....	80
2.2.5	Display of Faults.....	80

3	Installation and Initial Operation.....	81
3.1	Safety Notices.....	81
3.2	Staff Qualifications.....	82
3.3	Installation.....	83
3.3.1	Setup.....	83
3.3.2	Connection of Monitors.....	84
3.3.3	Connection of a Keyboard.....	84
3.3.4	Connection of Foot Switches.....	84
3.3.5	Connection of a Recording Device.....	85
3.3.6	Connection of a Light Source.....	86
3.3.7	Connection to Mains Power.....	86
3.3.8	Initial Operation and Function Check.....	87
4	Configuration.....	88
4.1	Basic Settings.....	88
4.1.1	System Language.....	88
4.1.2	Position of the On-Screen Display on the Monitor.....	88
4.1.3	Date and Time, Display Format.....	89
4.1.4	Image Format for Photos.....	89
4.1.5	Image Refresh Rate.....	90
4.1.6	On-Screen Display on the Monitor.....	91
4.1.7	Acoustic Signals during Photo Capture.....	92
4.1.8	Activation of Automatic Light Control.....	92
4.1.9	Shutter Speed.....	93
4.1.10	Foot Switch Configuration.....	94
4.1.11	System Information Display.....	95
4.2	Image Setting.....	95
4.2.1	General Information.....	95
4.2.2	Changing the Settings.....	99
4.2.3	User-Defined Image Setting.....	100
4.3	Reset to Default Settings.....	101
4.3.1	Full Reset.....	101
4.3.2	Reset of an Individual Specialty.....	102

5	Operation.....	103
5.1	Safety Notices.....	103
5.2	Staff Qualifications.....	105
5.3	General Notes on Operation.....	105
5.4	Technical Inspection Prior to Use.....	105
5.4.1	Visual Inspection.....	105
5.4.2	Function Check.....	106
5.5	Switching On and Off, Separation from the Power Supply.....	107
5.6	Preparation for Use.....	107
5.6.1	Connection of the Camera Head.....	107
5.6.2	Selection of a Specialty.....	108
5.6.3	Setting the White Balance.....	108
5.6.4	Entering Patient Data.....	109
5.7	Creation of Recordings.....	110
5.7.1	Creation of Photos.....	110
5.7.2	Creation of Videos.....	110
5.7.3	Creation of Recordings with an External Device.....	111
5.8	Storage and Transport.....	111
6	Cleaning and Disinfection.....	112
6.1	Safety Instructions.....	112
6.2	Performance.....	112
7	Maintenance and Repair.....	113
7.1	Safety Notices.....	113
7.2	Staff Qualifications.....	113
7.3	Safety Inspection.....	113
7.4	Changing of the Battery.....	114
7.5	Changing of the Fuse.....	115
7.6	Repair.....	116
7.7	Troubleshooting Table.....	117
8	Product Data.....	119
8.1	Technical Data.....	119
8.2	Ambient Conditions.....	119
8.3	Spare Parts and Accessories.....	120
8.4	Open-Source Licenses.....	120
9	Disposal.....	121

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0 Important Information about this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following products:

Item name: **Controller Full HD**

Item numbers: **95-3960** - Controller Full HD USB Image Capture

95-3970 - Controller Full HD USB Image/Video Capture

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, medical engineers and employees of central sterile supply who are entrusted with the installation, operation, maintenance, and reprocessing of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

0.4 Supplementary Documents

EMC The electromagnetic compatibility (EMC) requirements are described in the included separate brochure

Medical electrical equipment

pursuant to IEC 60601-1-2

The rules and information contained therein must be observed during installation and operation of this device.

Instructions for use of the camera head

To safely use this device, also observe the instructions for use of the camera head being used.

1 General Information about the Device

1.1 Scope of Delivery

The scope of delivery for the device includes:

- one controller
- one DVI cable
- two 3.5 mm jack plug cables (1.2 m)
- one cable for the MIS-Bus port (0.75 m)
- one power cord (item no. available on request)
- one set of instructions for use (TPA471-000-20)
- one brochure specifying the open-source licenses used (TPO003-121-01)
- one EMC brochure (TPI007-121-20)

Check scope of delivery!

Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Product Description

1.2.1 Performance Characteristics and Function

Medical endoscopy	<p>The present device is a controller for use in medical endoscopy.</p> <p>The combination of the controller and camera head result in a camera based on CMOS technology.</p>
Auto light control	<p>If a compatible light source with automatic light control is connected, the required light intensity can be controlled fully automatic.</p>
Full HD – 1080p	<p>When used with a compatible monitor, the camera delivers a native full HD image resolution using progressive scanning (1080p).</p>
Extensive configuration options	<p>In addition to the general setting options, the camera offers various endoscopic application profiles. These profiles are preconfigured to offer an optimized image display; however, the profiles can also be modified and customized for your specific needs. Ten additional profiles can also be configured.</p> <p>The general setting options include language selection for the configuration menu and monitor displays, and the camera head button settings. The default language setting is English.</p>
Captures	<p>The controller with the item number 95-3960 allows creating and saving photo captures onto a USB storage device, which can be inserted on the front side of the controller. The controller with the item number 95-3970 allows to record video captures in addition to photos.</p> <p>The controller can be used with various camera heads.</p>

1.2.2 Front View of the Controller


Figure 1-1: Front view of the controller.

Fig. 1-1 legend

- ① ON / standby button
- ② Connection for camera head
- ③ Photo capture button / upwards in menu
- ④ White balance button / downwards in menu
- ⑤ Open Menu button / Enter
- ⑥ USB storage device connection
- ⑦ Video capture start / stop button (optional)

1.2.3 Rear View of the Controller


Figure 1-2: Rear view of the controller.

Fig. 1-2 legend

- ① Appliance coupler
- ② Fuse holder
- ③ Main power switch
- ④ MIS-Bus
- ⑤ Recording device connection (Remote-Out)
- ⑥ Potential equalization terminal
- ⑦ Foot switch connection (Remote-In)
- ⑧ Port for service (covered)
- ⑨ 2x DVI output
- ⑩ Battery compartment (Battery VARTA CR2032)
- ⑪ 2x 3G-SDI output
- ⑫ Keyboard connection

1.2.4 Menu - Functions, Access and Navigation

Configuration	The menu is where you configure the device and define the basic settings and image settings.
Selecting the specialty	Furthermore, you select the specialty with the predefined, optimized profile settings and enter the patient data (optional).
Control functions	Perform image settings, white balance or captures before or during procedures, with the appropriate configuration, do so directly with the pre-defined camera head buttons.
Entry and navigation	Entry and navigation in the menu are done through the buttons on the front side of the device, with the keyboard or using the camera head buttons.

Navigation in the menu using the buttons on the front side of the device

	Entry into menu:	Press button briefly.
	Upwards in the menu:	Press button briefly.
	Downwards in the menu:	Press button briefly.
	Confirm selection in menu:	Press button briefly.
	Leave the menu:	Press and hold button.

Navigation in the menu using the keyboard

→	Entry into menu:	Press right arrow button.
↑	Upwards in the menu:	Press up arrow button.
↓	Downwards in the menu:	Press down arrow button.
→	Go one level deeper in menu:	Press right arrow button.
→	Confirm selection in menu:	Press right arrow button.
←	Go one level back in menu:	Press left arrow button.
→	Enter Setup menu:	Press and hold right arrow button.
→	Leave the menu:	Press and hold right arrow button.

1.3 Usage

1.3.1 Intended Use

The product is used to visualize anatomical structures of the human body.

Depending on the camera head connected, the product meets the type BF, CF or CF defibrillation-proof classification requirements according to IEC 60601-1 for protection against electric shock.

This degree of protection must be assured for every connection between the endoscope and connected devices.


Products that meet the requirements classified as type BF feature the adjacent symbol.

Products that meet the type CF classification requirements feature the adjacent symbol, offer an even higher level of protection against electric shock and are approved for direct conductive contact with the heart.

Products that meet the type CF defibrillation-proof classification requirements feature the adjacent symbol, offer an even higher level of protection against electric shock and are defibrillation-proof.

Products that are defibrillation-proof in addition to that are marked with one of the adjacent symbols, depending on the degree of insulation.

1.3.2 Indications

Minimally invasive application.

Surgery, general endoscopy.

1.3.3 Contraindications

The use of the camera components and its accessories is contraindicated if endoscopic procedures are contraindicated for any reason. As is appropriate with any surgical procedure, consideration must be given to patient size and workspace volume when using the camera components. Depending on the patient's disease, there may exist contraindications that rely on the general condition of the patient or the specific disease pattern. The decision to perform an endoscopic procedure rests with the responsible surgeon and should be made on the basis of an individual risk-benefit analysis.

1.4 Conformity

1.4.1 Standards and Directives

The device consisting of controller and camera head meets the requirements of the following applicable standards:

- **IEC 60601-1** Medical electrical equipment - Part 1:
General requirements for basic safety and essential performance
- **IEC 60601-1-2** Medical electrical equipment - Part 1-2:
Electromagnetic compatibility
- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:
Particular requirements for basic safety and essential performance of endoscopic equipment

1.4.2 Medical Device Classification

The product is a Class I medical device within the meaning of applicable **EU guidelines**.

1.5 Marking

1.5.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the product or packaging.

Adhere to the instructions for use

CE marking

Item number

Serial number

Medical Device

Manufacturer

Date of manufacture

Caution (IEC 60601-1 3rd edition) /

Take note of accompanying documents (IEC 60601-1 2nd edition)

MR unsafe

Potential equalization

Permissible storage and transport temperature

Permissible relative air humidity during storage and transport

Permissible atmospheric pressure during storage and transport

Caution: Federal (US) law restricts this device to sale by or on order of a physician

Separate collection for WEEE (waste of electrical and electronic equipment)

Fuse

Alternating current

1.5.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Dangerous electrical voltage warning sign

Biohazard warning sign, risk of infection

1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

1.7 Obligation to Report Serious Incidents

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 General Safety Information
2.1 Safety Messages in this Document
2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

⚠ WARNING

Nature and origin of the injury risk

Potential consequences of non-observance

- > Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

⚠ CAUTION

Nature and origin of the injury risk

Potential consequences of non-observance

- > Preventive action

The following kind of message refers to a risk that could lead to material damage:

NOTICE

Nature and origin of the risk of material damage

Potential consequences of non-observance

- > Preventive action

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Risks despite highest of quality standards

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

Keep the original packaging. Transport and store the product in its original packaging and use it to return goods if service support is required.

If you detect any faults or malfunctions, inform us immediately.

-  **WARNING!** Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.
-  **WARNING!** Operating the device in an open condition. Risk of electric shock. Do not open the device. Operate only when closed.
-  **WARNING!** Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.
-  **WARNING!** Risks from the arrangement, setup, combination, or properties of connected or surrounding devices or equipment. Follow the instructions for use of the respective devices. Perform a risk assessment.
-  **WARNING!** Risk of suffocation. Keep packaging material out of reach of children.
-  **WARNING!** Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.
-  **CAUTION!** Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. Do not kink, crush or strongly bend camera cable. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.
-  **WARNING!** Improper handling, maintenance, and use carries a risk to the patient, user and other persons, or can lead to premature wearing of the device.

2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for installation, operation, and maintenance of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

2.2.3 Electromagnetic Compatibility

Medical electrical equipment is subject to stringent electromagnetic compatibility (EMC) requirements.

Despite the device's high interference immunity and low emitted interference, compliance with EMC-related requirements is necessary with regard to installation, installation location, and ambient conditions.

It is therefore important that you observe the related safety messages in this document and the EMC brochure included with the device.

According to IEC/CISPR 11, medical electrical equipment is divided into groups and classes:

Group 1: The device uses or generates RF-energy only for its internal functioning.

Group 2: The device generates RF-energy in the frequency range 9 kHz to 400 GHz in the form of electromagnetic radiation and through inductive or capacitive coupling for purposes of the treatment of material, its examination or analysis.

Class A: The device **is not suitable** for use in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

Note: The emissions characteristics of this device make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this device might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the device.

The use of the device is only intended for professional healthcare environment.

Class B: The device **is suitable** for operation in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

WARNING! Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

WARNING! Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the device specified by the manufacturer and covered by the scope of this document, including cables. Otherwise, degradation of the performance of the device could result.

WARNING! Interference to the live image caused by external sources (for example portable devices, such as mobile phones). Risk to the patient.

- > Remove all sources of interference
- > Observe minimum distances
- > Ensure that the monitor is correctly displaying the live image

2.2.4 Combination with Medical Electrical Equipment

The device can be combined with components from other manufacturers provided that all components are compliant with the medical electrical equipment safety requirements according to IEC 60601-1.

It is the operator's responsibility to check and make sure that the system is and remains fully operational.

When using devices from different manufacturers and when operating an endoscope and/or endoscopic accessories with electromedical devices, it must be ensured that the applied part is properly insulated (in accordance with the type BF/CF or CF defibrillation-proof requirements defined in IEC 60601-1).

Contact our Service Department if you have any questions regarding compatibility.

2.2.5 Display of Faults

Overheating

The device features an internal temperature sensor.

In the event of controller overheating, the message "Excess temperature! Attention, Camera will switch off in two minutes!" appears on the screen.

If the menu is opened, it closes automatically when the message appears.

Incompatible USB storage device

Compatible USB storage devices must be formatted with the FAT32 file system, and may only contain image and video files.

The USB storage device must meet the USB 2.0 standard (or higher) and be capable of a minimum data transfer rate of 20 MB/s.

If an incompatible USB storage device is used when attempting to save a capture, the message "USB storage device in incorrect format" appears.

It is recommended to only use USB flash drives in conjunction with the devices listed under the scope of validity ([section 0.1](#)) of these instructions for use. USB hard drives (e.g. HDD or SSD) are not suitable.

To ensure compatibility of the USB storage device, it is recommended to make a test capture. Please note that the devices have no formatting capabilities.

3 Installation and Initial Operation

3.1 Safety Notices

WARNING

Electrical connections installed improperly

Risk of fire, short circuit, or electric shock

- > Make sure that the electrical connections are installed in accordance with the relevant national technical regulations

WARNING

Medical electrical system installation

Risk of fire, short circuit, or electric shock

- > Installation must only be performed by qualified personnel
- > If a number of electrical devices are being used in combination, observe IEC 60601-1
- > Non-medical devices that meet the applicable IEC safety standards can be connected only via a medical isolating transformer
- > Do not connect any additional non-medical devices to a medical electrical system
- > Signal lines from devices with a functional ground (FE) connection, which are connected to different branches of the mains power supply, must be galvanically isolated at both ends
- > Connect the devices only to a power supply with a protective ground (PE)
- > After installing a medical electrical system, it must be inspected according to IEC 62353

WARNING

Electromagnetic interactions

Possible impaired function and image quality

- > Observe the EMC brochure supplied
- > Make sure that neighboring devices are not causing reciprocal electromagnetic interference
- > Perform a post-installation function check

**⚠ WARNING****Use of multiple sockets**

Risk of fire, short circuit, electrical shock, reduced level of safety

- > Whenever possible, avoid the use of multiple sockets
- > If required, use medically approved multiple sockets
- > Never connect multiple sockets in series
- > Do not cover multiple sockets (leads to heat accumulation)
- > Do not place multiple sockets on the floor
- > Use traction relief
- > Devices should be connected to a shared multiple socket only if they are all intended to be used as part of the medical electrical system

3.2 Staff Qualifications

Qualified personnel Personnel responsible for setting up or installing medical electrical equipment must possess the necessary qualifications and be familiar with the safety regulations in place at the installation location and the applicable official regulations.

3.3 Installation

3.3.1 Setup

Setting up the device

Observe the safety notices provided at the start of this chapter and refer to the enclosed EMC brochure.

Proceed as follows:

WARNING! Risk from installation in potentially explosive locations. Increased risk of fire and explosion in oxygen-enriched atmospheres. Device must be installed outside of potentially explosive locations and not in the vicinity of flammable substances.

WARNING! Risk resulting from electromagnetism. Possible malfunctions and image interference. Use only the cables supplied. Ensure that all devices operated in the vicinity meet the EMC requirements. Perform a post-installation function check.

WARNING! Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other devices should be observed to verify that they are operating normally.

1. Install the device.

Install safely!

When installing the device, ensure that:

- it is positioned on a level, non-slip surface that is able to support its weight,
- it is installed outside of the patient environment and the sterile area on a sufficiently stable support,
- it is protected from water drips and splashes,
- it will not be subject to vibration during operation,
- the ventilation slots are not obstructed,
- the main switch on the rear of the device is accessible.

2. Connect the **potential equalization terminal** on the rear of the device to the potential equalization terminal at the installation site.

Use a potential equalization conductor in accordance with DIN 42801 with at least a 4 mm² cross-section.

- The device is properly installed.

3.3.2 Connection of Monitors

Connecting monitors

Multiple monitors can be connected. The main monitor should be connected to the best signal available. Please observe the instructions for use of the monitor.

Proceed as follows:

1. Connect one of the **DVI** ports on the rear of the device to the DVI input on the main monitor.
Use the supplied DVI cable.
 2. If necessary, connect any additional monitors to the remaining video outputs.
- Secure the cable Secure the connection cables at both ends to prevent them from being pulled out inadvertently.
- Route cables safely! 3. When routing the cables, make sure they do not present a trip or entanglement hazard.
- The monitors are connected properly.

3.3.3 Connection of a Keyboard

Connecting a keyboard

Proceed as follows:

1. Connect the keyboard connection cable to the USB connection labeled **USB-Keyboard (USB keyboard)** on the rear of the device.
 2. When routing the cables, make sure they do not present a trip or entanglement hazard.
- Route cables safely! The keyboard has been properly connected.

3.3.4 Connection of Foot Switches

Connecting foot switches

Up to two single-pedal foot switches can be connected. Connect foot switches if you want to control the camera by foot switch.

Note: Only the foot switch specified in [section 8.3](#) is compatible.

Proceed as follows:

1. Connect a foot switch to the **Remote-In 1** connection of the controller.
 2. Connect a second foot switch to the **Remote-In 2** connection of the controller.
 3. Configure the foot switches as described in [section 4.1.10](#).
 4. When routing the cables, make sure they do not present a trip or entanglement hazard.
- Configure the foot switches The foot switches are connected properly.
- Route cables safely!

3.3.5 Connection of a Recording Device

Connecting a recording device

Observe the instructions for use supplied with the recording device.

Proceed as follows:

1. Depending on the type of signal supported by the recording device, connect either the **DVI**, or an **3G-SDI** connection on the rear of the medical device to the corresponding input on the recording device.
2. If you would like to use the buttons on the camera head or a foot switch to make recordings, connect the Remote-Out ports on the rear of the medical device to the corresponding inputs on the recording device.

Use the supplied 3.5 mm jack plug cable to do this.

Remote-Out 1
for photos

Connect the **Remote-Out 1** port to the photo recording connection on the recording device.

Remote-Out 2
for video recordings

Connect the **Remote-Out 2** port to the video recording connection on the recording device.

3. If necessary, modify the button layout of the camera head.

Route cables safely!

4. When routing the cables, make sure they do not present a trip or entanglement hazard.

The recording device has been properly connected.

3.3.6 Connection of a Light Source

Connecting the light source

The connected light source can be controlled using the buttons on the camera head.

Note: Only compatible light sources can be controlled with the camera head buttons.

Proceed as follows:

1. Connect the MIS-Bus port on the controller to that on the light source. Use the supplied MIS-Bus cable.
 2. Ensure the connections are secured on both ends.
 3. If necessary, modify the button layout of the camera head so that they can be used to switch the light source on and off or to control its brightness.
- The light source has been properly connected.

Configure camera head buttons

3.3.7 Connection to Mains Power

Connecting the device to the mains

Do not connect the device to the mains until you have completed all the previous activities. Observe the safety notices provided at the start of this chapter.

The use of an uninterruptible power supply (UPS) is recommended to protect against input power interruptions.

NOTICE! Risk of condensation with extreme differences in temperature. A device stored in very cold conditions can result in condensation in the device from humidity in the air when there is a rapid change in temperature from cold to warm. Wait until the device has come to the ambient temperature to connect to power supply and switch on.

Proceed as follows:

1. Check that the supply voltage at the site of operation corresponds to that specified on the device nameplate.
2. Insert the power cord into the mains connection socket on the rear of the device.

WARNING! Sparks in a potentially explosive location. Risk of explosion. Ensure that the device is not situated in a potentially explosive location when connecting to the power supply.

3. Connect the power cord to the power supply.
 4. When routing the cable, make sure it does not present a trip or entanglement hazard.
- The device is properly connected to the power.

Route cables safely!

3.3.8 Initial Operation and Function Check**Performing a function check**

When operating the device for the first time, check the installation and make sure the device is in full working order.

Precondition: Installation must have been completed.

Proceed as follows:

- | | |
|---------------------|---|
| | 1. Switch ON at the main switch on the rear of the device. |
| | > <i>The device switches to standby mode.</i> |
| | 2. Switch on the device and all neighboring devices. |
| Connect camera head | 3. Connect a camera head as described in section 5.6.1 . |
| Set white balance | 4. Perform white balance as described in section 5.6.3 . |
| Focus the image | 5. Direct the camera head at an object in the room and focus the image. |
| | 6. Check that a good image is displayed on the monitor without vertical lines, color variations, or flickering. |

**Problems with the image displayed on screen?**

Electromagnetism can cause interference.

Try rearranging the different devices to increase the distance between them.
Contact our Service Department if you require further assistance.

7. Check that the buttons on the camera head can be operated properly.
- This concludes the function check.

4 Configuration

4.1 Basic Settings

4.1.1 System Language

Setting system language

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
 - > The main menu appears on the screen.
 2. Navigate to the **Options** menu entry.
 3. Press the **Menu** button for approx. 2 seconds.
 - > The setup menu appears on the screen.
 4. Navigate to **Configuration** and select.
 5. Navigate to **Language selection** and select.
 6. Select the desired language and confirm.
 - > Your selection is displayed in white type.
 - > The language is changed as desired.
- The system language is set.

Language selection

4.1.2 Position of the On-Screen Display on the Monitor

Changing the on-screen display position

To change the position of the on-screen display on the monitor, choose an aspect ratio of 16:9 or 16:10, depending on the monitor used.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
 - > The main menu appears on the screen.
2. Navigate to the **Options** menu entry.
3. Press the **Menu** button for approx. 2 seconds.
 - > The setup menu appears on the screen.
4. Navigate to **Configuration** and select.
5. Navigate to **16:9 / 16:10** and select.
 - > Your selection is displayed in white type.

16:9 / 16:10

6. Select the desired aspect ratio and confirm.
 The on-screen display position is changed.

4.1.3 Date and Time, Display Format

Setting date and time

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> The main menu appears on the screen.
2. Navigate to the **Options** menu entry.
3. Press the **Menu** button for approx. 2 seconds.
> The setup menu appears on the screen.
4. Navigate to **Configuration** and select.
5. Navigate to **Date/time** and select.
6. Select the desired display format and confirm.
> Your selection is displayed in white type.
7. Set the date and time and save in each case by pressing Enter.
 Date and time are set.

Date/time

4.1.4 Image Format for Photos

Defining image format (optional)

Define what file format should be used for photos taken. Select between BMP and JPEG.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> The main menu appears on the screen.
2. Navigate to the **Options** menu entry.
3. Press the **Menu** button for approx. 2 seconds.
> The setup menu appears on the screen.

4. Navigate to **Configuration** and select.
- Image format 5. Navigate to **Image Format** and select.
6. Select the desired file format and confirm.
Select BMP if you want to create uncompressed photos in higher quality. One image can be up to 6 MB in size.
Select JPG if you want to create compressed photos in good quality with a smaller file size of approx. 400 KB.
> Your selection is displayed in white type.
- The image format has been set.

4.1.5 Image Refresh Rate

Defining image refresh rate

Leave the factory-preset image refresh rate unchanged. This setting only needs to be changed in exceptional circumstances.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> *The main menu appears on the screen.*
2. Navigate to the **Options** menu entry.
3. Press the **Menu** button for approx. 2 seconds.
> *The setup menu appears on the screen.*
4. Navigate to **Configuration** and select.
- 50Hz / 60Hz 5. Navigate to **50Hz / 60Hz** and select.
6. Select the desired image refresh rate and confirm.
> Your selection is displayed in white type.
> A dialog window for confirmation opens up.
7. Confirm the dialog window using **Perform**.
- The image refresh rate is set.

4.1.6 On-Screen Display on the Monitor

Defining on-screen display on the monitor

Determine whether the date, time and if necessary patient data should be overlaid onto the screen. Also determine whether the configuration of the camera head buttons should be displayed temporarily after switching on the camera system.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
 - > *The main menu appears on the screen.*
2. Navigate to the **Options** menu entry.
3. Press the **Menu** button for approx. 2 seconds.
 - > *The setup menu appears on the screen.*
4. Navigate to **Configuration** and select.
- OSD 5. Navigate to **OSD** and select.
6. Select the desired on-screen displays and confirm.
 - > Your selection is displayed in white type.

When selected,

- the date and time are overlaid onto the screen,
 - the patient data are overlaid onto the screen,
 - the configuration of the camera head buttons is displayed after switching on the camera system for several seconds.
- The on-screen displays on the monitor are set.

4.1.7 Acoustic Signals during Photo Capture

Switching signal sounds on or off (optional)

Determine whether an acoustic signal should confirm when photos are taken.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> The main menu appears on the screen.

2. Navigate to the **Options** menu entry.

3. Press the **Menu** button for approx. 2 seconds.
> The setup menu appears on the screen.

4. Navigate to **Configuration** and select.

Beeper 5. Navigate to **Beeper** and select.

6. Activate or deactivate the signal sound and confirm.
> Your selection is displayed in white type.

The signal sound is set.

4.1.8 Activation of Automatic Light Control

Compatible light sources allow for the use of the auto light control feature. To this purpose, it is necessary to activate auto light control.

Activating automatic light control

Proceed as follows:

1. Using the MIS-Bus cable, connect the camera to a compatible light source. Doing so, observe the instructions for use supplied with the light source.

2. Activate auto light control via the camera menu.

> When auto light control is activated, this message appears on the light source display.

Auto light control is now activated.

4.1.9 Shutter Speed

Setting shutter speed

The shutter speed is normally set automatically. Setting a fixed shutter speed can be beneficial for special applications, e.g. for minimizing the effects of vibration during stroboscopy.

Set a fixed shutter speed if required.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> *The main menu appears on the screen.*

2. Navigate to the **Options** menu entry.

3. Press the **Menu** button for approx. 2 seconds.
> *The setup menu appears on the screen.*

4. Navigate to **Configuration** and select.

Shutter 5. Navigate to **Shutter** and select.

6. Select "Auto" or "Manual" and confirm.

If you choose "Manual", set the desired shutter speed and confirm.

> Your selection is displayed in white type.

The shutter speed is set.

4.1.10 Foot Switch Configuration

Configure foot switches

Assign a function to foot switches as desired.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
 - > The main menu appears on the screen.
2. Navigate to the **Options** menu entry.
3. Press the **Menu** button for approx. 2 seconds.
 - > The setup menu appears on the screen.
4. Navigate to **Configuration** and select.

Foot switch

5. Navigate to **Foot switch** and select.
6. Assign the foot switch function as desired.
Select one of the following functions:
 - Selective color enhancement (SCE)
 - Smoke reduction
 - Grid removal
 - Light source on/off
 - Photo [capture] (onto a USB storage device)
 - Video [capture] (onto a USB storage device)
 - Remote-Out 1 (capturing of photos through an external recording device)
 - Remote-Out 2 (capturing of videos through an external recording device)

Note: The selectable image settings are described in [section 4.2.1](#) in more detail.

- The foot switch is configured.

4.1.11 System Information Display

Displaying system information

The system information can be useful to the Service Department.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> The main menu appears on the screen.
2. Navigate to the **Options** menu entry.
3. Press the **Menu** button for approx. 2 seconds.
> The setup menu appears on the screen.

System information

4. Navigate to **System information** and select.
 5. The desired information can be displayed:
 - Software version
 - Hardware version
 - Serial number > Controller, Serial number > Camera head and Serial number > Light source
- The system information was displayed.

4.2 Image Setting

4.2.1 General Information

Optimized camera settings can be selected in the configuration menu for selected endoscopic specialties.

The image settings can be adjusted and saved according to your requirements.

Adjustment of image settings

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> The main menu appears on the screen.
2. To reach the image setting you wish to adjust, navigate to the menu item **Settings** and select.
> A list of the image settings described in the following opens up.

3. Navigate to the desired image setting and select.
 - > An additional window appears on the screen next to the main menu.
4. Make the desired adjustments using the arrow buttons.

Save your changes by pressing the Enter button each time.

- The image settings were changed for the desired specialty.

The following configurable parameters are available for each specialty:

Brightness

Brightness controls the exposure control for the camera.

Reduce brightness to decrease blooming and increase it to make the image brighter.

Contrast

Contrast changes the contrast level of the image.

Reduce the contrast to obtain a homogeneously illuminated image for endoscopes with poorly illuminated edges.

Note: Dimming on the edges of the image can also result from the window size adjustment being too small.

Edge enhancement

Edge enhancement amplifies the contours at the transitions of adjacent structures.

Increase the edge enhancement during laparoscopic and arthroscopic applications with rod lens systems to improving the view of vessels.

Reduce the edge enhancement with image bundle systems to decrease unwanted structures of the image bundle system.

Note: Edge enhancement can increase image noise.

Color saturation

Saturation changes color intensity.

Change the saturation according to your personal impression. Reduce the saturation to decrease unwanted color effects, such as moiré effects.

Window

Window size defines the image section used by the camera for controlling image brightness.

Choose a window that is similar in size or smaller than the round endoscope image.

Choose an even smaller window to achieve a better view in the distance in tubular cavities, such as the urethra.

Note: If the window setting is too large and exceeds the round endoscope image, blooming can result.

Light source +/-

This menu item sets the brightness of the light source, if it is compatible and appropriately connected.

Select a light intensity of 70 %-100 % with an LED light source.

Reduce the light intensity to reduce blooming. Increase the light intensity to reduce image noise.

Note: With LED light sources, too low a light intensity can change color representation.

Light source on/off

This menu item allows you to switch the light source on and off, if it is compatible and appropriately connected.

Auto light control

Compatible light sources allow for the use of the auto light control feature. To this purpose, it is necessary to activate auto light control ([section 4.1.8](#)).

Zoom

The digital zoom allows up to three-fold digital enlargement of the image displayed.

Note: The digital zoom reduces the image resolution and decreases the image quality. For this reason, use the optical zoom on the camera head preferably whenever possible.

Selective color enhancement (SCE)

With selective color enhancement, you can intensify the predominant color in the image.

Use selective color enhancement to better differentiate structures from one another in environments with little differences in color, for example, to accentuate the colors of vessels.

Note: Selective color enhancement changes the natural color representation. For this reason, the message “Attention: false color rendering” is displayed on the screen as long as this setting is activated.

Smoke reduction

Smoke reduction allows a clearer image in a smoke-filled environment (e.g. when there is coagulation with electrosurgical units).

Use smoke reduction to compensate for a poor view due to smoke formation.

Note: Smoke reduction creates a very contrast-rich image and can promote blooming.

Color Shift 1 / 2

Color Shift 1 / 2 changes the color of vessels and tissue structures.

Use Color Shift 1 / 2 for a better view of vessels and tissue structures.

Note: Color Shift 1 / 2 alters the natural color representation. For this reason, the message “Color Shift 1” or “Color Shift 2” is displayed on the screen as long as this setting is activated.

Grid removal

Grid removal reduces or removes the grid that is caused by use of image bundle systems.

Use grid removal to suppress interfering grids.

Note: Grid removal causes a slight reduction in image sharpness.

Noise reduction

Noise reduction reduces noise in the image.

Adjust noise reduction at the end or when a high edge enhancement has been set and further reduction of noise is desired.

Note: Noise reduction leads to a slight reduction in image sharpness starting at the middle setting. A higher setting may cause artifacts with quick movements (streaking).

Picture out of Picture (PoP)

Picture out of Picture displays the image adjusted by video algorithms on the right-hand side of the original image.

The menu item Picture out of Picture allows you to choose the video algorithm to be used for the display of both images next to one another for comparison purposes.

Note: The original image is displayed on the left, the adjusted image is displayed on the right. When selecting Selective color enhancement (SCE), Smoke reduction or Color Shift 1 / 2, both live images show the same image section. When selecting Zoom, the center of the original image is displayed in enlarged format.

4.2.2 Changing the Settings

Changing image settings for a specialty

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> The main menu appears on the screen.

2. To reach the specialty you wish to change, navigate to the menu item **Specialty** and select.

Navigate to the desired specialty and select.

> The specialty is changed.

3. Navigate to **Settings** and select.
4. Make the desired changes.

Save your changes by pressing the Enter button each time.

5. Set the white balance for a realistic view of the changes you have made.

The image settings were changed for the desired specialty.

4.2.3 User-Defined Image Setting

Making and saving user-defined image settings as a profile

Create a user-defined profile based on the existing reference profile and save it under a name of your choosing.

Note: A keyboard must be connected to enter a user-defined name.

User-defined profiles can also be saved without a keyboard. In this case, your settings are saved under the preset name, consisting of the camera head model and User 1-10 (e.g. FVHD2 User 1, FVHD2 User 2, ...).

Up to 10 user-defined profiles can be saved.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> The main menu appears on the screen.
2. Navigate to the **Specialty** menu item and select.
3. Navigate to the **Edit user profile** menu item and select.
4. Navigate to the **Rename** menu item and select.

Enter name

5. Press the right arrow button on the keyboard and enter a name for the user-defined profile.

Up to 30 characters are available for the name.

6. Confirm the entry with the Enter button on the keyboard.

> The user-defined profile is renamed with the name given.

7. Navigate back in the main menu to **Reference specialty** and select.

Select reference profile

8. Select the reference profile that you wish to base the user-defined profile on and confirm.

9. Navigate back in the main menu to **Settings** and select.

10. Make the desired changes.

Save your changes by pressing the Enter button each time.

- The user-defined profile is saved.

4.3 Reset to Default Settings

4.3.1 Full Reset

Reset the device completely to default settings

NOTICE! All settings are reset to the default settings. If necessary, only restore the default settings of an individual specialty.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
> *The main menu appears on the screen.*

2. Navigate to the **Options** menu entry.

3. Press the **Menu** button for approx. 2 seconds.

> *The setup menu appears on the monitor.*

4. Navigate to **Configuration** and select.

Default settings

5. Navigate to **Default settings** and select.

6. Select **Perform** and confirm.

> *A dialog window for confirmation is displayed on the screen. The option **Perform** is activated.*

*If you wish to cancel the process, select **Cancel** and confirm.*

> *The dialog window closes and the setup menu is shown.*

If you wish to reset the device to the default settings, confirm the selection.

> *The device is reset to the default settings and restarts.*

The device has been reset to the default settings.

4.3.2 Reset of an Individual Specialty

Resetting the image settings of a specialty

Follow the steps below if you want to reset the image settings of an individual specialty to the default settings.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Select the specialty you want to reset, as described in [section 5.6.2](#).

Settings

3. Navigate to **Settings** and select.
4. Navigate to **Reset?** and select.
5. Navigate to **Perform** and select.

The image settings of the specialty have been reset.

5 Operation

5.1 Safety Notices

WARNING

Modifications to the installation

Risk of fire, short circuit, and electric shock

- > Follow the safety notices in the Installation and Initial Operation chapter
- > Do not make unauthorized changes to the electrical installation
- > If the installation contains a multiple socket, do not connect additional devices to this without authorization
- > Never connect multiple sockets in series

WARNING

Current leakage when touching the patient

Risk of electric shock

- > Never touch the patient and the device housing or accessible contacts of electric devices at the same time
- > Never touch the patient and non-medical electric devices at the same time

WARNING

Positioning the equipment cart inside the sterile area

Risk of infection

- > Do not position the equipment cart inside the sterile area
- > Measure the cable lengths so there is enough room for movement without compromising sterility
- > The light guide must be anchored securely within the direct OR environment to avoid any risk of accidental slipping

WARNING

Displaying a recorded image instead of the live image

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera at all times

⚠ WARNING

Live image failure during operation of a defibrillator

Risk to the patient

- > Make sure that 3 seconds (recovery time) after operation of a defibrillator the live image is displayed correctly again
- > Restart the controller if the live image is incorrectly displayed

⚠ WARNING

Interactions between devices in simultaneous use (e.g. lasers, electro-surgery)

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the required type BF, CF or CF defibrillation-proof classification requirements according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active RF electrodes
- > Do not activate RF electrodes in the vicinity of flammable gases or liquids
- > Before using RF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

⚠ WARNING

Addition of leakage currents. When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive

Risk to the patient

- > Use endoscopes in conjunction with endotherapy devices only if both have matching type applied parts (CF or CF defibrillation-proof)

⚠ CAUTION

Electrostatic discharge at the connection field of the monitor

Temporary image failure, risk to patient

- > Do not touch the connection field and connections at the rear of the monitor during ongoing procedures

5.2 Staff Qualifications

This device may be operated only by physicians and medical assistants with the requisite training, advanced training, knowledge, and practical experience in the endoscopic discipline concerned as defined by the provisions in place at the site of operation.

5.3 General Notes on Operation

Start every procedure by creating a new patient. This ensures that the captures created during the procedure are saved in a separate file folder.

Enter the patient data for this as described in [section 5.6.4](#).

New patient Find the **“New patient”** entry in the menu. Confirming this entry deletes the existing patient data. If you do not wish to enter any patient data, confirm this entry anyway so that a new file folder is generated on the USB storage device for the subsequent captures.

Sufficient disk space? Ensure that a USB storage device is connected with sufficient disk space for your recordings.

Depending on the components used, navigation in the menu can be done with the buttons on the front side of the device, the camera head buttons or via a connected keyboard.

Operation on the front side of the device is described in the following.

5.4 Technical Inspection Prior to Use

5.4.1 Visual Inspection

Performing a visual inspection

*A visual inspection must be performed **prior to each medical procedure**.*

Proceed as follows:

- | | |
|----------------------------------|--|
| Properly installed? | <ol style="list-style-type: none"> 1. Check that the device is installed as described in section 3.3.1. Make sure in particular that the potential equalization is properly connected.

Do not use the device if it is not properly set up and connected. |
| Is the housing undamaged? | <ol style="list-style-type: none"> 2. Inspect the housing for external damage.

Do not use the device if there is any external damage to the housing. |
| Is the cable undamaged? | <ol style="list-style-type: none"> 3. Check that the mains cable is undamaged.

Do not use the device if the cable is damaged. |
| Is the socket dry and undamaged? | <ol style="list-style-type: none"> 4. Check that the contacts on the camera cable connection socket are free from moisture and dirt.

Do not use the device if there is any moisture in the socket or contacts are soiled. <p><input checked="" type="checkbox"/> This concludes the visual inspection.</p> |

5.4.2 Function Check

Performing a function check

*A function check must be performed **prior to each medical procedure.***

Proceed as follows:

- | | |
|----------------------|--|
| Connect camera head | 1. Connect a camera head as described in section 5.6.1 . |
| Switch on the system | 2. Switch on all the system components that you intend using for the procedure. |
| Set white balance | 3. Perform white balance as described in section 5.6.3 . |
| Focus the image | 4. Direct the camera head at an object in the room and focus the image.
Check that you are able to achieve a sharp, good-quality image.
Do not use the system if the image contains vertical lines, color variations, or flickers, or you are unable to achieve a sharp, good-quality image. |
| | 5. Adjust the monitor as required to achieve an optimum viewing angle for the surgeon. |
| | <input checked="" type="checkbox"/> This concludes the function check. |

5.5 Switching On and Off, Separation from the Power Supply

Switching on the device

Proceed as follows:

1. Switch on the device at the on/off button on the rear.
 - > *The device switches to standby mode.*

2. Press the ON button.
 - > *The ON switch illuminates in white.*

- The device is on.

Switching off the device

Proceed as follows:

1. Switch off the device at the front.
 - > *The device switches to standby mode.*

To switch off the device completely, switch off at the main switch on the rear.

- The device is off.

Separating the device from power

Proceed as follows:

1. To separate the device completely from its power supply, disconnect the power cord. The appliance coupler is considered as disconnecting device.

Make sure that the rear of the device is accessible at any time.

- The device is completely separated from the power supply.

5.6 Preparation for Use

5.6.1 Connection of the Camera Head

Connecting up the camera head

NOTICE! Damage to the camera cable by kinking, crushing or sharp bending. Do not kink, crush or strongly bend camera cable.

Proceed as follows:

1. Slide the connection plug fully into the connection socket on the device until it clicks into place.

> *The name of the connected camera head model is displayed on the screen.*

- The camera head is connected.

5.6.2 Selection of a Specialty

Selecting a specialty

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
 - > The main menu appears on the screen.

2. Navigate to the **Specialty** menu item and select.

The following optimized camera settings are available, for example, depending on the camera head model connected:

- Laparoscopy
- Urology
- URO Fiberscope
- Gynecology
- GYN Fiberscope
- Arthroscopy
- ENT rigid
- ENT flexible
- Neuroendoscopy

Navigate to the desired specialty and confirm.

> The parameters of the current specialty are taken over.

- The specialty is selected.

5.6.3 Setting the White Balance

Setting white balance

Set the white balance **prior to each medical procedure**.

Note: Keep in mind that ambient light that deviates significantly from the color temperature of the light source can affect the results of the white balance.

Proceed as follows:

1. Connect the endoscope that you intend on using for the procedure.
2. Connect the light guide and light source that you intend on using for the procedure.
3. Switch on all system components.
- 5 cm away 4. Direct the endoscope at a white object, a sterile wipe for example, that is approximately **5 cm** away and focus the image.

5. Press the **White balance** button on the front of the controller for about 2 seconds or use the designated button on the camera head.

> *The message **White balance OK** appears on the screen.*

Does the message "White balance failed" appear on the screen?

The white balance was not performed as described above.

Repeat the white balance.

- The white balance has been set.

5.6.4 Entering Patient Data

Entering a patient data

A keyboard must be connected to enter a patient data.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

Proceed as follows:

1. Press the **Menu** button on the front side of the device for approx. 2 seconds.
 - > *The main menu appears on the screen.*
2. Navigate to the **Patient** menu item and select.
 - > *The menu item "New patient" is activated.*
3. Press the right arrow button to delete the existing patient data.
4. Navigate to the first field for data entry using the down arrow button.
5. Enter the data. Jump from field to field with the **Enter** button.
6. Navigate to **Confirm selection** and save your entries using the right arrow button.
 - > *With the first recording, a new file folder is generated on the USB storage device, in which the captures are saved.*

- The patient data are entered.

5.7 Creation of Recordings

CAUTION! Recorded photos and videos can be compressed or changed by modified image settings. The photo and video documentation is not suitable for making a diagnosis.

5.7.1 Creation of Photos

Creating a photo

Proceed as follows:

1. Ensure that a USB storage device is connected.
2. Activate the release.

Dependent on model, configuration and connected accessories

- press the **Photo capture** button on the front of the device,
 - press the affected button on the camera head,
 - activate the affected foot switch,
 - release the capture in the menu.
- > *The indicator on the main screen shows the photo captured for approx. one half second (freeze function).*
- > *The photo is saved to the USB storage device.*

The photo is recorded.

5.7.2 Creation of Videos

Creating video recordings (95-3970 only)

Proceed as follows:

1. Ensure that a USB storage device is connected.
2. Activate the release to start the recording.

Dependent on model, configuration and connected accessories

- press the **Video capture** button on the front of the device,
 - press the affected button on the camera head,
 - activate the affected foot switch,
 - release the capture in the menu.
- > *The main screen shows a symbol indicating that the video capture has been started. The video is saved to the USB storage device.*

3. Activate the release once more to exit the recording.

The video is created.

5.7.3 Creation of Recordings with an External Device

Creating photo or video recordings

The recording and ending of videos can only be controlled with the camera head buttons or via foot switch. The camera head or foot switch must be configured appropriately. For video capture, a recording device must be connected.

Observe the instructions for use supplied with the recording device.

Proceed as follows:

1. Ensure that a recording device is connected and switched on.
2. Activate the release to start the recording.

Dependent on configuration and connected accessories

- press the affected button on the camera head,
- activate the affected foot switch.

> *The photo or video recording is saved to the external device.*

3. Video recording only: Activate the release once more to exit the recording.
 The capture is saved.

5.8 Storage and Transport

Store the device in a dry, well-ventilated, temperature-controlled room, where it is protected from dust.

Switch off the device prior to storage and remove the power cord and accessories.

When placing into storage or storing temporarily, make sure that the device cannot be damaged by its surroundings.

To avoid damage, keep out of direct sunlight and away from radioactivity and strong electromagnetic radiation.

The ambient conditions for transport, storage, and operation are described in [section 8.2](#).

6 Cleaning and Disinfection

6.1 Safety Instructions

⚠ WARNING

Risks associated with improper cleaning

Risk of fire, short circuit, and electric shock

- > Before cleaning, switch off the device at its rear main switch and disconnect from the mains power supply
- > Make sure no moisture penetrates the device
- > The device should not be connected to the mains power supply until it has completely dried

NOTICE

Improper cleaning

Risk of damage to the device

- > Do not use sponges or wipes that are dripping wet
- > Do not use abrasive cleaners or solvents

6.2 Performance

Cleaning the device

Proceed as follows:

Separate the device
from power

1. Switch off the device at its rear main switch and remove the power cord.

NOTICE! Pulling the power cable will damage it. Disconnect the power at the plug.

2. Clean the device.

Use slightly moistened
wipe only

Wipe the external surfaces with a soft wipe.

Only slightly moisten the wipe with a mild soap solution or isopropyl alcohol.

Make sure that the contacts of the plug do not get in touch with humidity.

3. Connect the device to the power supply and switch back on at the main switch.

This concludes the cleaning process.

7 Maintenance and Repair

7.1 Safety Notices

⚠ WARNING

Improper maintenance and repair

Risk of fire, short circuit, and electric shock

- > All maintenance and repair work must only be performed by qualified personnel
- > Do not make any modifications to the device

7.2 Staff Qualifications

Qualified personnel Personnel responsible for maintaining the device or carrying out safety inspections must possess the necessary qualifications and be familiar with the local safety regulations in place at the installation location and the applicable official regulations.

Repairs must be performed by authorized service technicians only.

7.3 Safety Inspection

We recommend to arrange for safety inspections to be performed every 12 months in accordance with Section 5 of IEC 62353.

The inspection comprises the following steps:

- visual inspection
- measuring the protective resistance (only for protection class I devices)
- measuring the leakage currents
- measuring the insulation resistance

Document results The results of the inspection should be documented in an inspection report in accordance with Section 6 and Annex G of IEC 62353.

We recommend that all components of a medical electrical system be subjected to documented technical safety inspections at regular intervals.

7.4 Changing of the Battery

Changing the battery

Proceed as follows:

1. Isolate the device from its power supply and disconnect the power cord.

NOTICE! Pulling the power cable will damage it. Disconnect the power at the plug.

2. Loosen the screws of the battery compartment.


Figure 7-1: Changing of the battery.

3. Pull out the battery compartment.
 4. Replace the battery (VARTA CR2032).
 5. Gently slide the battery compartment back in and secure it with the screws.
- This concludes the changing of the battery process.

7.5 Changing of the Fuse

Changing the fuse

Proceed as follows:

1. Isolate the device from its power supply and disconnect the power cord.

NOTICE! Pulling the power cable will damage it. Disconnect the power at the plug.

2. Release the fuse holder with a screwdriver.


Figure 7-2: Releasing the fuse holder with a screwdriver.

3. Carefully pull out the fuse holder.

WARNING! Unsuitable electrical fuse. Risk of fire. Insert the fuse specified in [section 8.3](#) only and never bypass the fuse.

4. Replace the faulty fuse.
 5. Carefully slide the fuse holder back in until it clicks into place.
- This concludes the fuse replacement.

7.6 Repair

Should you need to arrange a repair for the device, please contact one of our subsidiaries. Contact details can be found on the back of these instructions for use.

When you send in equipment, please enclose as accurate a **fault description** as possible, and record the item number and serial number of the product on the delivery note. You will find these details on the nameplate.

WARNING

Contaminated device

Risk of infection

- > Reprocess the device properly prior to shipping (*chapter 6*)
- > Clearly identify contaminated products as such

Send only thoroughly cleaned goods in for repair.

For shipping, please use the original packaging whenever possible. Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

7.7 Troubleshooting Table

Issue	Possible causes	Fixes
On-screen display set off	Aspect ratio set incorrectly	Correct the position of the on-screen display (section 4.1.2)
Image too dark, too little illumination	Glass surfaces on camera head are soiled	Clean glass surfaces
	Stubborn residue on the glass surfaces	Remove residue
	Light guide defective	Connect a new light guide. Send in for repair
	Camera brightness set too low	Increase camera brightness
	Endoscope optical system defective	If the image is too dark even without the camera, use a different endoscope and send in the endoscope for repair
	Light output set too dark	Increase the light output at the light source
No image on the monitor	No power supply	Connect and switch on devices (chapter 3)
	Fuse defective	Replace fuse (chapter 7)
	Controller incorrectly connected to monitor	Connect controller to monitor correctly (chapter 3)
	Video cable defective	Replace video cable
	Connection cable not connected or defective	Connect the camera head to the controller. Check the plug for moisture. If the cable is defective, send it in for repair
	Monitor not tuned in to the correct input signal	Set monitor for the correct input signal
Date / time not correct	Date / time not correctly set	Set the date / time correctly in the operating menu of the camera
	Battery empty or defective	Change the battery (section 7.4)
Light source switching on/off on its own	MIS-Bus fault in the camera	Disconnect camera and light source, send in camera for repair
	MIS-Bus fault in the light source	Disconnect camera and light source, send in light source for repair
Auto light control not possible	Light source does not support auto light control	Use a compatible light source model

Table 7-1: Troubleshooting table.

Message	Possible causes	Fixes
“Camera head not connected”	Camera head not connected to controller	Connect compatible camera head to controller. For combination options see the respective instructions for use of the camera heads
“Incompatible camera head”	Camera head not compatible with controller	
“Excess temperature! Attention, Camera will switch off in two minutes!”	Controller ventilation slots are covered	Uncover ventilation slots
	Controller fan defective	Send controller in for repair
“No light source connected / Light source incompatible”	Light source not compatible with controller	Connect compatible light source to controller, using a MIS-Bus cable
“No USB storage available”	No USB storage device connected to controller	Connect USB storage device to front side of the device
“USB storage device in incorrect format”	USB storage device with incompatible file system	Format USB storage device using the FAT32 file system
“USB device connected to wrong port”	USB device inserted on the rear of the device	Insert USB device on the front side of the controller
“USB device is read-only”	USB device defective	Replace USB device
“Unknown USB device error”		
“USB device removed while saving”	USB device removed while saving	Remove USB device only after saving has concluded
“White balance failed”	White balance not performed as described in the instructions for use	Repeat white balance (section 5.6.3)

Table 7-2: Troubleshooting table.

8 Product Data

8.1 Technical Data

Dimensions (W x H x D)	295 x 100 x 355 mm
Weight	approximately 5.7 kg
Current consumption	0.8 - 0.4A
Supply voltage	100-240V~, 50/60Hz
Protection class as per IEC 60601-1	Protection class I
Device fuse	T2.0AH 250V
Image refresh rate	50/60Hz (switchable)
Video outputs	2x DVI 1080p 2x 3G-SDI 1080p
Interfaces	1x MIS-Bus 2x Remote IN, 2x Remote OUT 2x USB
Degree of protection	IP 21
Classification as per IEC/CISPR 11	Group 1 Class B

8.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	50 kPa to 106 kPa

Operating conditions

Temperature	+10 °C to +40 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

CAUTION

Non-observance of the ambient conditions

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

> Observe operating conditions as well as transport and storage conditions

8.3 Spare Parts and Accessories

Use original spare parts and accessories only.

Designation	Item number
Micro-fuse T2.0AH 250V, 5x20 mm	A075305
MIS-Bus cable, 2.25 m	A057635
MIS-Bus cable, 0.75 m	A059584
DVI-cable	A090275
3.5 mm jack plug cable, 1.2 m	A056998
BNC cable, 2 m	BNC-BNC.2
Power cord (country-specific)	on request
USB keyboard, German layout	PS00510
USB keyboard, US layout	PS00511
Foot switch	6226-ACAA-BZZZ-002
Battery type VARTA CR2032	A052755

Table 8-1: Spare parts and accessories.

8.4 Open-Source Licenses

Information on open-source licenses

The device uses open-source licenses. For more information on the terms governing their use, please refer to the supplied brochure specifying the open-source licenses used.

9 Disposal

WARNING

Contaminated device

Risk of infection

> The device must be reprocessed prior to disposal (*chapter 6*)

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

The product is packed in a polyethylene (PE) bag and a PE foam inlay (0.50 kg), as well as a corrugated cardboard box (1.18 kg). Other packing materials used are a cardboard insert (0.22 kg) and a membrane pad made of PE film and corrugated cardboard (0.30 kg).

Any product carrying this symbol must be disposed of separately through dedicated electrical and electronic devices recycling. Within the EU, such disposal is taken care of by the manufacturer free of charge.

SCHÖLLY FIBEROPTIC GMBH
Robert-Bosch-Str. 1-3
79211 Denzlingen
Germany