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LED Lichtquelle

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Instructions For Use

LED Light Source

Light Source for Medical Endoscopy

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	Content	Page
0	Important Information about this Document.....	39
0.1	Scope of Validity, Identification, Purpose.....	39
0.2	Target Group.....	39
0.3	Using and Storing this Document.....	39
0.4	Supplementary Documents.....	39
1	General Information about the Device.....	40
1.1	Scope of Delivery.....	40
1.2	Product Description.....	40
1.2.1	Performance Characteristics and Function.....	40
1.2.2	Visual Overview.....	41
1.2.3	Compatible Light Guides.....	42
1.3	Usage.....	43
1.3.1	Intended Use.....	43
1.3.2	Indications.....	43
1.3.3	Contraindications.....	43
1.4	Conformity.....	43
1.4.1	Standards and Directives.....	43
1.4.2	Medical Device Classification.....	43
1.5	Marking.....	44
1.5.1	Pictograms and Information on the Device and Packaging.....	44
1.5.2	Pictograms in this Document.....	45
1.6	Service Department Contact Details.....	45
1.7	Obligation to Report Serious Incidents.....	45
2	General Safety Information.....	46
2.1	Safety Messages in this Document.....	46
2.1.1	Safety Messages at the Start of a Chapter.....	46
2.1.2	Safety Messages in the Body of the Text.....	47
2.2	Product Safety.....	47
2.2.1	Basic Safety Notices.....	47
2.2.2	Staff Qualifications.....	48
2.2.3	Electromagnetic Compatibility.....	49
2.2.4	Combination with Medical Electrical Equipment.....	50
2.2.5	Display of Faults.....	50

3	Installation and Initial Operation	51
3.1	Safety Notices	51
3.2	Staff Qualifications	52
3.3	Installation	53
3.3.1	Setting up the Device	53
3.3.2	Connecting the Device to the Mains	54
3.4	Initial Operation	54
3.4.1	Performing a Function Check	54
4	Operation	55
4.1	Safety Notices	55
4.2	Staff Qualifications	56
4.3	Technical Inspection Prior to Use	57
4.3.1	Performing a Visual Inspection	57
4.3.2	Performing a Function Check	57
4.4	Switching On and Off, Separation from the Power Supply	58
4.5	Connecting a Light Guide	58
4.6	Decoupling the Light Guide	59
4.7	Setting the Amount of Light	59
4.7.1	Manual Adjustment	59
4.7.2	Automatic Adjustment	59
4.8	Storage and Transport	60
5	Reprocessing	61
5.1	Safety Notices	61
5.2	Performing Cleaning and Disinfection	61
6	Maintenance and Repair	62
6.1	Safety Notices	62
6.2	Staff Qualifications	62
6.3	Safety Inspection	62
6.4	Changing the Fuse	63
6.5	Repair	64
6.6	Troubleshooting Table	65
7	Product Data	66
7.1	Technical Data	66
7.2	Ambient Conditions	66
7.3	Spare Parts and Accessories	67
8	Disposal	70

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0 Important Information about this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following product:

Item name: **LED Light Source**

Item number: **05-0760led**

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and medical engineers who are entrusted with the operation, cleaning, and maintenance of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

0.4 Supplementary Documents

The electromagnetic compatibility (EMC) requirements are described in the included separate brochure

Medical electrical equipment

pursuant to IEC 60601-1-2

The rules and information contained therein must be observed during installation and operation of this device.

1 General Information about the Device

1.1 Scope of Delivery

The scope of delivery for the device includes:

- 1x light source (item no. 05-0760led)
- 1x power cord (item no. available on request)
- 1x instructions for use
- 1x EMC brochure (TPI009-121-20)

Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Product Description

1.2.1 Performance Characteristics and Function

Medical endoscopy	The present product is a light source for use in medical endoscopy.
White light with a CRI >90	The light source provides white light with a color rendering index of CRI R _a >90. The high color rendering index allows for a particularly natural color representation of the endoscopic image.
Auto light control	If a compatible camera system with automatic light control is connected, the required light intensity can be controlled fully automatic.
Multi-light-guide adapter	The multi-light-guide adapter allows the light source to be connected with light guides manufactured by Storz, Olympus and Richard Wolf.
Remote controlled	With a compatible camera system, the light source functions can be remotely controlled, for instance using the camera head buttons.
Display	Operation instructions and other messages are shown on the device display.
Safe	After a power failure, the light source will restart in standby mode. This prevents any unintentional light emission.
Glare protection	Without a connected light guide, no light can be emitted.
Silent	The light source is designed to be inaudible, even from short distances. This contributes to improving OR working conditions and reducing stress.
Defibrillation-proof type CF	The light source meets the high defibrillation-proof type CF standard.
EMC Class B	The device is suitable for operation in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.
LED	The light is generated using the latest LED technology. The long life of the LED means that the device is maintenance free. Be sure to observe the contents of <i>Maintenance and Repair</i> .

1.2.2 Visual Overview

Front view


Figure 1-1: Front view of the light source.

- ① ON / standby button
- ② Multi-light-guide adapter
- ③ Display
- ④ Light intensity adjustment


Rear view

Figure 1-2: Rear view of the light source.

- ❶ Fuse holder
- ❷ Connection for power cord
- ❸ Main power switch
- ❹ Potential equalization terminal
- ❺ MIS-Bus ports
- ❻ Port for service (covered)

1.2.3 Compatible Light Guides

The product is compatible with cold light guides and high-power light guides supplied by Karl Storz, Richard Wolf, Olympus and our own company:

- Fiber bundles of 3.5 to 4.8 mm Ø
- Length up to 300 cm

1.3 Usage

1.3.1 Intended Use

This light source is intended for use in endoscopy systems for human medicine. In combination with a light guide, endoscope, and camera, its purpose is to illuminate the inside of the human body.

1.3.2 Indications

There are no medical indications associated with the use of the light source if separated from the endoscopic system.

The indications for an endoscopic application depend on the medical condition of the patient and the individual risk/benefit analysis by the surgeon.

1.3.3 Contraindications

The use of the endoscopic light source and its accessories is contraindicated in all cases where endoscopic procedures are contraindicated for any reason.

As is appropriate with any surgical procedure, consideration must be given to patient size and workspace volume when using the endoscopic light source.

Depending on the patient's medical condition, there may exist contraindications that rely on the general condition of the patient or the specific disease pattern.

The decision to perform an endoscopic procedure rests with the responsible surgeon and should be made on the basis of an individual risk-benefit analysis.

1.4 Conformity

1.4.1 Standards and Directives

The device meets the requirements of the following applicable standards:

- **IEC 60601-1** Medical electrical equipment - Part 1:
General requirements for basic safety and essential performance
- **IEC 60601-1-2** Medical electrical equipment - Part 1-2:
Electromagnetic compatibility

1.4.2 Medical Device Classification

The product is a Class I medical device within the meaning of applicable **EU guidelines**.


1.5 Marking

1.5.1 Pictograms and Information on the Device and Packaging

	Adhere to the instructions for use
	CE marking
	Item number
	Serial number
	Medical Device
	Manufacturer
	Date of manufacture
	Caution (IEC 60601-1 3rd and 4th edition) / Take note of accompanying documents (IEC 60601-1 3rd and 4th edition)
	MR unsafe
	Type BF applied part according to IEC 60601-1
	Type CF defibrillation-proof applied part according to IEC 60601-1
	Potential equalization
	Release
	ON / standby
	Permissible storage and transport temperature
	Permissible relative air humidity during storage and transport
	Permissible atmospheric pressure during storage and transport
R_x only	Caution: Federal (US) law restricts this device to sale by or on order of a physician
	Separate collection for WEEE (waste of electrical and electronic equipment)
	Fuse
	Alternating current

Hot temperature

Possibly hazardous optical radiation

1.5.2 Pictograms in this Document

General warning sign

Electricity hazard warning sign

Biohazard warning sign, risk of infection

Hot surfaces warning sign

1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

1.7 Obligation to Report Serious Incidents

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 General Safety Information
2.1 Safety Messages in this Document
2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

<p>⚠ WARNING</p> <p>Nature and origin of the injury risk</p> <p>Potential consequences of non-observance</p> <p>> Preventive action</p>
--

The following kind of message refers to a risk that could lead to minor or moderate injury:

<p>⚠ CAUTION</p> <p>Nature and origin of the injury risk</p> <p>Potential consequences of non-observance</p> <p>> Preventive action</p>
--

The following kind of message refers to a risk that could lead to material damage:

<p>NOTICE</p> <p>Nature and origin of the risk of material damage</p> <p>Potential consequences of non-observance</p> <p>> Preventive action</p>

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

Keep the original packaging. Transport and store the product in its original packaging and use it to return goods if service support is required.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

If you detect any faults or malfunctions, inform us immediately.

WARNING! Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

WARNING! Operating the device in an open condition. Risk of electric shock. Do not open the device. Operate only when closed.

WARNING! Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

WARNING! Risks from the arrangement, setup, combination, or properties of connected or surrounding devices or equipment. Follow the instructions for use of the respective devices. Perform a risk assessment.

WARNING! High-intensity light source. Risk of injury to eyes. Do not look directly into the open end of the light guide.

WARNING! Risk of suffocation. Keep packaging material out of reach of children.

WARNING! Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

WARNING! Using the device during operation of a defibrillator. Risk to persons. Before discharging a defibrillator, remove device from the surgical field.

CAUTION! Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

WARNING! Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

2.2.2 Staff Qualifications

Specific qualifications are required for installation, operation, and maintenance of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

2.2.3 Electromagnetic Compatibility

Medical electrical equipment is subject to stringent electromagnetic compatibility (EMC) requirements.

Despite the device's high interference immunity and low emitted interference, compliance with EMC-related requirements is necessary with regard to installation, installation location, and ambient conditions.

It is therefore important that you observe the related safety messages in this document and the EMC brochure included with the device.

According to IEC/CISPR 11, medical electrical equipment is divided into the subsequent groups and classes.

Refer to Technical Data (section 7.1) to find out which apply to the present device.

Group 1: The device uses or generates RF energy only for its internal functioning.

Group 2: The device generates HF energy in the frequency range 9 kHz to 400 GHz in the form of electromagnetic radiation and through inductive or capacitive coupling for purposes of the treatment of material, its examination or analysis.

Class A: The device **is not suitable** for use in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

Note: The emissions characteristics of this device make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this device might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the device.

The use of the device is only intended for professional healthcare environment.

Class B: The device **is suitable** for operation in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

WARNING! Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

WARNING! Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the device specified by the manufacturer and covered by the scope of this document, including cables. Otherwise, degradation of the performance of the device could result.

2.2.4 Combination with Medical Electrical Equipment

The device can be combined with components from other manufacturers provided that all components are compliant with the medical electrical equipment safety requirements according to IEC 60601-1.


It is the operator's responsibility to check and make sure that the system is and remains fully operational.

When using devices from different manufacturers and when operating an endoscope and/or endoscopic accessories with electromedical devices, it must be ensured that the applied part is properly insulated (in accordance with the requirements for type BF/CF applied parts defined in IEC 60601-1).

Contact our Service Department if you have any questions regarding compatibility.

2.2.5 Display of Faults

No light guide or light guide not attached correctly

If the light guide is not attached or is not attached correctly, this message appears on the light source display.

If this occurs, you should attach the light guide, attach it correctly, or switch the device off.

Warning with error code

In case of an operating error, this message (exemplary illustration) appears on the light source display.

Possible causes

LED light intensity cannot be controlled or LED is defective

Failure of housing fan

Housing temperature exceeds 60 °C. Light output is decreased. The device will not switch off, but is not guaranteed to remain fully functional

The device features two internal fans for cooling the LED. The first fan is defective

The device features two internal fans for cooling the LED. The second fan is defective

Fixes

Replace LED light source as quickly as possible by a second unit

Replace LED light source as quickly as possible by a second unit

Replace LED light source as quickly as possible by a second unit

Replace LED light source as quickly as possible by a second unit

Replace LED light source as quickly as possible by a second unit

3 Installation and Initial Operation

3.1 Safety Notices

WARNING

Electrical connections installed improperly

Risk of fire, short circuit, or electric shock

- > Make sure that the electrical connections are installed in accordance with the relevant national technical regulations

WARNING

Medical electrical system installation

Risk of fire, short circuit, or electric shock

- > Installation must only be performed by qualified personnel
- > If a number of electrical devices are being used in combination, observe IEC 60601-1
- > Non-medical devices that meet the applicable IEC safety standards can be connected only via a medical isolating transformer
- > Do not connect any additional non-medical devices to a medical electrical system
- > Signal lines from devices with a functional ground (FE) connection, which are connected to different branches of the mains power supply, must be galvanically isolated at both ends
- > Connect the devices only to a power supply with a protective ground (PE)
- > After installing a medical electrical system, it must be inspected according to IEC 62353

WARNING

Electromagnetic interactions

Possible impaired function and image quality

- > Observe the EMC brochure supplied
- > Make sure that neighboring devices are not causing reciprocal electromagnetic interference
- > Perform a post-installation function check

**⚠ WARNING****Use of multiple sockets**

Risk of fire, short circuit, electrical shock, reduced level of safety

- > Whenever possible, avoid the use of multiple sockets
- > If required, use medically approved multiple sockets
- > Never connect multiple sockets in series
- > Do not cover multiple sockets (leads to heat accumulation)
- > Do not place multiple sockets on the floor
- > Use traction relief
- > Devices should be connected to a shared multiple socket only if they are all intended to be used as part of the medical electrical system

3.2 Staff Qualifications

Personnel responsible for setting up or installing medical electrical equipment must possess the necessary qualifications and be familiar with the safety regulations in place at the installation location and the applicable official regulations.

3.3 Installation

3.3.1 Setting up the Device

Observe the safety notices provided at the start of this chapter and refer to the enclosed EMC brochure.

WARNING! Risk from installation in potentially explosive locations. Increased risk of fire and explosion in oxygen-enriched atmospheres. Device must be installed outside of potentially explosive locations and not in the vicinity of flammable substances.

WARNING! Risk resulting from electromagnetism. Possible malfunctions and image interference. Use only the cables supplied. Ensure that all devices operated in the vicinity meet the EMC requirements. Perform a post-installation function check.

WARNING! Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other devices should be observed to verify that they are operating normally.

1. Install the device.

Maintain minimum distances to other devices:

- at least 5 cm to the side,
- at least 6 cm from the back.

When installing the device, ensure that:

- it is positioned on an level, non-slip surface,
- it is located outside of the sterile area on a sufficiently stable support,
- it is protected from water drips and splashes,
- it will not be subject to vibration during operation,
- the ventilation slots are not obstructed,
- the main switch on the rear of the device is accessible.

2. Connect the **potential equalization terminal** on the rear of the device to the potential equalization terminal at the installation site or on the equipment cart.

Use a potential equalization conductor in accordance with DIN 42801 with at least a 4 mm² cross-section.

3. If you would like to control the light source using the buttons on the camera head or the camera controller, connect the camera to the light source at one of the ports labeled **MIS-BUS**.

Observe the instructions for use supplied with the camera.

3.3.2 Connecting the Device to the Mains

Do not connect the device to the mains until you have completed all the previous activities. Observe the safety notices provided at the start of this chapter.

1. Check that the supply voltage at the site of operation corresponds to that specified on the device nameplate.
2. Insert the power cord into the mains connection socket on the rear of the device.

WARNING! Sparks in a potentially explosive location. Risk of explosion. Ensure that the device is not situated in a potentially explosive location when connecting to the power supply.

3. Connect the power cord to the power supply.
4. When routing the cable, make sure it does not present a trip or entanglement hazard.

3.4 Initial Operation

3.4.1 Performing a Function Check

When operating the device for the first time, check the installation and make sure the device is in full working order.

Precondition: Installation must have been completed.

1. Switch the main switch on the rear of the device to ON.
 - > *The device switches to standby mode and the ON / standby button on the front is illuminated blue.*
2. Switch on the device.
 - > *After approximately 10 seconds, the device is ready for operation.*
3. Connect a light guide as described in [section 4.5](#).

WARNING! High-intensity light. Risk of injury to eyes. Do not look directly into the open end of the light guide.

4. Check that you are able to adjust the light intensity as described in [section 4.7](#).
5. Check that the LED can go into standby mode and be reactivated as described in [section 4.4](#).

4 Operation

4.1 Safety Notices

WARNING

Modifications to the installation

Risk of fire, short circuit, and electric shock

- > Follow the safety notices in the Installation and Initial Operation chapter
- > Do not make unauthorized changes to the electrical installation
- > If the installation contains a multiple socket, do not connect additional devices to this without authorization
- > Never connect multiple sockets in series

WARNING

Current leakage when touching the patient

Risk of electric shock

- > Never touch the patient and the device housing or accessible contacts of electric devices at the same time
- > Never touch the patient and non-medical electric devices at the same time

WARNING

Positioning the equipment cart inside the sterile area

Risk of infection

- > Do not position the equipment cart inside the sterile area
- > Measure the cable lengths so there is enough room for movement without compromising sterility
- > The light guide must be anchored securely within the direct OR environment to avoid any risk of accidental slipping

WARNING

Interactions between devices in simultaneous use (e.g. lasers, electro-surgery)

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the required type BF, CF or CF defibrillation-proof classification requirements according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active HF electrodes
- > Do not activate HF electrodes in the vicinity of flammable gases or liquids
- > Before using HF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

CAUTION

High temperatures in combination with light guides and endoscopes

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest possible light output setting that will allow you to illuminate the target area
- > Do not touch ends of the light guide or allow to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Do not touch light guide connection and distal end of the endoscope or allow to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Always use sheath with endoscope

4.2 Staff Qualifications

This device may be operated only by physicians and medical assistants with the requisite training, advanced training, knowledge, and practical experience in the endoscopic discipline concerned as defined by the provisions in place at the site of operation.

4.3 Technical Inspection Prior to Use

4.3.1 Performing a Visual Inspection

*A visual inspection must be performed **prior to each medical procedure.***

1. Inspect the housing of the light source for external damage.
Do not use the light source if there is any external damage to the housing.
2. Check that all electric cables are undamaged and positioned safely.
Do not use the device if the electric cables are damaged or present a trip or entanglement hazard.
3. Check that the light guide you intend to use is undamaged.
Do not use the light guide if it is broken or twisted.

4.3.2 Performing a Function Check

*A function check must be performed **prior to each medical procedure.***

1. Connect the light guide.
2. The light guide should be directed away from the eyes to avoid dazzling anybody.
3. Activate the light source.
Check that you are able to adjust the light intensity and that the system is functioning correctly.
Do not use the light source if you are unable to adjust the light intensity.

4.4 Switching On and Off, Separation from the Power Supply

Switching on the device

1. Switch on the device at the front.

Device does not switch on?

The main switch on the rear is switched off.

Switch ON at the main switch.

- > After approximately 10 seconds, the device is ready for operation.
- > The ON / standby button on the front is illuminated white.

Switching off the device

1. Switch off the device at the front.
 - > The device switches to standby mode.
 - > The ON / standby button on the front is illuminated blue.
2. To switch off the device completely, switch off at the main switch on the rear.

Separating the device from power

1. To separate the device completely from its power supply, disconnect the power cord. The appliance coupler is considered as disconnecting device.
Make sure that the rear of the device is accessible at any time.

4.5 Connecting a Light Guide

The multi-light-guide connector enables the connection of light guides supplied by different manufacturers.

Original light guides from Karl Storz, Richard Wolf and Olympus can be attached directly without adapters. Light guides of other manufacturers than that must be used with the corresponding light guide adapters.

WARNING! High-intensity light. Risk of injury to eyes. Do not look directly into the open end of the light guide.

1. Slide the light guide into the multi-light-guide adapter until it clicks into place, as far as it will go.
 - > The term **STANDBY** and the corresponding symbol appear on the display.

4.6 Decoupling the Light Guide


1. Rotate the multi-light-guide adapter clockwise as far as it will go.
 - > *The light guide detaches.*
 - > *The LED switches off, the device remains switched on.*
2. Pull out the light guide.
3. Detach possible adapters used, and arrange for reprocessing the light guide. Observe the instructions for use supplied with the light guide when doing so.

4.7 Setting the Amount of Light

4.7.1 Manual Adjustment

The light intensity can be set in 5 % increments.

1. Press the arrow buttons to set the light intensity.
 - Tap an arrow button briefly to set the light intensity in 5 % increments.
 - > *The regulated value of light intensity is shown on the display in “%”.*
 - Press and hold down an arrow button to continually increase or decrease the light intensity.

4.7.2 Automatic Adjustment

Compatible cameras allow for the use of the auto light control feature. To this purpose, it is necessary to activate auto light control. This can be done using the camera.

Activating automatic light control


1. Using the MIS-Bus cable, connect the light source to a compatible camera. Observe the instructions for use supplied with the camera.
2. Activate auto light control via the camera menu.

- > *When auto light control is activated, this message appears on the light source display.*

Deactivating automatic light control

1. Deactivate auto light control via the camera menu or pressing one of the light source arrow buttons.

- > *When auto light control is deactivated, the light intensity is shown on the display in “%”.*

4.8 Storage and Transport

To prevent the device from falling when being transported on the equipment cart, secure with a tension belt.


Figure 4-1: Device on equipment cart, secured with tension belt.

Store the device in a dry, well-ventilated, temperature-controlled room, where it is protected from dust.

Switch off the device prior to storage and remove the power cord and accessories.

When placing into storage or storing temporarily, make sure that the device cannot be damaged by its surroundings.

To avoid damage, keep out of direct sunlight and away from radioactivity and strong electromagnetic radiation.

The ambient conditions for transport, storage, and operation are described in [section 7.2](#).

5 Reprocessing

5.1 Safety Notices

WARNING

Risks associated with improper cleaning

Risk of fire, short circuit, and electric shock

- > Before cleaning, switch off the device at its rear main switch and disconnect from the mains power supply
- > Make sure no moisture penetrates the device
- > The device should not be reconnected to the mains power supply until it has completely dried

NOTICE

Improper cleaning

Risk of damage to the device

- > Do not use sponges or wipes that are dripping wet
- > Do not use abrasive cleaners or solvents

5.2 Performing Cleaning and Disinfection

1. Switch off the device at its rear main switch.

NOTICE! Pulling the power cable will damage it. Pull at the plug to disconnect device from power.

2. Separate the device from the power supply by pulling out the plug of the power cord from the rear of the device.
3. Use disinfectant wipes¹ to clean all contaminated areas of the device for at least **1 minute** and to completely remove visible soiling.
4. Wipe the product until no more soiling can be seen. Use additional wipes, if required.
5. Use a fresh wipe to thoroughly wet all contaminated surfaces, ensuring that the device stays visibly wet for **3 minutes**. Use additional wipes, if required.
6. Continuously wipe the device during the **3 minutes**. Check that all seams, creases, crevices, recessed areas, and mated surfaces come into contact with the disinfectant.
7. Wait until all cleaned parts have thoroughly dried.
8. Reconnect the device to the power supply.

¹ Formulation: quaternary ammonium, bactericidal, virucidal, suitable for use on equipment.
For cleaning and disinfection validation Sani-Cloth Plus disinfectant wipes manufactured by PDI® were used.

6 Maintenance and Repair

6.1 Safety Notices

WARNING

Improper maintenance and repair

Risk of fire, short circuit, and electric shock

- > All maintenance and repair work must only be performed by qualified personnel
- > Do not make any modifications to the device

6.2 Staff Qualifications

Personnel responsible for maintaining the device or carrying out safety inspections must possess the necessary qualifications and be familiar with the local safety regulations in place at the installation location and the applicable official regulations.

Repairs must be performed by our service personnel only.

6.3 Safety Inspection

We recommend to arrange for safety inspections to be performed every 12 months in accordance with Section 5 of IEC 62353.

The inspection comprises the following steps:

- visual inspection
- measuring the protective resistance (only for protection class I devices)
- measuring the leakage currents
- measuring the insulation resistance

The results of the inspection should be documented in an inspection report in accordance with Section 6 and Annex G of IEC 62353.

We recommend that all components of a medical electrical system be subjected to documented technical safety inspections at regular intervals.

6.4 Changing the Fuse

1. Switch OFF the device at the main switch on the rear and pull out the power cord from the rear of the device to disconnect the power.

NOTICE! Pulling the power cable will damage it. Disconnect the power at the plug.

2. Release the fuse holder with a screwdriver.


Figure 6-1: Releasing the fuse holder with a screwdriver.

3. Pull out the fuse holder carefully, without tilting.

WARNING! Unsuitable electrical fuse. Risk of fire and short circuit. Use the specified fuse only and never bypass the fuse.

4. Replace the faulty fuse.
5. Carefully slide the fuse holder back in until it clicks into place.

6.5 Repair

Should you need to arrange a repair for the device, please contact one of our subsidiaries. Contact details can be found on the back of these instructions for use.

When you send in equipment, please enclose as accurate a **fault description** as possible, and record the item number and serial number of the product on the delivery note. You will find these details on the nameplate.

WARNING

Contaminated device

Risk of infection

- > Reprocess the product prior to shipping (*chapter 5*)
- > Clearly identify contaminated products as such

Send only thoroughly cleaned goods in for repair.

For shipping, please use the original packaging whenever possible. Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

6.6 Troubleshooting Table

Issue	Possible causes	Fixes
No operation	No power supply	Check the power connection and connect the power if necessary Check the fuse and replace if necessary
	Malfunction in power supply unit	Send the light source in for repair
No light from the light source	Device is overheated	Allow device to cool down
The light guide cannot be coupled	Adapter on light guide connection missing	Screw on adapter
Not enough light	Light guide not connected correctly	Check that the light guide is seated properly; correct if necessary
	Light guide defective	Replace the light guide
	Adapter not screwed on correctly	Fully screw on adapter
	Light guide not completely inserted	Use a compatible light guide
Fans do not work	Malfunction in fans	Send the light source in for repair
Light source cannot be remotely controlled	MIS-Bus cable not properly plugged in or defective	Check connection or replace cable
Auto light control not possible	Camera does not support auto light control	Use a compatible camera model
Light does not switch off, despite light guide removed	Light guide was not locked	Rotate multi-light-guide adapter
	Incompatible light guide	Use a compatible light guide

Table 6-1: Troubleshooting table.

7 Product Data

7.1 Technical Data

Dimensions (W x H x D)	295 x 130 x 355 mm
Weight	8 kg
Current consumption	1.6 - 0.6A
Supply voltage	100-240V~, 50/60Hz
Protection class acc. to IEC 60601-1	Protection class I
Device fuse	T1.6AH 250V
Interfaces	2x MIS-bus
Degree of protection	IP 21
Applied part as per IEC 60601-1	Type CF defibrillation-proof
Classification as per IEC/CISPR 11	Group 1 Class B

7.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	50 kPa to 106 kPa

Operating conditions

Temperature	+10 °C to +40 °C
Relative air humidity	10 % to 90 %
Atmospheric pressure	70 kPa to 106 kPa

CAUTION

Non-observance of the ambient conditions

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Observe operating conditions as well as transport and storage conditions

7.3 Spare Parts and Accessories

Use original spare parts and accessories only.

Image	Designation	Item number
	Micro-fuse T1.6AH, 5x20 mm, UL CSA	A085442
	MIS-Bus cable, 2.25 m	A057635
	MIS-Bus cable, 0.75 m	A059584
	Power cord (country-specific)	on request

Table 7-1: Electrical spare parts and accessories.

Image	Designation	Item number
	High-power light guide, 3.5 x 1,800 mm (standard light guide)	05.0084I
	High-power light guide, 3.5 x 2,300 mm (standard light guide)	05.0088I
	High-power light guide, 3.5 x 3,000 mm (standard light guide)	05.0085I
	High-power light guide, 4.8 x 1,800 mm (standard light guide)	05.0086I
	High-power light guide, 4.8 x 2,300 mm (standard light guide)	05.0090I
	High-power light guide, 4.8 x 3,000 mm (standard light guide)	05.0087I
	High-power light guide, 3.5 x 2,300 mm (high-temperature light guide)	05.0088I.ht
	High-power light guide, 3.5 x 3,000 mm (high-temperature light guide)	05.0085I.ht
	High-power light guide, 4.8 x 2,300 mm (high-temperature light guide)	05.0090I.ht

Table 7-2: Light guides.


Image	Designation	Item number
To connect to the light guide (light-source-side):		
	Light guide adapter, light-source-side, Storz systems	05.0100z
	Light guide adapter, light-source-side, Wolf systems	05.0102b
	Light guide adapter, light-source-side, Olympus systems	05.0101o
To connect to the light guide (endoscope-side):		
	Light guide adapter, endoscope-side, Storz systems	05.0108z
	Light guide adapter, endoscope-side, Wolf systems	05.0110b
	Light guide adapter, endoscope-side, ACMI systems	05.0112c
To connect to the endoscope:		
	Endoscope light guide adapter, Storz systems	05.0114z
	Endoscope light guide adapter, Wolf systems	05.0116b

Table 7-3: Light guide adapters.

8 Disposal

WARNING

Contaminated device

Risk of infection

> Reprocess the product prior to disposal (*chapter 5*)

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

The product is packed in a polyethylene (PE) bag and a PE foam inlay (0.50 kg), as well as a corrugated cardboard box (1.18 kg). Other packing materials used are a cardboard insert (0.22 kg) and a membrane pad made of PE film and corrugated cardboard (0.30 kg).

Any product carrying this symbol must be disposed of separately through dedicated electrical and electronic devices recycling. Within the EU, such disposal is taken care of by the manufacturer free of charge.

SCHÖLLY FIBEROPTIC GMBH
Robert-Bosch-Str. 1-3
79211 Denzlingen
Germany