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
## **NIR FI CCU 4K UHD**

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# **Instructions For Use**

## **NIR FI CCU 4K UHD**

4K Camera Control Unit for NIR Fluorescence Imaging  
(Not valid for USA)

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---

	<b>Content</b>	<b>Page</b>
<b>0</b>	<b>Important Information about this Document.....</b>	<b>74</b>
0.1	Scope of Validity, Identification, Purpose.....	74
0.2	Target Group.....	74
0.3	Using and Storing this Document.....	74
0.4	Supplementary Documents.....	74
<b>1</b>	<b>General Information about the Device.....</b>	<b>75</b>
1.1	Scope of Delivery.....	75
1.2	Product Description.....	76
1.2.1	Performance Characteristics and Function.....	76
1.2.2	Front View of the Camera Control Unit.....	77
1.2.3	Rear View of the Camera Control Unit.....	78
1.2.4	Menu - Functions, Access and Navigation.....	79
1.3	Usage.....	80
1.3.1	Intended Use.....	80
1.3.2	Indications.....	80
1.3.3	Contraindications.....	80
1.3.4	User Profile.....	80
1.3.5	Compatibility with HF Surgical Devices.....	80
1.3.6	Components Required for Operation.....	81
1.4	Conformity.....	81
1.4.1	Standards and Directives.....	81
1.4.2	Degree of Protection against Electric Shock.....	82
1.4.3	Medical Device Classification.....	82
1.5	Marking.....	83
1.5.1	Pictograms and Information on the Device and Packaging.....	83
1.5.2	Pictograms in this Document.....	84
1.6	Service Department Contact Details.....	84
1.7	Obligation to Report Serious Incidents.....	84

---

<b>2</b>	<b>General Safety Information.....</b>	<b>85</b>
2.1	Safety Messages in this Document.....	85
2.1.1	Safety Messages at the Start of a Chapter.....	85
2.1.2	Safety Messages in the Body of the Text.....	86
2.2	Product Safety.....	86
2.2.1	Basic Safety Notices.....	86
2.2.2	Staff Qualifications.....	87
2.2.3	Electromagnetic Compatibility.....	87
2.2.4	Combination with Medical Electrical Equipment.....	89
2.2.5	Display of Faults.....	89
<b>3</b>	<b>Installation and Initial Operation.....</b>	<b>90</b>
3.1	Safety Notices.....	90
3.2	Staff Qualifications.....	91
3.3	Installation.....	91
3.3.1	Setting up the Device.....	91
3.3.2	Connecting Monitors.....	92
3.3.3	Connecting a Keyboard.....	93
3.3.4	Connecting Foot Switches.....	93
3.3.5	Connecting a Recording Device.....	94
3.3.6	Connecting and controlling a Light Source.....	94
3.3.7	Connecting the Device to the Mains.....	95
3.3.8	Performing a Function Check.....	95

---

<b>4</b>	<b>Configuration.....</b>	<b>96</b>
4.1	Basic Settings.....	97
4.1.1	Setting System Language.....	97
4.1.2	Changing the On-Screen Display Position.....	97
4.1.3	Setting Date and Time.....	98
4.1.4	Defining File Format for Photos (optional).....	98
4.1.5	Defining Image Refresh Rate.....	99
4.1.6	Defining On-Screen Displays.....	99
4.1.7	Switching Signal Sounds for taking Photos on or off (optional).....	100
4.1.8	Controlling the Amount of Light.....	100
4.1.9	Setting Shutter Speed.....	101
4.1.10	Configuring Foot Switches.....	102
4.1.11	Displaying System Information.....	103
4.2	Image Settings.....	104
4.2.1	General Information.....	104
4.2.2	Advanced Views.....	104
4.2.3	Customizing the Image.....	105
4.2.4	Changing the Image Settings of Specialties and Color Presets.....	107
4.2.5	Making and saving User-Specific Image Settings as a Profile.....	108
4.3	Reset to Default Settings.....	109
4.3.1	Resetting the Device to Default Settings.....	109
4.3.2	Resetting the Image Settings of a Specialty.....	109

---

<b>5</b>	<b>Operation</b>	<b>110</b>
5.1	Safety Notices	110
5.2	Staff Qualifications	112
5.3	General Notes on Operation	112
5.4	Technical Inspection Prior to Use	112
5.4.1	Performing a Visual Inspection	112
5.4.2	Performing a Function Check	113
5.5	Switching On and Off, Separation from the Power Supply	113
5.6	Preparation for Use	113
5.6.1	Connecting up a Camera Head	113
5.6.2	Selecting a Specialty or Color Preset	114
5.6.3	Setting the White Balance	115
5.6.4	Entering Patient Data	115
5.7	NIR Fluorescence Imaging Operation	116
5.7.1	Toggling between NIR FI and White Light Imaging Mode	116
5.7.2	Adjusting the Image in NIR FI Mode	117
5.8	Creating Recordings	118
5.8.1	Recording a Photo	118
5.8.2	Recording a Video	118
5.8.3	Recording a Photo or Video using an External Device	119
5.9	Storage and Transport	119
<b>6</b>	<b>Reprocessing</b>	<b>120</b>
6.1	Safety Instructions	120
6.2	Performing Cleaning and Disinfection	121
<b>7</b>	<b>Maintenance and Repair</b>	<b>122</b>
7.1	Safety Notices	122
7.2	Staff Qualifications	122
7.3	Safety Inspection	122
7.4	Changing the Battery	123
7.5	Changing the Fuse	124
7.6	Repair	125
7.7	Troubleshooting	126
<b>8</b>	<b>Product Data</b>	<b>129</b>
8.1	Technical Data	129
8.2	Ambient Conditions	129
8.3	Spare Parts and Accessories	130
8.4	Information on Open-Source Licenses	130
<b>9</b>	<b>Disposal</b>	<b>131</b>

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## **0 Important Information about this Document**

### **0.1 Scope of Validity, Identification, Purpose**

These instructions for use apply to the following product:

Item name: **NIR FI CCU 4K UHD**

Item number: **95-3985**

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

### **0.2 Target Group**

These instructions for use are intended for physicians, medical assistants, medical engineers and employees of central sterile supply who are entrusted with the installation, operation, maintenance, and reprocessing of the device.

### **0.3 Using and Storing this Document**

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

### **0.4 Supplementary Documents**

EMC The electromagnetic compatibility (EMC) requirements are described in the included separate brochure

**Medical electrical equipment**

pursuant to IEC 60601-1-2

The rules and information contained therein must be observed during installation and operation of this device.

Instructions for use of camera head and light source

To safely use this device, also observe the instructions for use of the camera head and light source being used.

# 1 General Information about the Device

## 1.1 Scope of Delivery

The scope of delivery for the device includes:

- one camera control unit
- four BNC cables (2 m)
- one DVI cable (2 m)
- two 3.5 mm jack plug cables (1.2 m)
- one cable for the MIS-Bus port (0.75 m)
- one power cord (item no. available on request)
- one set of instructions for use
- one brochure specifying the open-source licenses used (TPO003-121-01)
- one EMC brochure (TPI014-121-20)

Check scope of  
delivery!

Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

## 1.2 Product Description

### 1.2.1 Performance Characteristics and Function

CMOS camera	The combination of the camera control unit and the camera head result in a camera based on CMOS technology that is used to provide endoscopic real-time video for display on a monitor.
4K/UHD	When used with a compatible monitor, the camera offers a 16:9 4K/UHD image resolution using progressive scanning (3840 x 2160p).
NIR fluorescence imaging (NIR FI)	The camera control unit can be used to visualize near-infrared light for fluorescence imaging applications. For this purpose, the system is designed to work with a light source that provides NIR illumination in addition to visible white light, a camera head and endoscope sensitive in the visible and NIR spectrum, and indocyanine green (ICG) which promotes the fluorescing appearance of the camera image.
Auto light control	If a compatible light source with automatic light control is connected, the required light intensity can be controlled fully automatic.
Extensive configuration options	<p>In addition to the general setting options, the camera offers various endoscopic application profiles. These profiles are preconfigured to offer an optimized image display; however, the profiles can also be modified and customized for your specific needs. Ten additional profiles can also be configured.</p> <p>The general setting options include language selection for the configuration menu and monitor displays, and the camera head button settings. The default language setting is English.</p>
Captures	<p>The present camera control unit allows creating and saving video and photo captures onto a USB storage device, which can be inserted on the front side of the camera control unit.</p> <p>The camera control unit can be used with various camera heads.</p>

### 1.2.2 Front View of the Camera Control Unit


Figure 1-1: Front view of the camera control unit.

#### Fig. 1-1 legend

- ① ON / standby button
- ② Connection for camera head
- ③ Photo capture button / upwards in menu
- ④ White balance button / downwards in menu
- ⑤ Open Menu button / Enter
- ⑥ USB storage device connection
- ⑦ Video capture start / stop button

#### Indicator light and control button color coding:

- ① **ON / standby button -**  
blue: device is in standby mode / white: device is ready for use
- ③ **Photo capture button / upwards in menu -**  
white: device is ready for use / green: laser is on (NIR FI mode)
- ④ **White balance button / downwards in menu -**  
white: device is ready for use / green: laser is on (NIR FI mode)
- ⑤ **Open Menu button / Enter -**  
white: device is ready for use / green: laser is on (NIR FI mode)
- ⑦ **Video capture start / stop button -**  
white: device is ready for use / green: laser is on (NIR FI mode)

### 1.2.3 Rear View of the Camera Control Unit


Figure 1-2: Rear view of the camera control unit.

#### Fig. 1-2 legend

- ① Appliance coupler
- ② Potential equalization terminal
- ③ Main power switch
- ④ MIS-Bus (MIS-BUS)
- ⑤ Recording device connection (Remote-Out)
- ⑥ Port for service (covered)
- ⑦ Foot switch connection (Remote-In)
- ⑧ 4x 3G SDI output for 4K/UHD (4K/UHD 3840 x 2160 SDI 1 - SDI 4)
- ⑨ Fuse holder
- ⑩ 2x DVI output (DVI 1 / DVI 2)
- ⑪ Battery compartment (Battery CR2032)
- ⑫ 2x 3G-SDI output (3G-SDI 1 / 2)
- ⑬ Keyboard connection (USB-keyboard)

### 1.2.4 Menu - Functions, Access and Navigation

Configuration	The menu is where you configure the device and define the basic settings and image settings.
Selecting the specialty	Furthermore, you select the specialty with the predefined, optimized profile settings and enter the patient data (optional).
Control functions	Perform image settings, white balance or captures before or during procedures, with the appropriate configuration, do so directly with the pre-defined camera head buttons.
Entry and navigation	Entry and navigation in the menu are done through the buttons on the front side of the device, with the keyboard or using the camera head buttons.

#### Navigation in the menu using the buttons on the front side of the device

	Entry into menu:	Press button briefly.
	Upwards in the menu:	Press button briefly.
	Downwards in the menu:	Press button briefly.
	Confirm selection in menu:	Press button briefly.
	Leave the menu:	Press and hold button.

#### Navigation in the menu using the keyboard

→	Entry into menu:	Press right arrow button.
↑	Upwards in the menu:	Press up arrow button.
↓	Downwards in the menu:	Press down arrow button.
→	Go one level deeper in menu:	Press right arrow button.
→	Confirm selection in menu:	Press right arrow button.
←	Go one level back in menu:	Press left arrow button.
→	Enter Setup menu:	Press and hold right arrow button.
→	Leave the menu:	Press and hold right arrow button.

## **1.3 Usage**

### **1.3.1 Intended Use**

The camera control unit is used in combination with a camera head or video endoscope to provide real-time video and documentation during diagnostic or therapeutic endoscopic procedures.

The camera control unit is also used to provide intraoperative fluorescence imaging.

### **1.3.2 Indications**

The camera control unit is indicated for use in surgical endoscopic procedures that physicians deem appropriate for the patient wherever a camera is required for intraoperative imaging.

The decision to perform an endoscopic procedure depends on the condition of the patient and remains with the responsible physician.

In order to assure a safe application, standards/guidelines released by official organizations such as scientific societies of endoscopic surgery have to be considered.

### **1.3.3 Contraindications**

The use is contraindicated in all cases where endoscopic procedures are contraindicated for any reason.

Medical contraindications that are listed for suitable fluorescence dyes that may be used apply.

For any interventions the responsible physician must decide applying professional standards and based on the general condition of the patient and an individual benefit-risk analysis whether the intended application of the device is admissible.

### **1.3.4 User Profile**

Intended to be used by trained healthcare professionals in a professional healthcare environment (hospital).

### **1.3.5 Compatibility with HF Surgical Devices**

The device is approved for use in conjunction with HF surgical devices. Applying HF surgery does not alter the conditions of the intended use.

### 1.3.6 Components Required for Operation

For NIR fluorescence imaging, the camera control unit is intended for use with the following components:

- NIR FI Camera Head Full HD Zoom (95-3908)
- NIR FI Light Source (05-0761nir)

For the complete system requirements for achieving the intended results when using the camera for NIR fluorescence imaging, observe the instructions for use supplied with the light source (05-0761nir).

Before using NIR FI mode, be familiar with the warnings and instructions in the instructions for use supplied with the compatible NIR FI light source.

For white light imaging, the camera control unit is intended for use with the following components:

#### Camera heads

- Camera Head 4K UHD Zoom (95-3905)
- Camera Head Full HD Zoom (95-3901)
- Camera Head Full HD 90° angled (95-3903)
- Camera Head Full HD autoclavable (95-3904)

#### Light source

- LED Light Source (05-0760led)

In order to guarantee full functionality, the joint use of the above-mentioned components is recommended plus:

- Endoscope
- Light guide

## 1.4 Conformity

### 1.4.1 Standards and Directives

The camera system consisting of camera control unit and camera head meets the requirements of the following applicable standards:

- **IEC 60601-1** Medical electrical equipment - Part 1:  
General requirements for basic safety and essential performance
- **IEC 60601-1-2** Medical electrical equipment - Part 1-2:  
Electromagnetic disturbances – Requirements and tests
- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:  
Particular requirements for basic safety and essential performance of endoscopic equipment

### 1.4.2 Degree of Protection against Electric Shock

Depending on the camera head connected, the product meets the type BF, CF or CF defibrillation-proof classification requirements according to IEC 60601-1 for protection against electric shock.

The camera control unit can be configured for cardiac applications in combination with applied parts intended for direct conductive contact with the heart only when TYPE CF applied parts (camera heads, video endoscopes and light guides) labelled as CF defibrillation-proof are used.

Classification types are clearly stated in the corresponding instructions for use and as as labeling on the device:


Products that meet the requirements classified as type BF feature the adjacent symbol.

Products that meet the type CF classification requirements feature the adjacent symbol, offer an even higher level of protection against electric shock and are approved for direct conductive contact with the heart.

Products that meet the type CF defibrillation-proof classification requirements feature the adjacent symbol, offer an even higher level of protection against electric shock and are defibrillation-proof.

Products that are defibrillation-proof in addition to that are marked with one of the adjacent symbols, depending on the degree of insulation.

### 1.4.3 Medical Device Classification

The product is a Class I medical device within the meaning of applicable **EU guidelines**.

## 1.5 Marking

### 1.5.1 Pictograms and Information on the Device and Packaging

Adhere to the instructions for use

CE marking

Item number

Serial number

Medical Device

Manufacturer

Date of manufacture

Caution (IEC 60601-1 3rd edition) /

Take note of accompanying documents (IEC 60601-1 2nd edition)

MR unsafe

Potential equalization

ON / standby

Permissible storage and transport temperature

Permissible relative air humidity during storage and transport

Permissible atmospheric pressure during storage and transport

**R<sub>x</sub>** only

Caution: Federal (US) law restricts this device to sale by or on order of a physician

Separate collection for WEEE (waste of electrical and electronic equipment)

Fuse

Alternating current

NIR FI Green imaging mode

NIR FI Contrast imaging mode

NIR FI Intensity imaging mode

Video capture

Laser on

### 1.5.2 Pictograms in this Document

General warning sign

Electricity hazard warning sign

Biohazard warning sign, risk of infection

Laser beam warning sign

## 1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

## 1.7 Obligation to Report Serious Incidents

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**2 General Safety Information**  
**2.1 Safety Messages in this Document**  
**2.1.1 Safety Messages at the Start of a Chapter**

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

**⚠ WARNING**

**Nature and origin of the injury risk**

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

**⚠ CAUTION**

**Nature and origin of the injury risk**

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

**NOTICE**

**Nature and origin of the risk of material damage**

Potential consequences of non-observance

> Preventive action

### 2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

**WARNING!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

**CAUTION!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

**NOTICE!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

## 2.2 Product Safety

### 2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Risks despite highest of quality standards

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

Keep the original packaging. Transport and store the product in its original packaging and use it to return goods if service support is required.

If you detect any faults or malfunctions, inform us immediately.

**WARNING!** Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

**WARNING!** Operating the device in an open condition. Risk of electric shock. Do not open the device. Operate only when closed.

**WARNING!** Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

**WARNING!** Risks from the arrangement, setup, combination, or properties of connected or surrounding devices or equipment. Follow the instructions for use of the respective devices. Perform a risk assessment.

**WARNING!** Risk of suffocation. Keep packaging material out of reach of children.

**WARNING!** Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

**CAUTION!** Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. Do not kink, crush or strongly bend camera cable. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

**WARNING!** Improper handling, maintenance, and use carries a risk to the patient, user and other persons, or can lead to premature wearing of the device.

### 2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for installation, operation, and maintenance of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

### 2.2.3 Electromagnetic Compatibility

Medical electrical equipment is subject to stringent electromagnetic compatibility (EMC) requirements.

Despite the device's high interference immunity and low emitted interference, compliance with EMC-related requirements is necessary with regard to installation, installation location, and ambient conditions.

It is therefore important that you observe the related safety messages in this document and the EMC brochure included with the device.

According to IEC/CISPR 11, medical electrical equipment is divided into the subsequent groups and classes.

Refer to Technical Data (section 8.1) to find out which apply to the present device.

**Group 1:** The device uses or generates RF energy only for its internal functioning.

**Group 2:** The device generates HF energy in the frequency range 9 kHz to 400 GHz in the form of electromagnetic radiation and through inductive or capacitive coupling for purposes of the treatment of material, its examination or analysis.

**Class A:** The device is **not suitable** for use in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

Note: The emissions characteristics of this device make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this device might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the device.

The device is only intended for use in professional healthcare environment (hospital).

**Class B:** The device **is suitable** for operation in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

**WARNING!** Use of accessories, transducers and cables not specified or supplied by the device manufacturer. Increased electromagnetic emissions/decreased electromagnetic immunity, improper operation. Use only accessories, transducers and cables specified or supplied by the device manufacturer.

**WARNING!** Nonobservance of the 30 cm (12 in) minimum distance to any part of the device (including cables) specified by the manufacturer and covered by the scope of this document when using portable RF communications equipment (such as mobile phones [GSM 800/900, iDEN 820, LTE] or portable radio [TETRA, GM-RS/FRS 460]), including peripherals such as antenna cables and external antennas). Degradation of the device performance. Observe the specified minimum distance for portable RF communications equipment.

**WARNING!** Interference to the live image caused by external sources (such as mobile phones [GSM 800/900, iDEN 820, LTE] or portable radio [TETRA, GM-RS/FRS 460]). Risk to the patient.

- > Remove all sources of interference
- > Observe a minimum distance of 30 cm (12 in)
- > Ensure that the monitor is correctly displaying the live image

**WARNING!** Interference to the live image caused by electromagnetic disturbances (white light and NIR fluorescence imaging mode). Risk to the patient due to altered image representation, system locked in FI mode, live image failure (e.g. delay, frozen image), image interference (e.g. flickering, stripes), altered image orientation.

- > Remove all sources of interference
- > Observe the recommended minimum distances according to the Information on Electromagnetic Compatibility
- > Ensure that the monitor is correctly displaying the live image
- > Ensure that the live image is displayed in the correct orientation (not mirrored)

### 2.2.4 Combination with Medical Electrical Equipment

The device can be combined with components from other manufacturers provided that all components are compliant with the medical electrical equipment safety requirements according to IEC 60601-1.

It is the operator's responsibility to check and make sure that the system is and remains fully operational.

When using devices from different manufacturers and when operating an endoscope and/or endoscopic accessories with electromedical devices, it must be ensured that the applied part is properly insulated (in accordance with the requirements for type BF/CF applied parts defined in IEC 60601-1).

Contact our Service Department if you have any questions regarding compatibility.

### 2.2.5 Display of Faults

#### Overheating

The device features an internal temperature sensor.

In the event of camera control unit overheating, the message "Excess temperature! Attention, Camera will switch off in two minutes!" appears on the screen.

If the menu is opened, it closes automatically when the message appears.

#### Incompatible USB storage device

Compatible USB storage devices must be formatted with the FAT32 file system, and may only contain image and video files.

The USB storage device must meet the USB 2.0 standard (or higher) and be capable of a minimum data transfer rate of 20 MB/s.

If an incompatible USB storage device is used when attempting to save a capture, the message "USB storage device in incorrect format" appears.

It is recommended to only use USB flash drives in conjunction with the devices listed under the scope of validity (*section 0.1*) of these instructions for use. USB hard drives (e.g. HDD or SSD) are not suitable.

To ensure compatibility of the USB storage device, it is recommended to make a test capture. Please note that the devices have no formatting capabilities.

#### Incompatible hardware for NIR FI mode

If components are connected that do not support NIR FI mode, the message „NIR FI mode not available: Connect compatible ...“ appears, when attempting to enter NIR FI mode.

## 3 Installation and Initial Operation

### 3.1 Safety Notices

#### **WARNING**

##### **Electrical connections installed improperly**

Risk of fire, short circuit, or electric shock

- > Make sure that the electrical connections are installed in accordance with the relevant national technical regulations

#### **WARNING**

##### **Medical electrical system installation**

Risk of fire, short circuit, or electric shock

- > Installation must only be performed by qualified personnel
- > If a number of electrical devices are being used in combination, observe IEC 60601-1
- > Non-medical devices that meet the applicable IEC safety standards can be connected only via a medical isolating transformer
- > Do not connect any additional non-medical devices to a medical electrical system
- > Signal lines from devices with a functional ground (FE) connection, which are connected to different branches of the mains power supply, must be galvanically isolated at both ends
- > Connect the devices only to a power supply with a protective ground (PE)
- > After installing a medical electrical system, it must be inspected according to IEC 62353

#### **WARNING**

##### **Electromagnetic interactions**

Possible impaired function and image quality

- > Observe the EMC brochure supplied
- > Make sure that neighboring devices are not causing reciprocal electromagnetic interference
- > Perform a post-installation function check

## **⚠ WARNING**

### **Use of multiple sockets**

Risk of fire, short circuit, electrical shock, reduced level of safety

- > Whenever possible, avoid the use of multiple sockets
- > If required, use medically approved multiple sockets
- > Never connect multiple sockets in series
- > Do not cover multiple sockets (leads to heat accumulation)
- > Do not place multiple sockets on the floor
- > Use traction relief
- > Devices should be connected to a shared multiple socket only if they are all intended to be used as part of the medical electrical system

## **3.2 Staff Qualifications**

**Qualified personnel** Personnel responsible for setting up or installing medical electrical equipment must possess the necessary qualifications and be familiar with the safety regulations in place at the installation location and the applicable official regulations.

## **3.3 Installation**

### **3.3.1 Setting up the Device**

*Observe the safety notices provided at the start of this chapter and refer to the enclosed EMC brochure.*

**WARNING!** Risk from installation in potentially explosive locations. Increased risk of fire and explosion in oxygen-enriched atmospheres. Device must be installed outside of potentially explosive locations and not in the vicinity of flammable substances.

**WARNING!** Risk resulting from electromagnetism. Possible malfunctions and image interference. Use only the cables supplied. Ensure that all devices operated in the vicinity meet the EMC requirements. Perform a post-installation function check.

**WARNING!** Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other devices should be observed to verify that they are operating normally.

1. Install the device.

**Install safely!**

When installing the device, ensure that:

- it is positioned on a level, non-slip surface that is able to support its weight,
- it is installed outside of the patient environment and the sterile area on a sufficiently stable support,

- it is protected from water drips and splashes,
- it will not be subject to vibration during operation,
- the ventilation slots are not obstructed,
- the main switch on the rear of the device is accessible.

2. Connect the **potential equalization terminal** on the rear of the device to the potential equalization terminal at the installation site.

Use a potential equalization conductor in accordance with DIN 42801 with at least a 4 mm<sup>2</sup> cross-section.

### 3.3.2 Connecting Monitors

*Multiple monitors can be connected. The main monitor should be connected to the best signal available. Please observe the instructions for use of the monitor.*

#### Connecting a 4K/UHD monitor

1. Connect the four **4K/UHD** ports on the rear of the device to the four 4K/UHD ports on the 4K/UHD monitor.

Be sure to observe the correct order of the ports.

Use the supplied BNC cables to do this.

Secure the cable

Secure the connection cables at both ends to prevent them from being pulled out inadvertently.

Route cables safely!

2. When routing the cables, make sure they do not present a trip or entanglement hazard.

Perform 4K/UHD settings

3. For the 4K/UHD signal to be correctly displayed, perform the following monitor menu settings:

Input: 4x 3G-SDI

Mode: Quad-View, square division

Color space: BT.2020

#### Connecting HD monitors

1. Connect one of the **DVI** ports on the rear of the device to the DVI input on the main monitor.

Use the supplied DVI cable.

2. If necessary, connect any additional monitors to the remaining video outputs.

Secure the cable

Secure the connection cables at both ends to prevent them from being pulled out inadvertently.

Route cables safely!

3. When routing the cables, make sure they do not present a trip or entanglement hazard.

### 3.3.3 Connecting a Keyboard

1. Connect the keyboard connection cable to the USB connection labeled **USB-Keyboard (USB keyboard)** on the rear of the device.
  2. When routing the cables, make sure they do not present a trip or entanglement hazard.
- Route cables safely!

### 3.3.4 Connecting Foot Switches

*Up to two single-pedal foot switches can be connected. Connect foot switches if you want to control the camera by foot switch.*

Note: Only the foot switch specified in [section 8.3](#) is compatible.

1. Connect a foot switch to the **Remote-In 1** connection of the controller.
  2. Connect a second foot switch to the **Remote-In 2** connection of the controller.
  3. Configure the foot switches as described in [section 4.1.10](#).
  4. When routing the cables, make sure they do not present a trip or entanglement hazard.
- Remote-In 1 input
- Remote-In 2 input
- Configure foot switches
- Route cables safely!

### 3.3.5 Connecting a Recording Device

*Observe the instructions for use supplied with the recording device.*

1. Depending on the type of signal supported by the recording device, connect either the **4x 3G-SDI** connections, a **DVI**, or a **3G-SDI** connection on the rear of the medical device to the corresponding inputs or input on the recording device.
2. If you would like to use the buttons on the camera head or a foot switch to launch recordings, connect the **Remote-Out** ports on the rear of the medical device to the corresponding inputs on the recording device.

Use the supplied 3.5 mm jack plug cable to do this.

Remote-Out 1 output  
for photos

Connect the **Remote-Out 1** port to the photo recording connection on the recording device.

Remote-Out 2 output  
for video recordings

Connect the **Remote-Out 2** port to the video recording connection on the recording device.

3. If necessary, modify the button function of the camera head.

Route cables safely!

4. When routing the cables, make sure they do not present a trip or entanglement hazard.

### 3.3.6 Connecting and controlling a Light Source

*Connect the light source and control it using the camera menu or the camera head buttons.*

Note: Only compatible light sources can be controlled using the camera menu or the camera head buttons.

1. Connect the MIS-Bus port on the camera control control unit to that on the light source. Use the supplied MIS-Bus cable.
2. Ensure the connections are secured on both ends.
3. Switch the light source on and off using the menu item **Light on/standby** in the camera main menu or control its brightness (*section 4.1.8*).

Configure camera  
head buttons

If you wish to do so using the camera head buttons, modify the camera head button function, if necessary.

### 3.3.7 Connecting the Device to the Mains

*Do not connect the device to the mains until you have completed all the previous activities. Observe the safety notices provided at the start of this chapter.*

*The use of an uninterruptible power supply (UPS) is recommended to protect against input power interruptions.*

**NOTICE!** Risk of condensation with extreme differences in temperature. A device stored in very cold conditions can result in condensation in the device from humidity in the air when there is a rapid change in temperature from cold to warm. Wait until the device has come to the ambient temperature to connect to power supply and switch on.

1. Check that the supply voltage at the site of operation corresponds to that specified on the device nameplate.
2. Insert the power cord into the mains connection socket on the rear of the device.

**WARNING!** Sparks in a potentially explosive location. Risk of explosion. Ensure that the device is not situated in a potentially explosive location when connecting to the power supply.

3. Connect the power cord to the power supply.

Route cables safely!

4. When routing the cable, make sure it does not present a trip or entanglement hazard.

### 3.3.8 Performing a Function Check

*When operating the device for the first time, check the installation and make sure the device is in full working order.*

Precondition: Installation must have been completed.

1. Switch ON at the main switch on the rear of the device.
  - > *The device switches to standby mode.*
2. Switch on the device and all neighboring devices.
3. Connect a camera head as described in [section 5.6.1](#).
4. Perform white balance as described in [section 5.6.3](#).
5. Direct the camera head at an object in the room and focus the image.
6. Check that a good image is displayed on the monitor without vertical lines, color variations, or flickering.

Connect camera head

Set white balance

Focus the image

#### **Problems with the image displayed on screen?**

Electromagnetism can cause interference.

Try rearranging the different devices to increase the distance between them. Contact our Service Department if you require further assistance.

7. Check that the buttons on the camera head can be operated properly.

## 4 Configuration

### Examples of actual menu paths


Figure 4-1: Setting brightness.


Figure 4-2: Setting the system language.


Figure 4-3: Entering patient data.

## 4.1 Basic Settings

### 4.1.1 Setting System Language

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **Options** menu item and select.

3. Navigate to **System setup** and select.

> *The setup menu appears on the screen.*

4. Navigate to **Configuration** and select.

Language selection

5. Navigate to **Language selection** and select.

6. Select the desired language and confirm.

> *Your selection is displayed in white type.*

> *The language is changed as desired.*

### 4.1.2 Changing the On-Screen Display Position

*To change the position of the on-screen display on the monitor, choose an aspect ratio of 16:9 or 16:10, depending on the monitor used.*

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **Options** menu item and select.

3. Navigate to **System setup** and select.

> *The setup menu appears on the screen.*

4. Navigate to **Configuration** and select.

16:9 / 16:10

5. Navigate to **16:9 / 16:10** and select.

> *Your selection is displayed in white type.*

6. Select the desired aspect ratio and confirm.

### 4.1.3 Setting Date and Time

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > *The main menu appears on the screen.*
2. Navigate to the **Options** menu item and select.
3. Navigate to **System setup** and select.
  - > *The setup menu appears on the screen.*
4. Navigate to **Configuration** and select.
5. Navigate to **Date/time** and select.
6. Select the desired display format and confirm.
  - > *Your selection is displayed in white type.*
7. Set the date and time and save in each case by pressing Enter.

Date/time

### 4.1.4 Defining File Format for Photos (optional)

Define what file format should be used for photos taken. Select between BMP and JPEG.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > *The main menu appears on the screen.*
2. Navigate to the **Options** menu item and select.
3. Navigate to **System setup** and select.
  - > *The setup menu appears on the screen.*
4. Navigate to **Configuration** and select.
5. Navigate to **Image Format** and select.
6. Select the desired file format and confirm.
  - Select BMP if you want to create uncompressed photos in higher quality. One image can be up to 6 MB in size.
  - Select JPG if you want to create compressed photos in good quality with a smaller file size of approx. 400 KB.
  - > *Your selection is displayed in white type.*

Image format

#### 4.1.5 Defining Image Refresh Rate

*Leave the factory-preset image refresh rate unchanged. This setting only needs to be changed in exceptional circumstances.*

*Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.*

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **Options** menu item and select.

3. Navigate to **System setup** and select.

> *The setup menu appears on the screen.*

4. Navigate to **Configuration** and select.

50Hz / 60Hz

5. Navigate to **50Hz / 60Hz** and select.

6. Select the desired image refresh rate and confirm.

> *Your selection is displayed in white type.*

> *A dialog window for confirmation opens up.*

7. Confirm the dialog window using **Perform**.

#### 4.1.6 Defining On-Screen Displays

*Determine whether the date, time and if necessary patient data should be overlaid onto the screen. Also determine whether the configuration of the camera head buttons should be displayed temporarily after switching on the camera system.*

*Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.*

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **Options** menu item and select.

3. Navigate to **System setup** and select.

> *The setup menu appears on the screen.*

4. Navigate to **Configuration** and select.

OSD

5. Navigate to **OSD** and select.

6. Select the desired on-screen displays and confirm.
  - > *Your selection is displayed in white type.*

When selected,

- the date and time are overlaid onto the screen,
- the patient data are overlaid onto the screen,
- the configuration of the camera head buttons is displayed after switching on the camera system for several seconds.

#### 4.1.7 Switching Signal Sounds for taking Photos on or off (optional)

*Determine whether an acoustic signal should confirm when photos are taken.*

*Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.*

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > *The main menu appears on the screen.*

2. Navigate to the **Options** menu item and select.

3. Navigate to **System setup** and select.

- > *The setup menu appears on the screen.*

4. Navigate to **Configuration** and select.

Beeper

5. Navigate to **Beeper** and select.

6. Activate or deactivate the signal sound and confirm.

- > *Your selection is displayed in white type.*

#### 4.1.8 Controlling the Amount of Light

Compatible light sources allow for the use of the auto light control feature. To this purpose, auto light control must be activated.

##### Activate automatic light control

1. Using the MIS-Bus cable, connect the camera to a compatible light source. Doing so, observe the instructions for use supplied with the light source.

2. Activate auto light control via the camera menu.

- > *When auto light control is activated, this message appears on the light source display.*

### Setting the light intensity manually

If you wish to control the amount of light manually and the light source used is compatible and appropriately connected, you can do so using the camera menu.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > The main menu appears on the screen.
2. Navigate to **Light source +/-** menu item and select.
3. Select a light intensity of 70 %-100 % with an LED light source and confirm.  
Reduce the light intensity to reduce blooming. Increase the light intensity to reduce image noise.

Note: With LED light sources, too low a light intensity can change color representation.

#### 4.1.9 Setting Shutter Speed

The shutter speed is set automatically. Setting a fixed shutter speed can be beneficial for special applications, e.g. for minimizing the effects of vibration during stroboscopy.

Set a fixed shutter speed if required.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > The main menu appears on the screen.
2. Navigate to the **Options** menu item and select.
3. Navigate to **System setup** and select.
  - > The setup menu appears on the screen.
4. Navigate to **Configuration** and select.
5. Navigate to **Shutter** and select.
6. Select "Auto" or "Manual" and confirm.  
If you choose "Manual", set the desired shutter speed and confirm.
  - > Your selection is displayed in white type.

#### 4.1.10 Configuring Foot Switches

*Assign a function to foot switches as desired.*

*Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.*

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **Options** menu item and select.

Foot switch

3. Navigate to **Foot switch** and select.

4. Assign the foot switch function as desired.

Select one of the following functions:

- NIR FI mode
- Selective color enhancement (SCE)
- Smoke reduction
- Light on/standby
- Photo [capture] (onto a USB storage device)
- Video [capture] (onto a USB storage device)
- Color Shift 1 / 2
- Remote Out 1 (output Remote-Out 1)  
(capturing of photos through an external recording device)
- Remote Out 2 (output Remote-Out 2)  
(capturing of videos through an external recording device)

Note: The selectable image settings are described in [section 4.2.1](#) in more detail.

#### 4.1.11 Displaying System Information

*The system information can be useful to the Service Department.*

*Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.*

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **Options** menu item and select.

3. Navigate to **System setup** and select.

> *The setup menu appears on the screen.*

System information

4. Navigate to **System information** and select.

5. The desired information can be displayed:

- Software version
- Hardware version
- Serial number > Controller, Serial number > Camera head and Serial number > Light source

## 4.2 Image Settings

### 4.2.1 General Information

Optimized camera settings can be selected in the configuration menu for selected endoscopic specialties.

The image settings can be adjusted and saved according to your requirements.

#### Adjustment of image settings

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. To reach the image setting you wish to adjust, navigate to the menu item **Image settings** and select.

3. Select **Advanced Views** or **Customize image**.

> *A list of the image settings described in the following opens up.*

Select **Image presets** if you would like to access predefined parameter sets for certain specialties or color modes ([section 5.6.2](#)) or work with user-defined profiles ([section 4.2.5](#)).

4. Navigate to the desired image setting and select.

> *An additional window appears on the screen next to the main menu.*

5. Make the desired adjustments using the arrow buttons.

Save your changes by pressing the Enter button each time.

### 4.2.2 Advanced Views

Note: The menu item **Advanced Views** is not available when in NIR FI mode.

The following configurable parameters are available in the **Advanced Views** menu:

#### Smoke reduction

Smoke reduction allows a clearer image in a smoke-filled environment (e.g. when there is coagulation with electrosurgical units).

Use smoke reduction to compensate for a poor view due to smoke formation.

Note: Smoke reduction creates a very contrast-rich image and can promote blooming.

### **Selective color enhancement (SCE)**

With Selective color enhancement, you can intensify the predominant color in the image.

Use Selective color enhancement to better differentiate structures from one another in environments with little differences in color, for example, to accentuate the colors of vessels.

Note: Selective color enhancement changes the natural color representation. For this reason, the message “Attention: false color rendering” is displayed on the screen as long as this setting is activated.

### **Color Shift 1 / 2**

Color Shift 1 / 2 changes the color of vessels and tissue structures.

Use Color Shift 1 / 2 for a better view of vessels and tissue structures.

Note: Color Shift 1 / 2 alters the natural color representation. For this reason, the message “Color Shift 1” or “Color Shift 2” is displayed on the screen as long as this setting is activated.

### **Grid removal**

Grid removal reduces or removes the grid that is caused by use of image bundle systems.

Use grid removal to suppress interfering grids.

Note: Grid removal causes a slight reduction in image sharpness.

## **4.2.3 Customizing the Image**

The following other parameters are available in the **Customize image** menu:

Note: The following parameters affect only the white light image and do not work in NIR FI mode:

- Brightness
- Contrast
- Edge enhancement
- Color saturation
- Noise reduction

### **Zoom**

The digital zoom allows up to three-fold digital enlargement of the image displayed.

Note: The digital zoom reduces the image resolution and decreases the image quality. For this reason, use the optical zoom on the camera head preferably whenever possible.

### **Brightness**

Brightness controls the exposure control for the camera.

Reduce brightness to decrease blooming and increase it to make the image brighter.

### **Contrast**

Contrast changes the contrast level of the image.

Reduce the contrast to obtain a homogeneously illuminated image for endoscopes with poorly illuminated edges.

Note: Dimming on the edges of the image can also result from the window size adjustment being too small.

### **Edge enhancement**

Edge enhancement amplifies the contours at the transitions of adjacent structures.

Increase the edge enhancement during laparoscopic and arthroscopic applications with rod lens systems to improving the view of vessels.

Reduce the edge enhancement with image bundle systems to decrease unwanted structures of the image bundle system.

Note: Edge enhancement can increase image noise.

### **Color saturation**

Saturation changes color intensity.

Change the saturation according to your personal impression. Reduce the saturation to decrease unwanted color effects, such as moiré effects.

### **Window**

Window size defines the image section used by the camera for controlling image brightness.

Choose a window that is similar in size or smaller than the round endoscope image.

Choose an even smaller window to achieve a better view in the distance in tubular cavities, such as the urethra.

Note: If the window setting is too large and exceeds the round endoscope image, blooming can result.

### **NIR FI Window**

Use NIR FI Window to define the image section used by the camera for controlling image brightness in NIR FI mode.

### Noise reduction

Noise reduction reduces noise in the image.

Adjust noise reduction at the end or when a high edge enhancement has been set and further reduction of noise is desired.

Note: Noise reduction leads to a slight reduction in image sharpness starting at the middle setting. A higher setting may cause artifacts with quick movements (streaking).

### Picture out of Picture (PoP)

Picture out of Picture displays the image adjusted by video algorithms on the right-hand side of the original image.

The menu item Picture out of Picture allows you to choose the video algorithm to be used for the display of both images next to one another for comparison purposes.

Picture out of Picture can also display different NIR FI modes simultaneously for comparison purposes.

Note: The original image is displayed on the left, the adjusted image is displayed on the right. When selecting Advanced Views or NIR FI modes, both live images show the same image section. When selecting Zoom, the center of the original image is displayed in enlarged format.

#### 4.2.4 Changing the Image Settings of Specialties and Color Presets

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. To reach the specialty or color mode you wish to change, navigate to the **Image settings** menu item and select.

3. Navigate to the **Image presets** menu item and select.

4. Navigate to the **Specialty** menu item and select.

Then navigate to the desired specialty or color mode and select.

> *The specialty or color mode is changed.*

5. Navigate to **Image settings** and select.

6. Navigate to **Advanced views** or **Customize image** and select.

7. Make the desired changes.

Save your changes by pressing the Enter button each time.

8. Set the white balance for a realistic view of the changes you have made.

#### 4.2.5 Making and saving User-Specific Image Settings as a Profile

Create a user-specific profile based on the existing reference profile and save it under a name of your choosing.

Note: A keyboard must be connected to enter a user-specific name.

User-specific profiles can also be saved without a keyboard. In this case, your settings are saved under the preset name, consisting of the camera head model and Surgeon 1-10 (e.g. FVHD2 Surgeon 1, FVHD2 Surgeon 2, ...).

Up to 10 user-specific profiles can be saved.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> The main menu appears on the screen.

2. Navigate to the **Image settings** menu item and select.

3. Navigate to the **Image presets** menu item and select.

4. Navigate to the **Saved users** menu item and select.

5. Navigate to the **Create User Profile** or **Edit User Profile** menu item and select.

Only when creating a user profile for the first time: Navigate to the **Enter profile name:** menu item and select.

6. Navigate to the **Rename** menu item and select.

Enter name

7. Press the right arrow button on the keyboard and enter a name for the user-specific profile.

Up to 30 characters are available for the name.

8. Confirm the entry with the Enter button on the keyboard.

> The user-specific profile is renamed with the name given.

9. Navigate back in the main menu to **Reference specialty** and select.

Select reference profile

10. Select the reference profile that you wish to base the user-specific profile on and confirm.

11. Navigate back in the main menu to **Image settings** and select.

12. Make the desired changes.

Save your changes by pressing the Enter button each time.

## 4.3 Reset to Default Settings

### 4.3.1 Resetting the Device to Default Settings

**NOTICE!** All settings are reset to the default settings. If necessary, only restore the default settings of an individual specialty.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> The main menu appears on the screen.

2. Navigate to the **Options** menu item and select.

3. Navigate to **System setup** and select.

> The setup menu appears on the monitor.

4. Navigate to **Configuration** and select.

Default settings

5. Navigate to **Default settings** and select.

6. Select **Perform** and confirm.

> A dialog window for confirmation is displayed on the screen. The option **Perform** is activated.

If you wish to cancel the process, select **Cancel** and confirm.

> The dialog window closes and the setup menu is shown.

If you wish to reset the device to the default settings, confirm the selection.

> The device is reset to the default settings and restarts.

### 4.3.2 Resetting the Image Settings of a Specialty

Follow the steps below if you want to reset the image settings of an individual specialty to the default settings.

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> The main menu appears on the screen.

2. Select the specialty you want to reset, as described in [section 5.6.2](#).

Settings

3. Navigate to **Image settings** and select.

4. Navigate to **Customize image** and select.

5. Navigate to **Reset?** and select.

6. Navigate to **Perform** and select.

## 5 Operation

### 5.1 Safety Notices

#### **WARNING**

##### **Modifications to the installation**

Risk of fire, short circuit, and electric shock

- > Follow the safety notices in the Installation and Initial Operation chapter
- > Do not make unauthorized changes to the electrical installation
- > If the installation contains a multiple socket, do not connect additional devices to this without authorization
- > Never connect multiple sockets in series

#### **WARNING**

##### **Current leakage when touching the patient**

Risk of electric shock

- > Never touch the patient and the device housing or accessible contacts of electric devices at the same time
- > Never touch the patient and non-medical electric devices at the same time

#### **WARNING**

##### **Positioning the equipment cart inside the sterile area**

Risk of infection

- > Do not position the equipment cart inside the sterile area
- > Measure the cable lengths so there is enough room for movement without compromising sterility
- > The light guide must be anchored securely within the direct OR environment to avoid any risk of accidental slipping

#### **WARNING**

##### **Displaying a recorded image instead of the live image**

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera at all times

**⚠ WARNING**

**Live image failure during operation of a defibrillator**

Risk to the patient

- > Make sure that 3 seconds (recovery time) after operation of a defibrillator the live image is displayed correctly again
- > Restart the controller if the live image is incorrectly displayed

**⚠ WARNING**

**Interactions between devices in simultaneous use (e.g. lasers, electro-surgery)**

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the required type BF, CF or CF defibrillation-proof classification requirements according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active HF electrodes
- > Do not activate HF electrodes in the vicinity of flammable gases or liquids
- > Before using HF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

**⚠ WARNING**

**Addition of leakage currents. When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive**

Risk to the patient

- > Use endoscopes in conjunction with endotherapy devices only if both have matching type applied parts (CF or CF defibrillation-proof)

**⚠ CAUTION**

**Electrostatic discharge at the connection field of the monitor**

Temporary image failure, risk to patient

- > Do not touch the connection field and connections at the rear of the monitor during ongoing procedures

## 5.2 Staff Qualifications

This device may be operated only by physicians and medical assistants with the requisite training, advanced training, knowledge, and practical experience in the endoscopic discipline concerned as defined by the provisions in place at the site of operation.

## 5.3 General Notes on Operation

Start every procedure by creating a new patient. This ensures that the captures created during the procedure are saved in a separate file folder.

Enter the patient data for this as described in [section 5.6.4](#).

**New patient** Find the **“New patient”** entry in the menu. Confirming this entry deletes the existing patient data. If you do not wish to enter any patient data, confirm this entry anyway so that a new file folder is generated on the USB storage device for the subsequent captures.

**Sufficient disk space?** Ensure that a USB storage device is connected with sufficient disk space for your recordings.

Depending on the components used, navigation in the menu can be done with the buttons on the front side of the device, the camera head buttons or via a connected keyboard.

Operation on the front side of the device is described in the following.

## 5.4 Technical Inspection Prior to Use

### 5.4.1 Performing a Visual Inspection

*A visual inspection must be performed **prior to each medical procedure**.*

**Properly installed?** 1. Check that the device is installed as described in [section 3.3.1](#). Make sure in particular that the potential equalization is properly connected.

Do not use the device if it is not properly set up and connected.

**Is the housing undamaged?** 2. Inspect the housing for external damage.

Do not use the device if there is any external damage to the housing.

**Is the cable undamaged?** 3. Check that the mains cable is undamaged.

Do not use the device if the cable is damaged.

**Is the socket dry and undamaged?** 4. Check that the contacts on the camera cable connection socket are free from moisture and dirt.

Do not use the device if there is any moisture in the socket or contacts are soiled.

### 5.4.2 Performing a Function Check

*A function check must be performed prior to each medical procedure.*

- |                      |  |
|----------------------|--|
| Connect camera head  | 1. Connect a camera head as described in <a href="#">section 5.6.1</a> .   |
| Switch on the system | 2. Switch on all the system components that you intend using for the procedure.  |
| Set white balance    | 3. Perform white balance as described in <a href="#">section 5.6.3</a> .   |
| Focus the image      | 4. Direct the camera head at an object in the room and focus the image.<br>Check that you are able to achieve a sharp, good-quality image.<br>Do not use the system if the image contains vertical lines, color variations, or flickers, or you are unable to achieve a sharp, good-quality image. |
|                      | 5. Adjust the monitor as required to achieve an optimum viewing angle for the surgeon.   |

## 5.5 Switching On and Off, Separation from the Power Supply

### Switching on the device

- Switch on the device at the on/off button on the rear.  
> *The device switches to standby mode.*

- Press the ON button.  
> *The ON switch illuminates in white.*

### Switching off the device

- Switch off the device at the front.  
> *The device switches to standby mode.*

To switch off the device completely, switch off at the main switch on the rear.

### Separating the device from power

- To separate the device completely from its power supply, disconnect the power cord. The appliance coupler is considered as disconnecting device.

Make sure that the rear of the device is accessible at any time.

## 5.6 Preparation for Use

### 5.6.1 Connecting up a Camera Head

- Slide the connection plug fully into the connection socket on the device until it clicks into place.  
> *The name of the connected camera head model is displayed on the screen.*

### 5.6.2 Selecting a Specialty or Color Preset

Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **Image settings** menu item and select.

3. Navigate to the **Image presets** menu item and select.

4. Navigate to the **Specialty** menu item and select.

> *A list of the following specialties and color modes opens up:*

Select appropriate settings

- **Laparoscopy 1** (good visualization of small structures and fine details on medium-sized monitors up to 32 inch diagonal)

- **Laparoscopy 2** (very finely detailed, natural image. Ideal for the use of particularly large 4K/UHD monitors)

Select one of the 4K/UHD Laparoscopy settings according to your personal preference.

Note: For the most optimal image, use 4K/UHD endoscopes.

- Color mode 1

- Color mode 2

- Color mode 3

- Urology

- URO Fiberscope

- Gynecology

- GYN Fiberscope

- Arthroscopy

- ENT rigid

- ENT flexible

- Neuroendoscopy

Navigate to the desired specialty or color mode and select.

> *The parameters of the current specialty or color mode are taken over.*

### 5.6.3 Setting the White Balance

*Set the white balance **prior to each medical procedure.***

Note: Keep in mind that ambient light that deviates significantly from the color temperature of the light source can affect the results of the white balance.

1. Connect the endoscope that you intend on using for the procedure.
2. Connect the light guide and light source that you intend on using for the procedure.
3. Switch on all system components.
4. Direct the endoscope at a white object, a sterile wipe for example, that is approximately **5 cm** away and focus the image.

5 cm away

5. Press the **White balance** button on the front of the camera control unit for about 2 seconds or use the designated button on the camera head.

> *The message **White balance OK** appears on the screen.*

---

#### **Does the message “White balance failed” appear on the screen?**

The white balance was not performed as described above.

Repeat the white balance.

---

Note: White balance cannot be set when NIR FI mode is activated.

### 5.6.4 Entering Patient Data

*A keyboard must be connected to enter a patient data.*

*Navigate in the menu with the arrow buttons and confirm your selection each time using the Enter button.*

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > *The main menu appears on the screen.*
2. Navigate to the **Patient information** menu item and select.
  - > *The menu item “New patient” is activated.*
3. Press the right arrow button to delete the existing patient data.
4. Navigate to the first field for data entry using the down arrow button.
5. Enter the data. Jump from field to field with the **Enter** button.
6. Navigate to **Confirm selection** and save your entries using the right arrow button.
  - > *With the first recording, a new file folder is generated on the USB storage device, in which the captures are saved.*

## 5.7 NIR Fluorescence Imaging Operation

### 5.7.1 Toggling between NIR FI and White Light Imaging Mode

For NIR fluorescence visualization, press the affected button on the camera head or on the affected foot switch to toggle between the imaging modes in the following order:

#### White light imaging:

A standard white light reflectance image is displayed in natural color reproduction. No NIR fluorescence can be seen.

#### NIR FI Green:

A white light image is displayed in grayscale with NIR fluorescence appearing overlaid in green.

#### NIR FI Contrast:

NIR fluorescence is displayed in white and all other areas of the image appear dark.

#### NIR FI Intensity:

A white light image is displayed in levels of gray with NIR fluorescence appearing overlaid in green based on a multi-color intensity scale.

**Note:** The NIR fluorescence color appearance in NIR FI intensity mode depends on the position and distance of the endoscope tip relative to the tissue and does not provide quantitative information with regard to the NIR fluorescence signal strength.

### **WARNING**

#### **Anatomic background is not visible or inadequately displayed**

Risk of injury to the patient

Do not manipulate tissue if NIR FI Contrast mode is activated or if backlight intensity is reduced

#### Activating NIR FI mode in the menu

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.

> *The main menu appears on the screen.*

2. Navigate to the **NIR-FI mode** menu item and select.

> *The menu item „NIR FI mode“ is activated.*

3. Navigate to the desired imaging mode and select.

> *The main screen shows a symbol indicating that the NIR FI mode is on, next to the “LASER ON” symbol indicating that the NIR laser is on.*

### Deactivating NIR FI mode in the menu

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > *The main menu appears on the screen.*
2. Navigate to the **NIR-FI mode** menu item and select.
  - > *The menu item „NIR FI mode“ is activated.*
3. Navigate to the **OFF** menu item and select.
  - > *The white light reflectance image appears on the main screen.*

### 5.7.2 Adjusting the Image in NIR FI Mode

Note: In NIR FI mode the signal intensity varies depending on the working distance between endoscope and tissue. The closer the endoscope is to the tissue, the more intense the fluorescence signal appearance.

Note: If low levels of the fluorescence dye are not readily visible, reducing the amount of the LED white light may increase visibility.

1. Press the **Menu** button on the front side of the device for approximately 2 seconds.
  - > *The main menu appears on the screen.*
2. Navigate to the **NIR-FI mode** menu item and select.
  - > *The menu item „NIR FI mode“ is activated.*
3. Navigate to the desired image setting and make the desired adjustments using the arrow buttons.

Save your changes by pressing the Enter button each time.

#### NIR FI Gain

NIR FI Gain affects the intensity and discernability of the fluorescence signal appearance.

To improve the visibility of a weak fluorescence signal, increase NIR FI Gain, making the fluorescing appearance more intense.

Reduce NIR FI Gain to make the fluorescing appearance less intense.

#### NIR FI Brightness

NIR FI Brightness regulates the exposure control for the camera and affects the brightness appearance of the white light image.

## 5.8 Creating Recordings

**CAUTION!** Recorded photos and videos can be compressed or changed by modified image settings. The photo and video documentation is not suitable for making a diagnosis.

### 5.8.1 Recording a Photo

1. Ensure that a USB storage device is connected.
2. Activate the release.

Dependent on model, configuration and connected accessories

- press the **Photo capture** button on the front side of the device,
- press the affected button on the camera head,
- activate the affected foot switch,
- release the capture in the menu.

- > *The indicator on the main screen shows the photo captured for approx. one half second (freeze function).*
- > *The photo is saved to the USB storage device.*

### 5.8.2 Recording a Video

1. Ensure that a USB storage device is connected.
2. Activate the release to start the recording.

Dependent on model, configuration and connected accessories

- press the **Video capture start / stop** button on the front side of the device,
- press the affected button on the camera head,
- activate the affected foot switch,
- release the capture in the menu.

- > *The main screen shows a symbol indicating that the video capture has been started.*
- > *The video is saved to the USB storage device.*

3. Activate the release once more to exit the recording.

### 5.8.3 Recording a Photo or Video using an External Device

*Starting and stopping video recordings can only be controlled with the camera head buttons or via foot switch. The camera head or foot switch must be configured appropriately. For video capture, a recording device must be connected.*

*Observe the instructions for use supplied with the recording device.*

1. Ensure that a recording device is connected and switched on.
2. Activate the release to start the recording.

Dependent on configuration and connected accessories

- press the affected button on the camera head,
- activate the affected foot switch.

> *The photo or video recording is saved to the external device.*

3. Video recording only: Activate the release once more to exit the recording.

## 5.9 Storage and Transport

Store the device in a dry, well-ventilated, temperature-controlled room, where it is protected from dust.

Switch off the device prior to storage and remove the power cord and accessories.

When placing into storage or storing temporarily, make sure that the device cannot be damaged by its surroundings.

To avoid damage, keep out of direct sunlight and away from radioactivity and strong electromagnetic radiation.

The ambient conditions for transport, storage, and operation are described in [section 8.2](#).

**6 Reprocessing**  
**6.1 Safety Instructions**

**⚠ WARNING**

**Risks associated with improper cleaning**

Risk of fire, short circuit, and electric shock

- > Before cleaning, switch off the device at its rear main switch and disconnect from the mains power supply
- > Make sure no moisture penetrates the device
- > The device should not be connected to the mains power supply until it has completely dried

***NOTICE***

**Improper cleaning**

Risk of damage to the device

- > Do not use sponges or wipes that are dripping wet
- > Do not use abrasive cleaners or solvents

## 6.2 Performing Cleaning and Disinfection

1. Switch off the device at its rear main switch.

**NOTICE!** Pulling the power cable will damage it. Pull at the plug to disconnect device from power.

2. Separate the device from the power supply by pulling out the plug of the power cord from the rear of the device.
3. Use disinfectant wipes<sup>2</sup> to clean all contaminated areas of the device for at least **1 minute** and to completely remove visible soiling.
4. Wipe the product until no more soiling can be seen. Use additional wipes, if required.
5. Use a fresh wipe to thoroughly wet all contaminated surfaces, ensuring that the device stays visibly wet for **3 minutes**. Use additional wipes, if required.
6. Make sure you thoroughly wipe the sealing ring and the protective cap of the USB connection. Disinfecting this area can be challenging.
7. Open the protective cap of the USB connection.
8. Disinfect all parts of the protective cap and the sealing ring by gently wiping it.
9. Continuously wipe the device during the **3 minutes**. Check that all seams, creases, crevices, recessed areas, and mated surfaces come into contact with the disinfectant.
10. Wait until all cleaned parts have thoroughly dried.
11. Reconnect the device to the power supply.

<sup>2</sup> Formulation: quaternary ammonium, bactericidal, virucidal, suitable for use on equipment.  
For cleaning and disinfection validation Sani-Cloth Plus disinfectant wipes manufactured by PDI® were used.

## 7 Maintenance and Repair

### 7.1 Safety Notices

#### **⚠ WARNING**

##### **Improper maintenance and repair**

Risk of fire, short circuit, and electric shock

- > All maintenance and repair work must only be performed by qualified personnel
- > Do not make any modifications to the device

### 7.2 Staff Qualifications

**Qualified personnel** Personnel responsible for maintaining the device or carrying out safety inspections must possess the necessary qualifications and be familiar with the local safety regulations in place at the installation location and the applicable official regulations.

Repairs must be performed by authorized service technicians only.

### 7.3 Safety Inspection

We recommend to arrange for safety inspections to be performed every 12 months in accordance with Section 5 of IEC 62353.

The inspection comprises the following steps:

- visual inspection
- measuring the protective resistance (only for protection class I devices)
- measuring the leakage currents
- measuring the insulation resistance

**Document results** The results of the inspection should be documented in an inspection report in accordance with Section 6 and Annex G of IEC 62353.

We recommend that all components of a medical electrical system be subjected to documented technical safety inspections at regular intervals.

## 7.4 Changing the Battery

1. Isolate the device from its power supply and disconnect the power cord.

**NOTICE!** Pulling the power cable will damage it. Disconnect the power at the plug.

2. Loosen the screws of the battery compartment.


Figure 7-1: Changing of the battery.


3. Pull out the battery compartment.
4. Replace the battery (VARTA CR2032).
5. Gently slide the battery compartment back in and secure it with the screws.

## 7.5 Changing the Fuse

1. Isolate the device from its power supply and disconnect the power cord.

**NOTICE!** Pulling the power cable will damage it. Disconnect the power at the plug.

2. Release the fuse holder with a screwdriver.

*Figure 7-2: Releasing the fuse holder with a screwdriver.*

3. Carefully pull out the fuse holder.

**WARNING!** Unsuitable electrical fuse. Risk of fire. Insert the fuse specified in [section 8.3](#) only and never bypass the fuse.

4. Replace the faulty fuse.
5. Carefully slide the fuse holder back in until it clicks into place.

## 7.6 Repair

Should you need to arrange a repair for the device, please contact one of our subsidiaries. Contact details can be found on the back of these instructions for use.

When you send in equipment, please enclose as accurate a **fault description** as possible, and record the item number and serial number of the product on the delivery note. You will find these details on the nameplate.

### **⚠ WARNING**

#### **Contaminated device**

Risk of infection

- > Reprocess the device prior to shipping (*chapter 6*)
- > Clearly identify contaminated products as such

Send only thoroughly cleaned goods in for repair.

For shipping, please use the original packaging whenever possible. Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

## 7.7 Troubleshooting

Issue	Possible causes	Fixes
On-screen display set off	Aspect ratio set incorrectly	Correct the position of the on-screen display ( <a href="#">section 4.1.2</a> )
Image too dark, too little illumination	Glass surfaces on camera head are soiled	Clean glass surfaces
	Stubborn residue on the glass surfaces	Remove residue
	Light guide defective	Connect a new light guide. Send in for repair
	Camera brightness set too low	Increase camera brightness
	Endoscope optical system defective	If the image is too dark even without the camera, use a different endoscope and send in the endoscope for repair
	Light output set too dark	Increase the light output at the light source
4K/UHD lacks color	4K/UHD monitor set incorrectly	Set monitor correctly ( <a href="#">section 3.3.2</a> )
4K/UHD image is flawed	4K/UHD monitor set incorrectly	Set monitor correctly ( <a href="#">section 3.3.2</a> )
No image on the monitor	No power supply	Connect and switch on devices ( <a href="#">chapter 3</a> )
	Fuse defective	Replace fuse ( <a href="#">chapter 7</a> )
	Camera control unit connected to monitor incorrectly	Connect camera control unit to monitor correctly ( <a href="#">chapter 3</a> )
	Video cable defective	Replace video cable
	Connection cable not connected or defective	Connect the camera head to the camera control unit. Check the plug for moisture. If the cable is defective, send it in for repair
	Monitor not tuned in to the correct input signal	Set monitor for the correct input signal
Date / time not correct	Date / time not correctly set	Set the date / time correctly in the operating menu of the camera
	Battery empty or defective	Change the battery ( <a href="#">section 7.4</a> )

Issue	Possible causes	Fixes
Light source switching on/off on its own	MIS-Bus fault in the camera	Disconnect camera and light source, send in camera for repair
	MIS-Bus fault in the light source	Disconnect camera and light source, send in light source for repair
Auto light control not possible	Light source does not support auto light control	Use a compatible light source model
NIR FI mode not possible	Light source, camera head or light guide does not support NIR FI mode	Use a compatible light source, camera head or light guide

Table 7-1: Troubleshooting table.

Message	Possible causes	Fixes
“Camera head not connected” “Incompatible camera head”	Camera head not connected to camera control unit Camera head not compatible with camera control unit	Connect compatible camera head to camera control unit. For combination options see the respective instructions for use of the camera heads
“Excess temperature! Attention, Camera will switch off in two minutes!”	Camera control unit ventilation slots are covered	Uncover ventilation slots
	Camera control unit fan defective	Send camera control unit in for repair
“No light source connected / Light source incompatible”	Light source not compatible with camera control unit	Connect compatible light source to camera control unit, using a MIS-Bus cable
“No USB storage available”	No USB storage device connected to camera control unit	Connect USB storage device to front side of the device
“USB storage device in incorrect format”	USB storage device with incompatible file system	Format USB storage device using the FAT32 file system
“USB device connected to wrong port”	USB device inserted on the rear of the device	Insert USB device on the front side of the controller
“USB device is read-only” “Unknown USB device error”	USB device defective	Replace USB device

Message	Possible causes	Fixes
“USB device removed while saving”	USB device removed while saving	Remove USB device only after saving has concluded
“White balance failed”	White balance not performed as described in the instructions for use	Repeat white balance ( <a href="#">section 5.6.3</a> )
„Function not available: Leave NIR FI mode before repeating your command“	Function activated that is only available in white light imaging mode	Repeat function in white light imaging mode
„NIR FI mode not available: Connect compatible light source“	Light source guide does not support NIR FI mode	Connect compatible light source
	MIS-Bus connection fault in the light source or camera	Check MIS-Bus connection
„NIR FI mode not available: Connect compatible camera head“	Camera head does not support NIR FI mode	Connect compatible camera head
„NIR FI mode not available: Connect compatible light guide“	Light guide does not support NIR FI mode	Connect compatible light guide

Table 7-2: Troubleshooting table.

## 8 Product Data

### 8.1 Technical Data

Dimensions (W x H x D)	295 x 129.5 x 355 mm
Weight	approx. 7 kg
Current consumption	1.2 - 0.5 A
Supply voltage	100-240 V~, 50/60 Hz
Appliance class as per IEC 60601-1	Class I
Device fuse	T2.0AH 250V
Image refresh rate	50/60 Hz (switchable)
Video outputs	1x 4K/UHD 3840 x 2160p (4x 3G-SDI) 2x DVI 1080p 2x 3G-SDI 1080p
Interfaces	1x MIS-Bus 2x Remote IN, 2x Remote OUT 2x USB
Degree of protection	IP 21
Classification as per IEC/CISPR 11	Group 1 Class A

### 8.2 Ambient Conditions

#### Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	50 kPa to 106 kPa

#### Operating conditions

Temperature	+10 °C to +40 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

#### **CAUTION**

##### **Non-observance of the ambient conditions**

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

> Observe operating conditions as well as transport and storage conditions

### 8.3 Spare Parts and Accessories

**CAUTION!** Incompatible spare parts or accessories. Risk to the patient resulting from possible malfunctions. Use original spare parts and accessories only.

Designation	Item number
Micro-fuse T2.0AH 250V, 5x20 mm	A075305
MIS-Bus cable, 2.25 m	A057635
MIS-Bus cable, 0.75 m	A059584
DVI cable, 2 m	A090275
3.5 mm jack plug cable, 1.2 m	A056998
BNC cable (number required for a 4K/UHD connection: 4)	BNC-BNC.2
Power cord (country-specific)	on request
USB keyboard, German layout	PS00510
USB keyboard, US layout	PS00511
Foot switch, 1-pole momentary, 3.5 mm jack plug, IP X7, 2 m	6226-ACAA-BZZZ-002
Battery VARTA CR2032	A052755

Table 8-1: Spare parts and accessories.

### 8.4 Information on Open-Source Licenses

The device uses open-source licenses. For more information on the terms governing their use, please refer to the supplied brochure specifying the open-source licenses used.

## 9 Disposal

### **WARNING**

#### **Contaminated device**

Risk of infection

> The device must be reprocessed prior to disposal (*chapter 6*)

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

The product is packed in a polyethylene (PE) bag and a PE foam inlay (0.50 kg), as well as a corrugated cardboard box (1.18 kg). Other packing materials used are a cardboard insert (0.22 kg) and a membrane pad made of PE film and corrugated cardboard (0.30 kg).

Any product carrying this symbol must be disposed of separately through dedicated electrical and electronic devices recycling. Within the EU, such disposal is taken care of by the manufacturer free of charge.

SCHÖLLY FIBEROPTIC GMBH  
Robert-Bosch-Str. 1-3  
79211 Denzlingen  
Germany