Gebrauchsanweisung FLEXILUX 200 LED

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Instructions For Use FLEXILUX 200 LED

Light Source for Medical Endoscopy

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	Content	Page
0	Important Information about this Document.....	43
0.1	Scope of Validity, Identification, Purpose.....	43
0.2	Target Group.....	43
0.3	Using and Storing this Document.....	43
0.4	Supplementary Documents.....	43
1	General Information about the Device.....	44
1.1	Scope of Delivery.....	44
1.2	Product Description.....	44
1.2.1	Performance Characteristics and Function.....	44
1.2.2	Visual Overview.....	45
1.2.3	Compatible Light Guides.....	46
1.3	Usage.....	47
1.3.1	Intended Use.....	47
1.3.2	Indications.....	47
1.3.3	Contraindications.....	47
1.4	Conformity.....	47
1.4.1	Standards and Directives.....	47
1.4.2	Medical Device Classification.....	47
1.5	Marking.....	48
1.5.1	Pictograms and Information on the Device and Packaging.....	48
1.5.2	Pictograms in this Document.....	49
1.6	Service Department Contact Details.....	49
1.7	Obligation to Report Serious Incidents.....	49
2	General Safety Information.....	50
2.1	Safety Messages in this Document.....	50
2.1.1	Safety Messages at the Start of a Chapter.....	50
2.1.2	Safety Messages in the Body of the Text.....	51
2.2	Product Safety.....	51
2.2.1	Basic Safety Notices.....	51
2.2.2	Staff Qualifications.....	52
2.2.3	Electromagnetic Compatibility.....	53
2.2.4	Combination with Medical Electrical Equipment.....	54
2.2.5	Display of Faults.....	54

3	Installation and Initial Operation.....	56
3.1	Safety Notices.....	56
3.2	Staff Qualifications.....	57
3.3	Installation.....	58
3.3.1	Setup.....	58
3.3.2	Connection to Mains Power.....	59
3.4	Initial Operation.....	59
3.4.1	Initial Operation and Function Check.....	59
4	Operation.....	60
4.1	Safety Notices.....	60
4.2	Staff Qualifications.....	61
4.3	Technical Inspection Prior to Use.....	62
4.3.1	Visual Inspection.....	62
4.3.2	Function Check.....	62
4.4	Switching On and Off, Separation from the Power Supply.....	63
4.5	Connection of a Light Guide.....	64
4.6	Decoupling the Light Guide.....	64
4.7	Setting the Brightness.....	64
4.8	Using the LED Standby Function.....	65
4.9	Storage and Transport.....	65
5	Cleaning.....	66
5.1	Safety Notices.....	66
5.2	Performance.....	66
6	Maintenance and Repair.....	67
6.1	Safety Notices.....	67
6.2	Staff Qualifications.....	67
6.3	Safety Inspection.....	67
6.4	Changing of the Fuse.....	68
6.5	Repair.....	69
6.6	Troubleshooting Table.....	70
7	Product Data.....	71
7.1	Technical Data.....	71
7.2	Ambient Conditions.....	71
7.3	Spare Parts and Accessories.....	72
8	Disposal.....	75

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0 Important Information about this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following product:

Item name: **FLEXILUX 200 LED**

Item number: **05.0740LED**

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and medical engineers who are entrusted with the operation, cleaning, and maintenance of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

0.4 Supplementary Documents

EMC The electromagnetic compatibility (EMC) requirements are described in the included separate brochure

Medical electrical equipment

pursuant to IEC 60601-1-2

The rules and information contained therein must be observed during installation and operation of this device.

1 General Information about the Device

1.1 Scope of Delivery

The scope of delivery for the device includes:

- one light source (item no. 05.0740LED)
- one power cord (item no. available on request)
- one set of instructions for use (TPA305-000-1208-20_H)
- one EMC brochure (TPI001-121-0803-20)

Check scope of delivery!

Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Product Description

1.2.1 Performance Characteristics and Function

Medical endoscopy	The present product is a light source for use in medical endoscopy.
380-700 nm	<p>The light source provides white light in the spectral range of 380-700 nm and is ready for use seconds after it is switched on.</p> <p>When using the FLEXISCOPE C1 HD / FLEXIVISION camera, the light source can be controlled with the buttons on the camera head.</p>
Standby function	The standby function allows the light to be handled even more flexibly and without delay. The light can be controlled either directly on the device or using one of the convenient, fully programmable camera head buttons.
LED with an adjustable luminous flux of up to 1400 lumens	The light is generated using the latest LED technology. The brightness can be adjusted continuously and in 5 % increments up to a maximum luminous flux of 1400 lumens. The long life of the LED means that the device is maintenance free.
Multi-light-guide adapter	The multi-light-guide adapter allows the light source to be connected with fiber optic light guides manufactured by Karl Storz, Richard Wolf and Olympus up to a diameter of 4.8 mm as well as proprietary light guides.
Automatic shut-off and power regulation	The multi-light-guide adapter has been designed to cause the LED to automatically switch to standby mode as soon as the light guide is detached. If the light output has been set to 80 % or above, it automatically drops to 60 % the next time the LED is switched on. This helps to preserve the LED, save energy, and prevent dazzle.
Overheating protection	If the temperature of the device starts to increase, which can occur due to a prolonged high ambient temperature, for example, a warning will appear on the display. The internal cooling system can compensate for a temperature increase for approximately one hour; after this time, the device will switch off automatically to prevent damage (see also section 2.2.5 Display of Faults).

1.2.2 Visual Overview

Front view


Figure 1-1: Front view of the light source.

- 1 ON / standby button
- 2 Multi-light-guide adapter
- 3 Standby LED button
- 4 Brightness controls


Rear view

Figure 1-2: Rear view of the light source.

- 5 Fuse holder
- 6 Connection for power cord
- 7 Main power switch
- 8 Potential equalization terminal
- 9 MIS-Bus ports
- 10 Port for service (covered)

1.2.3 Compatible Light Guides

The product is compatible with cold light guides and high-power light guides supplied by Karl Storz, Richard Wolf, Olympus and our own company:

- Fiber bundles of 3.5 to 4.8 mm Ø
- Length up to 300 cm

We recommend the use of high-power light guides.

1.3 Usage

1.3.1 Intended Use

This light source is intended for use in endoscopy systems for human medicine. In combination with a light guide, endoscope, and camera, its purpose is to illuminate the inside of the human body.

In combination with a suitable light guide the device meets the BF classification requirements for protection from electric shock according to IEC 60601-1 and is approved for use in combination with applied parts. This approval does not include use of the device in combination with applied parts intended for direct conductive contact with the heart.

Light guides that meet the requirements for applied parts classified as type BF feature the adjacent symbol.

Light guides that meet the requirements for applied parts classified as type CF offer an even higher level of protection. They feature the adjacent symbol.

1.3.2 Indications

There are no medical indications associated with the use of the light source if separated from the endoscopic system.

The indications for an endoscopic application depend on the medical condition of the patient and the individual risk/benefit analysis by the surgeon.

1.3.3 Contraindications

There are no medical contraindications associated with use of the light source itself separate from the endoscopic system.

Depending on the medical condition of the patient, contraindications may exist in respect to the patient's general condition or specific disease pattern.

The decision to perform an endoscopic procedure remains with the responsible surgeon and should be made on the basis of an individual risk/benefit analysis.

1.4 Conformity

1.4.1 Standards and Directives

The device meets the requirements of the following applicable standards:

- **IEC 60601-1** Medical electrical equipment - Part 1:
General requirements for basic safety and essential performance
- **IEC 60601-1-2** Medical electrical equipment - Part 1-2:
Electromagnetic compatibility

1.4.2 Medical Device Classification

The product is a Class I medical device within the meaning of applicable **EU guidelines**.

1.5 Marking

1.5.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the product or packaging.

Adhere to the instructions for use

CE marking

Item number

Serial number

Medical Device

Manufacturer

Date of manufacture

Caution (IEC 60601-1 3rd edition) /

Take note of accompanying documents (IEC 60601-1 2nd edition)

MR unsafe

Type BF applied part according to IEC 60601-1

Type CF applied part according to IEC 60601-1

Potential equalization

Release

ON / standby

Standby LED

Permissible storage and transport temperature

Permissible relative air humidity during storage and transport

Permissible atmospheric pressure during storage and transport

Caution: Federal (US) law restricts this device to sale by or on order of a physician

Separate collection for WEEE (waste of electrical and electronic equipment)

Fuse

Alternating current

Hot temperature

1.5.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Dangerous electrical voltage warning sign

Biohazard warning sign, risk of infection

1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

1.7 Obligation to Report Serious Incidents

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 General Safety Information
2.1 Safety Messages in this Document
2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

⚠ WARNING

Nature and origin of the injury risk

Potential consequences of non-observance

- > Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

⚠ CAUTION

Nature and origin of the injury risk

Potential consequences of non-observance

- > Preventive action

The following kind of message refers to a risk that could lead to material damage:

NOTICE

Nature and origin of the risk of material damage

Potential consequences of non-observance

- > Preventive action

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Notices

Risks despite highest of quality standards

Our products are developed and manufactured to the highest quality standards.

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

Keep the original packaging. Transport and store the product in its original packaging and use it to return goods if service support is required.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

If you detect any faults or malfunctions, inform us immediately.

WARNING! Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

WARNING! Operating the device in an open condition. Risk of electric shock. Do not open the device. Operate only when closed.

WARNING! Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

WARNING! Risks from the arrangement, setup, combination, or properties of connected or surrounding devices or equipment. Follow the instructions for use of the respective devices. Perform a risk assessment.

WARNING! High-intensity light source. Risk of injury to eyes. Do not look directly into the open end of the light guide.

WARNING! Risk of suffocation. Keep packaging material out of reach of children.

WARNING! Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

WARNING! Using the device during operation of a defibrillator. Risk to persons. Before discharging a defibrillator, remove device from the surgical field.

CAUTION! Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

WARNING! Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for installation, operation, and maintenance of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

2.2.3 Electromagnetic Compatibility

Medical electrical equipment is subject to stringent electromagnetic compatibility (EMC) requirements.

Despite the device's high interference immunity and low emitted interference, compliance with EMC-related requirements is necessary with regard to installation, installation location, and ambient conditions.

It is therefore important that you observe the related safety messages in this document and the EMC brochure included with the device.

According to IEC/CISPR 11, medical electrical equipment is divided into groups and classes:

Group 1: The device uses or generates RF-energy only for its internal functioning.

Group 2: The device generates RF-energy in the frequency range 9 kHz to 400 GHz in the form of electromagnetic radiation and through inductive or capacitive coupling for purposes of the treatment of material, its examination or analysis.

Class A: The device **is not suitable** for use in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

Note: The emissions characteristics of this device make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this device might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the device.

Class B: The device **is suitable** for operation in residential environments and such places that are connected immediately to a low voltage power supply network, that (also) provides for residential buildings.

WARNING! Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

WARNING! Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of the device specified by the manufacturer and covered by the scope of this document, including cables. Otherwise, degradation of the performance of the device could result.

2.2.4 Combination with Medical Electrical Equipment

The device can be combined with components from other manufacturers provided that all components are compliant with the medical electrical equipment safety requirements according to IEC 60601-1.

It is the operator's responsibility to check and make sure that the system is and remains fully operational.

When using devices from different manufacturers and when operating an endoscope and/or endoscopic accessories with electromedical devices, it must be ensured that the applied part is properly insulated (in accordance with the requirements for type BF/CF applied parts defined in IEC 60601-1).

Contact our Service Department if you have any questions regarding compatibility.

2.2.5 Display of Faults

No light guide or light guide attached incorrectly

If the light guide is not attached or is not attached correctly, you will see the following on the display:

Figure 2-1: Symbols displayed if light guide is not attached.

If this occurs, you should attach the light guide, attach it correctly, or switch the device off.

Operating error in LED module

In case of an operating error in the LED module, the display shows the following:

Figure 2-2: Message displayed on operating error in LED module.

If this occurs, you should replace the light source unit as quickly as possible.

Failure of housing fan

If the housing fan fails, the following message is displayed:


Figure 2-3: Message displayed when housing fan fails.

When housing fan fails, the temperature inside the device climbs in relation to the ambient temperature. If the ambient temperature is high, the device may overheat after 30 to 45 minutes.

If this occurs, you should replace the light source unit as quickly as possible.

Overheating

The device features an internal temperature sensor.

In very rare cases, for instance if the fan in the device is obstructed, the internal temperature of the device could rise above 60 °C.

If the temperature rises above 60 °C, you will see the following warning on the display:

Figure 2-4: Warning displayed when temperature exceeds 60 °C.

In order to protect the LED, the device will, in the event of overheating, limit the maximum light output to 70 %. The device will not switch off, but is not guaranteed to remain fully functional.

If this occurs, you should replace the light source unit as quickly as possible.

3 Installation and Initial Operation

3.1 Safety Notices

WARNING

Electrical connections installed improperly

Risk of fire, short circuit, or electric shock

- > Make sure that the electrical connections are installed in accordance with the relevant national technical regulations

WARNING

Medical electrical system installation

Risk of fire, short circuit, or electric shock

- > Installation must only be performed by qualified personnel
- > If a number of electrical devices are being used in combination, observe IEC 60601-1
- > Non-medical devices that meet the applicable IEC safety standards can be connected only via a medical isolating transformer
- > Do not connect any additional non-medical devices to a medical electrical system
- > Signal lines from devices with a functional ground (FE) connection, which are connected to different branches of the mains power supply, must be galvanically isolated at both ends
- > Connect the devices only to a power supply with a protective ground (PE)
- > After installing a medical electrical system, it must be inspected according to IEC 62353

**⚠ WARNING****Use of multiple sockets**

Risk of fire, short circuit, electrical shock, reduced level of safety

- > Whenever possible, avoid the use of multiple sockets
- > If required, use medically approved multiple sockets
- > Never connect multiple sockets in series
- > Do not cover multiple sockets (leads to heat accumulation)
- > Do not place multiple sockets on the floor
- > Use traction relief
- > Devices should be connected to a shared multiple socket only if they are all intended to be used as part of the medical electrical system

**⚠ WARNING****Electromagnetic interactions**

Possible impaired function and image quality

- > Observe the EMC brochure supplied
- > Make sure that neighboring devices are not causing reciprocal electromagnetic interference
- > Perform a post-installation function check

3.2 Staff Qualifications

Qualified personnel Personnel responsible for setting up or installing medical electrical equipment must possess the necessary qualifications and be familiar with the safety regulations in place at the installation location and the applicable official regulations.

3.3 Installation

3.3.1 Setup

Setting up the device

Observe the safety notices provided at the start of this chapter and refer to the enclosed EMC brochure.

Proceed as follows:

WARNING! Risk from installation in potentially explosive locations. Increased risk of fire and explosion in oxygen-enriched atmospheres. Device must be installed outside of potentially explosive locations and not in the vicinity of flammable substances.

WARNING! Risk resulting from electromagnetism. Possible malfunctions and image interference. Use only the cables supplied. Ensure that all devices operated in the vicinity meet the EMC requirements. Perform a post-installation function check.

WARNING! Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other devices should be observed to verify that they are operating normally.

Maintain distance!

1. Install the device.

Maintain minimum distances to other devices:

- at least 5 cm to the side,
- at least 6 cm from the back.

Install safely!

When installing the device, ensure that:

- it is positioned on an level, non-slip surface,
- it is located outside of the sterile area on a sufficiently stable support,
- it is protected from water drips and splashes,
- it will not be subject to vibration during operation,
- the ventilation slots are not obstructed,
- the main switch on the rear of the device is accessible.

2. Connect the **potential equalization terminal** on the rear of the device to the potential equalization terminal at the installation site or on the equipment cart. Use a potential equalization conductor in accordance with DIN 42801 with at least a 4 mm² cross-section.

3. If you would like to control the light source using the buttons on the camera head, connect the camera to the light source at one of the ports labeled **MIS-BUS**.

Observe the instructions for use supplied with the camera.

- The device is properly installed.

3.3.2 Connection to Mains Power

Connecting the device to the mains

Do not connect the device to the mains until you have completed all the previous activities. Observe the safety notices provided at the start of this chapter.

Proceed as follows:

1. Check that the supply voltage at the site of operation corresponds to that specified on the device nameplate.
2. Insert the power cord into the mains connection socket on the rear of the device.

WARNING! Sparks in a potentially explosive location. Risk of explosion. Ensure that the device is not situated in a potentially explosive location when connecting to the power supply.

3. Connect the power cord to the power supply.

Route cables safely!

4. When routing the cable, make sure it does not present a trip or entanglement hazard.

- The device is properly connected to the power.

3.4 Initial Operation

3.4.1 Initial Operation and Function Check

Performing a function check

When operating the device for the first time, check the installation and make sure the device is in full working order.

Precondition: Installation must have been completed.

Proceed as follows:

1. Switch ON at the main switch on the rear of the device.
 - > *The device switches to standby mode and the ON / standby button on the front of the device is illuminated green.*
2. Switch on the device.
 - > *After approximately 10 seconds, the device is ready for operation.*
3. Connect a light guide as described in [section 4.5](#).

WARNING! High-intensity light. Risk of injury to eyes. Do not look directly into the open end of the light guide.

4. Check that you are able to adjust the light intensity as described in [section 4.7](#).
5. Check that the LED can go into standby mode and be reactivated as described in [section 4.8](#).

- This concludes the function check.

4 Operation

4.1 Safety Notices

WARNING

Modifications to the installation

Risk of fire, short circuit, and electric shock

- > Follow the safety notices in the Installation and Initial Operation chapter
- > Do not make unauthorized changes to the installation
- > If the installation contains a multiple socket, do not connect additional devices to this without authorization
- > Never connect multiple sockets in series

WARNING

Current leakage when touching the patient

Risk of electric shock

- > Never touch the patient and the device housing or accessible contacts of electric devices at the same time
- > Never touch the patient and non-medical electric devices at the same time

WARNING

Positioning the equipment cart inside the sterile area

Risk of infection

- > Do not position the equipment cart inside the sterile area
- > Measure the cable lengths so there is enough room for movement without compromising sterility
- > The light guide must be anchored securely within the direct OR environment to avoid any risk of accidental slipping

⚠ WARNING

Interactions between devices in simultaneous use (e.g. lasers, electro-surgery)

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the required type BF or type CF classification requirements according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active RF electrodes
- > Do not activate RF electrodes in the vicinity of flammable gases or liquids
- > Before using RF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

4.2 Staff Qualifications

This device may be operated only by physicians and medical assistants with the requisite training, advanced training, knowledge, and practical experience in the endoscopic discipline concerned as defined by the provisions in place at the site of operation.

4.3 Technical Inspection Prior to Use

4.3.1 Visual Inspection

Performing a visual inspection

*A visual inspection must be performed **prior to each medical procedure.***

Proceed as follows:

- | | |
|------------------------------------|--|
| Is the housing undamaged? | 1. Inspect the housing of the light source for external damage.
Do not use the light source if there is any external damage to the housing. |
| Are the electric cables undamaged? | 2. Check that all electric cables are undamaged and positioned safely.
Do not use the device if the electric cables are damaged or present a trip or entanglement hazard. |
| Is the light guide undamaged? | 3. Check that the light guide you intend to use is undamaged.
Do not use the light guide if it is broken or twisted. |
- This concludes the visual inspection.

4.3.2 Function Check

Performing a function check

*A function check must be performed **prior to each medical procedure.***

Proceed as follows:

1. Check that the light guide you intend to use meets the applicable electrical safety requirements for the procedure to be performed.
Do not use the light guide if its electrical insulation does not correspond to that of the applied part being used. Depending on the procedure being performed, the light guide must be compatible either with an applied part classified as type BF or one classified as type CF according to IEC 60601-1.
 2. Connect the light guide.
 3. The light guide should be directed away from the eyes to avoid dazzling anybody.
 4. Switch on the light source.
Check that you are able to adjust the brightness and that the system is functioning correctly.
Do not use the light source if you are unable to adjust the brightness.
- This concludes the function check.

4.4 Switching On and Off, Separation from the Power Supply

Switching on the device

Proceed as follows:

1. Switch on the device at the front.

Device does not switch on?

The main switch on the rear is switched off.

Switch ON at the main switch.

- > After approximately 10 seconds, the device is ready for operation.
- > The ON / standby button on the front is illuminated white.

- The device is on.

Switching off the device

Proceed as follows:

1. Switch off the device at the front.
 - > The device switches to standby mode.
 - > The ON / standby button on the front is illuminated green.

When the device switches back on, the light output depends on the setting prior to the switch to standby mode:

If the light output was 85 % or higher when in standby mode, it will be reduced to 80 % when switched back on. If lower than 85 %, the light output will be the same as prior to the switch to standby mode.

2. To switch off the device completely, switch off at the main switch on the rear.

- The device is off.

Separating the device from power

Proceed as follows:

1. To separate the device completely from its power supply, disconnect the power cord. The appliance coupler is considered as disconnecting device.

Make sure that the rear of the device is accessible at any time.

- The device is completely separated from the power supply.

4.5 Connection of a Light Guide

Connecting the light guide

The multi-light-guide connector enables the connection of light guides supplied by different manufacturers.

Original light guides from Karl Storz, Richard Wolf and Olympus can be attached directly without adapters. Proprietary light guides must be used with the corresponding light guide adapters.

Proceed as follows:

WARNING! Use of an unsuitable light guide. Risk of burns caused by insufficient electrical insulation. Use a light guide that meets at minimum the type BF classification requirements according to IEC 60601-1.

WARNING! High-intensity light. Risk of injury to eyes. Do not look directly into the open end of the light guide.

1. Slide the light guide into the multi-light-guide adapter until it clicks into place.
 - > *The light exits the distal end of the light guide.*

The light guide is connected.

4.6 Decoupling the Light Guide

Decoupling the light guide

Proceed as follows:

1. Rotate the multi-light-guide adapter clockwise as far as it will go.
 - > *The light guide detaches.*
 - > *The LED switches off, the device remains switched on.*

2. Pull out the light guide.
3. Detach the adapters used, and arrange for reprocessing the light guide. Observe the instructions for use supplied with the light guide.

The light guide has been decoupled.

4.7 Setting the Brightness

Setting the brightness

The brightness is continuously adjustable in 5 % increments.

Proceed as follows:

1. Press the arrow buttons on the device to set the brightness.

Tap briefly Tap an arrow button briefly to set the brightness in 5 % increments.
 > *The device emits a brief sound with each tap.*

Press and hold down Hold down an arrow button to continually increase or decrease the brightness.

The brightness is now set.

4.8 Using the LED Standby Function

Using the LED standby function

Use the LED standby function if you want to briefly dim the light to the minimum level.

Proceed as follows:

1. To set the light supply at 5 %, tap the LED standby button.

> The LED output is reduced to 5 %; the device remains switched on.

2. To set the light supply to its previous value, tap the LED standby button once more.

> The LED switches back on immediately and supplies light at its previous setting.

The LED standby function has been used.

4.9 Storage and Transport

To prevent the device from falling when being transported on the equipment cart, secure with a tension belt.


Figure 4-1: Device on equipment cart, secured with tension belt.

Store the device in a dry, well-ventilated, temperature-controlled room, where it is protected from dust.

Switch off the device prior to storage and remove the power cord and accessories.

When placing into storage or storing temporarily, make sure that the device cannot be damaged by its surroundings.

To avoid damage, keep out of direct sunlight and away from radioactivity and strong electromagnetic radiation.

The ambient conditions for transport, storage, and operation are described in [section 7.2](#).

5 Cleaning

5.1 Safety Notices

⚠ WARNING

Risks associated with improper cleaning

Risk of fire, short circuit, and electric shock

- > Before cleaning, switch off the device at its rear main switch and disconnect from the mains power supply
- > Make sure no moisture penetrates the device
- > The device should not be reconnected to the mains power supply until it has completely dried

NOTICE

Improper cleaning

Risk of damage to the device

- > Do not use sponges or wipes that are dripping wet
- > Do not use abrasive cleaners or solvents

5.2 Performance

Cleaning the device

Proceed as follows:

Separate the device from power

1. Switch OFF the device at the main switch on the rear and pull out the power cord from the rear of the device to disconnect the power.
2. Clean the device.

Use slightly moistened wipe only

Wipe the external surfaces with a soft wipe. Only slightly moisten the wipe with a mild soap solution or isopropyl alcohol.

Stubborn stains can be removed with a mild, ammonia-based detergent.

Do not reconnect the device to the power until all the cleaned parts are completely dry.

3. Connect the device to the power supply and switch back on at the main switch.

This concludes the cleaning process.

6 Maintenance and Repair

6.1 Safety Notices

⚠ WARNING

Improper maintenance and repair

Risk of fire, short circuit, and electric shock

- > All maintenance and repair work must only be performed by qualified personnel
- > Do not make any modifications to the device

6.2 Staff Qualifications

Qualified personnel Personnel responsible for maintaining the device or carrying out safety inspections must possess the necessary qualifications and be familiar with the local safety regulations in place at the installation location and the applicable official regulations.

Repairs must be performed by our service personnel only.

6.3 Safety Inspection

We recommend to arrange for safety inspections to be performed every 24 months in accordance with Section 5 of IEC 62353.

The inspection comprises the following steps:

- visual inspection
- measuring the protective resistance (only for protection class I devices)
- measuring the leakage currents
- measuring the insulation resistance

Document results The results of the inspection should be documented in an inspection report in accordance with Section 6 and Annex G of IEC 62353.

We recommend that all components of a medical electrical system be subjected to documented technical safety inspections at regular intervals.

6.4 Changing of the Fuse

Changing the fuse

Proceed as follows:

1. Switch OFF the device at the main switch on the rear and pull out the power cord from the rear of the device to disconnect the power.

NOTICE! Pulling the power cable will damage it. Disconnect the power at the plug.

2. Release the fuse holder with a screwdriver.


Figure 6-1: Releasing the fuse holder with a screwdriver.

3. Pull out the fuse holder carefully, without tilting.

WARNING! Unsuitable electrical fuse. Risk of fire and short circuit. Use the specified fuse only and never bypass the fuse.

4. Replace the faulty fuse.
 5. Carefully slide the fuse holder back in until it clicks into place.
- This concludes the fuse replacement.

6.5 Repair

Should you need to arrange a repair for the device, please contact one of our subsidiaries. Contact details can be found on the back of these instructions for use.

When you send in equipment, please enclose as accurate a **fault description** as possible, and record the item number and serial number of the product on the delivery note. You will find these details on the nameplate.

WARNING

Contaminated device

Risk of infection

- > Reprocess the device properly prior to shipping (*chapter 5*)
- > Clearly identify contaminated products as such

For shipping, please use the original packaging whenever possible. Affix a label to the outer packaging indicating the disinfection status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

6.6 Troubleshooting Table

Issue	Possible causes	Fixes
No operation	No power supply	Check the mains connection and connect the power if necessary Check the fuse and replace if necessary
	Malfunction in power supply unit	Send the light source in for repair
Light doesn't switch on	Device is overheated	Allow device to cool down
	Shutter closed/light off	Open the shutter/switch on light
The light guide cannot be coupled	Adapter on light guide connector missing	Screw on adapter
Not enough light	Light guide not connected correctly	Check that the light guide is positioned properly; correct if necessary
	Light guide defective	Replace the light guide
	Adapter not screwed on correctly	Screw on the adapter
Fan doesn't work	Malfunction in fan	Send the light source in for repair
Brightness significantly too low	Light guide not completely inserted	Use the correct light guide connection
Light source cannot be controlled	Interface cable not correctly plugged in or defective	Check connection or replace cable

Table 6-1: Troubleshooting table.

7 Product Data

7.1 Technical Data

Dimensions (W x H x D)	295 x 130 x 355 mm
Weight	5.8 kg
Current consumption	1.0 - 0.42A
Supply voltage	100-240V~, 50/60Hz
Protection class acc. to IEC 60601-1	Protection class I
Device fuse	T1.6AL 250V
Interfaces	2x MIS-Bus
Degree of protection	IP X0
Applied part as per IEC 60601-1	BF
Classification as per IEC/CISPR 11	Group 1 Class B

7.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

Operating conditions

Temperature	+10 °C to +40 °C
Relative air humidity	10 % to 90 %
Atmospheric pressure	70 kPa to 106 kPa

CAUTION

Non-observance of the ambient conditions

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

> Observe operating conditions as well as transport and storage conditions

7.3 Spare Parts and Accessories

Use original spare parts and accessories only.

Image	Designation	Item number
	Micro-fuse T1.6AL 250V, 5x20 mm	94600.
	MIS-Bus cable, 2.25 m	A057635
	MIS-Bus cable, 0.75 m	A059584
	Power cord (country-specific)	on request

Table 7-1: Spare parts and accessories.

Image	Designation	Item number
	High-power light guide, 3.5 x 1,800 mm (standard light guide)	05.0084I
	High-power light guide, 3.5 x 2,300 mm (standard light guide)	05.0088I
	High-power light guide, 3.5 x 3,000 mm (standard light guide)	05.0085I
	High-power light guide, 4.8 x 1,800 mm (standard light guide)	05.0086I
	High-power light guide, 4.8 x 2,300 mm (standard light guide)	05.0090I
	High-power light guide, 4.8 x 3,000 mm (standard light guide)	05.0087I
	High-power light guide, 3.5 x 2,300 mm (high-temperature light guide)	05.0088I.ht
	High-power light guide, 3.5 x 3,000 mm (high-temperature light guide)	05.0085I.ht
	High-power light guide, 4.8 x 2,300 mm (high-temperature light guide)	05.0090I.ht

Table 7-2: Spare parts and accessories.


Image	Designation	Item number
Light guide adapter (to connect to the light guide, light source side) for:		
	Storz systems	05.0100z
	Wolf systems	05.0102b
	Olympus systems	05.0101o
Light guide adapter (to connect to the light guide, endoscope side) for:		
	Storz endoscopes	05.0108z
	Wolf endoscopes	05.0110b
	ACMI endoscopes	05.0112c
Light guide adapter (to connect to the endoscope) for:		
	Storz light guides	05.0114z
	Wolf light guides	05.0116b

Table 7-3: Spare parts and accessories.

8 Disposal

WARNING

Contaminated device

Risk of infection

> The device must be reprocessed prior to disposal (*chapter 5*)

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

The product is packed in a polyethylene (PE) bag and a PE foam inlay (0.50 kg), as well as a corrugated cardboard box (1.18 kg). Other packing materials used are a cardboard insert (0.22 kg) and a membrane pad made of PE film and corrugated cardboard (0.30 kg).

Any product carrying this symbol must be disposed of separately through dedicated electrical and electronic devices recycling. Within the EU, such disposal is taken care of by the manufacturer free of charge.

SCHÖLLY FIBEROPTIC GMBH
Robert-Bosch-Str. 1-3
79211 Denzlingen
Germany