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## **Instructions For Use**

Arthroscopes  
Endoscopes for Orthopedic Procedures  
(not valid for USA)

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| <b>Table of Contents</b>                                     |  | <b>Page</b> |
|--|--|-------------|
| <b>0 Important Information about this Document</b>           |  | <b>29</b>   |
| 0.1 Scope of Validity, Identification, Purpose               |  | 29          |
| 0.2 Target Group   |  | 29          |
| 0.3 Using and Storing this Document                          |  | 29          |
| <b>1 General Information about the Device</b>                |  | <b>30</b>   |
| 1.1 Scope of Delivery  |  | 30          |
| 1.2 Product Description                                      |  | 30          |
| 1.2.1 Performance Characteristics and Function               |  | 30          |
| 1.2.2 Visual Overview  |  | 30          |
| 1.2.3 Required Accessories                                   |  | 31          |
| 1.2.4 Compatible Light Guides                                |  | 31          |
| 1.3 Usage  |  | 31          |
| 1.3.1 Intended Use   |  | 31          |
| 1.3.2 Indications  |  | 31          |
| 1.3.3 Contraindications                                      |  | 31          |
| 1.4 Conformity   |  | 32          |
| 1.4.1 Standards and Directives                               |  | 32          |
| 1.4.2 Medical Device Classification                          |  | 32          |
| 1.5 Marking  |  | 32          |
| 1.6 Service Department Contact Details                       |  | 32          |
| 1.7 Obligation to Report Serious Incidents                   |  | 32          |
| <b>2 General Safety Information</b>                          |  | <b>33</b>   |
| 2.1 Safety Messages in this Document                         |  | 33          |
| 2.1.1 Safety Messages at the Start of a Chapter              |  | 33          |
| 2.1.2 Safety Messages in the Body of the Text                |  | 33          |
| 2.2 Product Safety   |  | 34          |
| 2.2.1 Basic Safety Notices                                   |  | 34          |
| 2.2.2 Staff Qualifications                                   |  | 34          |
| <b>3 Operation</b>   |  | <b>35</b>   |
| 3.1 Safety Notices   |  | 35          |
| 3.2 Staff Qualifications                                     |  | 36          |
| 3.3 Technical Inspection Prior to Use                        |  | 37          |
| 3.3.1 Performing a Visual Inspection                         |  | 37          |
| 3.3.2 Performing a Function Check                            |  | 37          |
| 3.4 Connecting a Light Guide                                 |  | 38          |
| 3.5 Stopping Usage and completing Pre-Cleaning               |  | 38          |
| <b>4 Reprocessing</b>  |  | <b>39</b>   |
| 4.1 Safety Notices   |  | 39          |
| 4.2 Staff Qualifications                                     |  | 40          |
| 4.3 Validated Procedures                                     |  | 40          |
| 4.4 The Process Flow   |  | 41          |
| 4.5 Cleaning and Disinfecting Agents                         |  | 41          |
| 4.6 Manual Cleaning and Disinfection                         |  | 42          |
| 4.6.1 Overview   |  | 42          |
| 4.6.2 Performing Manual Cleaning and Disinfection            |  | 42          |
| 4.7 Automated Cleaning and Thermal Disinfection              |  | 43          |
| 4.7.1 Overview   |  | 43          |
| 4.7.2 Performing Automated Cleaning and Thermal Disinfection |  | 44          |
| 4.8 Sterilization  |  | 45          |
| 4.8.1 Performing Steam Sterilization                         |  | 45          |
| 4.8.2 STERRAD® Sterilization Process 100S, NX, 100NX         |  | 46          |
| <b>5 Maintenance and Repairs</b>                             |  | <b>47</b>   |
| 5.1 Troubleshooting  |  | 47          |
| 5.2 Repair   |  | 48          |
| <b>6 Product Data</b>  |  | <b>49</b>   |
| 6.1 Technical Data   |  | 49          |
| 6.2 Ambient Conditions                                       |  | 50          |
| 6.3 Spare Parts and Accessories                              |  | 50          |
| <b>7 Disposal</b>  |  | <b>51</b>   |

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## 0 Important Information about this Document

### 0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following products:

**Item numbers**

21.0003a

21.0008a

21.0011a

21.0101a

21.0617a

21.0619a

21.0621a

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

These instructions for use do not describe endoscopic procedures or techniques during surgery.

### 0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and employees of central sterile supply who are entrusted with the operation, handling, and reprocessing of the device.

### 0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

## 1 General Information about the Device

### 1.1 Scope of Delivery

The scope of delivery for the device includes:

- 1x endoscope, comprising  
(pre-installed on the endoscope):
  - endoscope light guide adapter, Wolf systems
  - endoscope light guide adapter, Storz systems
- 1x instructions for use

Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

### 1.2 Product Description

#### 1.2.1 Performance Characteristics and Function

Rigid endoscopes permit the visualization of the inside of the body and body cavities.

Depending on the endoscopic discipline and anatomical region, the approach to the body cavity may be through a natural orifice or a surgically created opening.

Rigid endoscopes consist of a fiber optic cable and sensitive image transmission system with eyepiece.

The fiber optic cable is used to illuminate the site inside the body. The connector for connecting the light guide to the light source is situated at the proximal end of the endoscope. The endoscope light guide adapters required to connect light guides are included in the scope of delivery.

The distal end of the endoscope features an objective lens that captures the image from inside the body. The image is sent through the image transmission system to the eyepiece.

The eyepiece cap connects the eyepiece to an endocoupler, which transmits the image to a camera. The camera's controller converts the signal so that the image can be displayed on a monitor.

Special instruments are required for endoscopic procedures.

It is the surgeon's responsibility to make sure the instruments required for a particular endoscopic procedure are assembled and available for use.

#### 1.2.2 Visual Overview

An arthroscope design illustration is provided below.


Figure 1-1: Arthroscope illustration.

- ① Objective
- ② Eyepiece
- ③ Light guide connector (ACMI)
- ④ Endoscope light guide adapter, Wolf systems
- ⑤ Endoscope light guide adapter, Storz systems

### 1.2.3 Required Accessories

Additional instruments are required for endoscopic procedures. It is the surgeon's responsibility to make sure the instruments required are assembled and available for use.

There is no guarantee that instruments selected solely using maximum insertion portion width and working length will be compatible in combination.

Depending on the arthroscope used ([section 6.1](#)) standard trocars with a diameter suitable for 2,7 mm or 4 mm endoscopes are required for the application.

### 1.2.4 Compatible Light Guides

The specifications for compatible light guides are:

- Fiber bundles of 3.5 mm Ø
- Length up to 300 cm

The endoscope light guide adapters included in the scope of delivery ([section 1.1](#)) allow for the connection of light guides from different manufacturers.

## 1.3 Usage

### 1.3.1 Intended Use

Rigid endoscopes without a working channel permit visualization of the inside of the body and body cavities.

Arthroscopes are intended for use in minimally invasive orthopedic procedures. They are used to visualize and illuminate particular regions of the body during diagnostic and therapeutic procedures. They are used in particular for procedures on the knee, shoulder, jaw, wrist, ankle and elbow, for example.

The rigid endoscopes described herein are intended to be used in combination with medical electrical devices that meet at **minimum** the type **BF** classification requirements according to IEC 60601-1 for elevated protection from electric shock.

This insulation barrier can be realized by the medical electrical devices themselves or by the connection cables to the endoscope and must be assured for every connection between the endoscope and connected devices.

Devices or connection cables that meet the type BF classification requirements for applied parts feature the adjacent symbol.

Devices and connection cables that meet the type CF classification requirements for applied parts offer an even higher level of protection. They bear the adjacent heart symbol and are approved for direct conductive contact with the heart.

Defibrillation-proof devices or connection cables are also compatible. These will be marked with one of the adjacent symbols depending on the degree of protection against electric shock offered.

### 1.3.2 Indications

The indications for an endoscopic application depend on the medical condition of the patient and the individual risk/benefit analysis by the surgeon.

### 1.3.3 Contraindications

Contraindications may exist depending on the patient's general condition or specific symptoms.

The decision to perform an endoscopic procedure remains with the responsible surgeon and must be made on the basis of an individual risk/benefit analysis.

## 1.4 Conformity

### 1.4.1 Standards and Directives

The device meets the requirements of the following applicable standard:

- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:  
Particular requirements for basic safety and essential performance of endoscopic equipment

### 1.4.2 Medical Device Classification

The product is a Class IIa medical device within the meaning of applicable **EU guidelines**.

## 1.5 Marking

Observe the symbols used on product and packaging.

In addition to the international standardized symbols, we use the following symbols:

 The device is suitable for autoclave sterilization

 Fully refurbished device ("refurbished device")

 Federal (US) law restricts this device to sale by or on order of a physician

Find a list of all symbols used on the product, packaging and in documents online:

<https://ifu.schoelly.de>

## 1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

## 1.7 Obligation to Report Serious Incidents

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## 2 General Safety Information

### 2.1 Safety Messages in this Document

#### 2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

#### **WARNING**

##### **Nature and origin of the injury risk**

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

#### **CAUTION**

##### **Nature and origin of the injury risk**

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

#### **NOTICE**

##### **Nature and origin of the risk of material damage**

Potential consequences of non-observance

> Preventive action

#### 2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

**WARNING!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

**CAUTION!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

**NOTICE!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

## 2.2 Product Safety

### 2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

The original packaging should be retained and reused when returning the device for servicing.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

**WARNING!** This device is supplied in a non-sterile condition. Risk of infection. Carry out reprocessing of the device and its accessories prior to first use.

**WARNING!** Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

**WARNING!** Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

**WARNING!** High-intensity light source. Risk of injury to eyes. Do not look directly into the open end of the light guide or light exit of the endoscope.

**WARNING!** Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

**WARNING!** Using the device during operation of a defibrillator. Risk to persons. Before discharging a defibrillator, remove device from the surgical field.

**CAUTION!** Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

**WARNING!** Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

**WARNING!** When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. This is particularly important if a type CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

Note: The product may consist of metal alloys containing cobalt. Refer to the corresponding pictogram on the packaging label.

### 2.2.2 Staff Qualifications

Specific qualifications are required for operation and reprocessing of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

## 3 Operation

### 3.1 Safety Notices

#### **WARNING**

##### **Interactions between devices in simultaneous use (e.g. lasers, electrosurgery)**

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the requirements classified as type BF according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active HF electrodes
- > Do not activate HF electrodes in the vicinity of flammable gases or liquids
- > Before using HF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

#### **WARNING**

##### **Displaying a recorded image instead of the live image or an altered image orientation**

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera
- > Ensure that the live image is displayed in the correct orientation (not mirrored)

#### **WARNING**

##### **Use of non-sterile parts**

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Always carry out a visual inspection prior to use
- > Anchor the light guide in the surgical field and secure against slipping, allowing some slack for maneuvering

#### **WARNING**

##### **Recontamination resulting from improper handling**

Infection risk for the patient

- > Adhere to the infection control practices that are in place

**CAUTION**

**High temperatures resulting from light source usage**

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage, temperature on the distal tip > 41 °C

- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest possible light output setting that will allow you to illuminate the target area
- > Do not touch the light source in the vicinity of the lamp
- > Do not touch light guide connections
- > Do not allow the distal end of the endoscope to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Do not exceed the maximum permitted operating conditions for the endoscope

**CAUTION**

**Coupling laser beams in the endoscope**

Damage to the eyes caused by looking directly into the eyepiece

- > Wear laser safety glasses

**NOTICE**

**Improper handling**

Risk of damage to the product

- > Set down carefully without bending or knocking
- > Avoid lever forces
- > Avoid scratching the surface
- > Take care when removing instruments from the surgical field

**3.2 Staff Qualifications**

The device may be operated only by physicians and medical assistants who have received instruction in how to use the device and who have the requisite training or advanced training, knowledge, and practical experience in the endoscopic procedure to be applied as defined by the provisions in place at the site of operation.

### 3.3 Technical Inspection Prior to Use

#### 3.3.1 Performing a Visual Inspection

*A visual inspection must be performed **prior to each use**.*

1. Check the endoscope and all the components to be used for external damage.

Do not use the product if it has any sharp edges or corners, bulges, or rough surfaces that might cause injury to the patient.

2. Check that the fiber optics in the endoscope are in full working order.

To do this, hold the distal end toward a bright light (not a cold light source) and the light guide connection toward your eyes. Gently move the endoscope from side to side, taking note of the brightness of the fibers.

If more than 20 % of the fibers remain dark, the device must not be used. Submit the endoscope to repair.

3. Check that the glass surfaces and fiber optic end faces of the endoscope are clean and smooth.

Do not use the endoscope if the surfaces are soiled or scratched. Set aside the endoscope and follow the steps described in [section 5.1](#).

4. Make sure you can see a clear, bright, complete image.

Look through the eyepiece and assess the image quality.

Do not use the endoscope if the image is yellowed, dark, speckled, or cropped. Set aside the endoscope and follow the steps described in [section 5.1](#).

5. Check that all the equipment to be used for the procedure is mutually compatible and can be interlocked, if applicable.

Only use the OR equipment if the isolation barrier for protecting against electric shock meets **at minimum** the type **BF** classification requirements according to IEC 60601-1 (see [section 1.3.1](#)).

#### 3.3.2 Performing a Function Check

*A function check must be performed **prior to each medical procedure**.*

1. Connect a camera.
2. Check that camera head and endoscope are firmly seated and securely connected to each other.
3. Switch on all the system components that you intend using for the procedure.

4. Direct the camera head at a nearby object in the room and focus the image.

Check that you are able to achieve a sharp, bright, good quality image.

Do not use the system if the image contains vertical lines or color variations, flickers, or you are unable to achieve a sharp, bright, good quality image.

### 3.4 Connecting a Light Guide

*Different light guides require different light guide adapters to connect to the endoscope.*


Figure 3-1: Mounting endoscope light guide adapters.

1. Screw the suitable light guide adapter to the light guide connector on the endoscope and, if applicable, to the light guide.
2. Connect the light guide to the light guide connector.

### 3.5 Stopping Usage and completing Pre-Cleaning

*Complete the pre-cleaning process **immediately after use.***

1. Remove and discard the sterile drape, if present.
2. Uncouple the endoscope from the endocoupler.
3. Detach the light guide and all detachable parts (such as light guide adapters).

**NOTICE!** The eyepiece cap is non-detachable.

4. Pre-clean the product at the conclusion of the surgery by wiping it with a lint-free and wetted, but not dripping wipe until no more residue can be seen.
5. Dry the product using a soft, lint-free wipe.
6. Arrange for reprocessing.

Ensure that the product is reprocessed within **6 hours.**

## 4 Reprocessing

### 4.1 Safety Notices

#### **WARNING**

**If the product is potentially contaminated with Creutzfeldt-Jakob pathogens, it is not suitable for reprocessing**

Risk of cross-contamination during use and reprocessing

- > Do not reprocess products that are potentially contaminated
- > Products that are potentially contaminated must be disposed of

#### **CAUTION**

**Improper cleaning and disinfection**

Risk to patient as a result of inadequate cleaning and disinfection, damage to the device

- > Use an automatic washer/disinfector that meets the performance requirements of ISO 15883-1
- > Reprocess all detachable parts (such as light guide adapters)
- > Load the device with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Connect instruments with lumens and channels directly to the dedicated rinsing ports on the device
- > Disassemble all valves on instruments
- > Ensure that the washer/disinfector is properly maintained
- > Only use cleaning and disinfecting agents that are approved for the product
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing cleaning and disinfecting agents (active ingredients: Aldehyde, alcohol)

#### **NOTICE**

**Contact with chloride-based solutions**

Corrosion and destruction of the device

- > Avoid contact with chloride-based solutions, which can be found, for example, in surgical residue, lotions, drugs, saline solutions, and cleaning and disinfecting agents
- > After contact with chloride-based solutions, rinse products sufficiently with deionized water and dry completely

#### **NOTICE**

**Use of ultrasonic bath**

Damage to the endoscope

- > Do not expose endoscope to ultrasonic bath

### **NOTICE**

#### **Unsuitable cleaning agents and process chemicals**

Corrosion damage, premature aging, and visible material changes

- > Only use cleaning and disinfecting agents that are approved for the product
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials

## **4.2 Staff Qualifications**

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must always fall under the responsibility of qualified personnel who have the necessary knowledge and expertise.

This knowledge and expertise can be acquired by completing further training in this field or as a result of completing a dedicated qualification with practical experience, supplemented by appropriate further training measures, where required.

## **4.3 Validated Procedures**

The effectiveness of the following procedures:

- Manual cleaning and disinfection
- Automated cleaning and thermal disinfection
- Steam sterilization
- Low temperature plasma sterilization (STERRAD® 100S, NX, 100NX)

as described in this document, has been fully validated.

It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

#### 4.4 The Process Flow

The product has to be prepared for reprocessing immediately after use by pre-cleaning.

Refer to [section 3.5](#) for pre-cleaning instructions.

The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Cleaning and disinfection (manual or automated)
- Sterilization

During reprocessing, staff should wear personal protective clothing.

The product must be cleaned thoroughly at the beginning of reprocessing. It is essential that the sterilization medium reaches all parts of the product.

The optimum and most reliable reprocessing results are achieved using automated cleaning and disinfection with subsequent steam sterilization using the pre-vacuum steam procedure.

Adhere to national statutory regulations, national and international standards and directives, and the infection control practices that are in place at your institution for reprocessing.

New and returned products following repair have to undergo complete reprocessing before first use.

The alternating use of different reprocessing methods can cause premature aging of the product.

Please see [www.a-k-i.org](http://www.a-k-i.org) for more detailed information about hygienically safe, value-preserving reprocessing methods that also protect the materials.

#### 4.5 Cleaning and Disinfecting Agents

The **manual cleaning** process was validated using:

- Cidezyme®/Enzol® (Johnson & Johnson)

The **manual disinfection** process was validated using:

- Cidex® OPA (Johnson & Johnson)

The **automated cleaning** process was validated using:

- needisher® MediClean forte 0.5 %  
(Chemische Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning and disinfecting agents only.

Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

## 4.6 Manual Cleaning and Disinfection

### 4.6.1 Overview

| Stage | Work step            | Temperature (°C/°F)     | Time (min) | Water quality  | Cleaning solution                 |
|-------|----------------------|-------------------------|------------|----------------|-----------------------------------|
| I     | Cleaning             | 35-45 / 95-113          | 5          | Drinking water | Enzymatic cleaning solution       |
| II    | Intermediate rinsing | Room temperature (cold) | 3 x 1      | Drinking water | ---                               |
| III   | Disinfection         | 20-25 / 68-77           | 12         | ---            | Disinfectant solution (undiluted) |
| IV    | Final rinsing        | Room temperature (cold) | 3 x 2      | DI water*      | ---                               |
| V     | Drying               | Room temperature        | ---        | ---            | ---                               |

Table 4-1: Overview of manual cleaning and disinfection.

\* DI water = Deionized water (demineralized, low-germ, max. 10 germs/ml, as well as low endotoxin contaminated, max. 0.25 endotoxin units/ml)

### 4.6.2 Performing Manual Cleaning and Disinfection

#### Complete stage I: Cleaning

**NOTICE!** Damage resulting from rough handling. Handle the product with care. Do not knock or bend. Set down carefully.

**NOTICE!** Scratch-sensitive surface. Risk of corrosion. Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads soaked with cleaning solution only.

1. Completely immerse the fully disassembled product in the cleaning solution.  
Make sure that the cleaning solution reaches all accessible surfaces. When placing hollow components into the bath, ensure that all air escapes completely from the lumens.
2. While the product soaks in the solution, use a soft wipe or suitable brush to remove all visible residue from its exterior surfaces.  
Brush all surfaces of products with hidden crevices, lumens with working channels, or complex geometries, for at least **1 minute** or until no more residue can be removed.
3. These surfaces should then be rinsed at least **5 times** thoroughly with the cleaning solution. Use a disposable 20 ml syringe.

#### Complete stage II: Intermediate rinsing

1. Fully rinse all accessible surfaces **3 times** for at least **1 minute** each.  
Use fresh water for each rinse.
2. Allow residual water to drip off for a sufficient length of time.

**Complete stage III: Disinfection**

1. Completely immerse the product in the disinfectant solution.

Make sure that the disinfectant solution reaches all accessible surfaces.

Remove all adherent air bubbles from the component surfaces. When placing hollow components into the bath, ensure that all air escapes completely from the lumens.

2. Thoroughly rinse all surfaces of products with hidden crevices, lumens with working channels, or complex geometries at least **5 times**. Use a disposable 20 ml syringe.

**Complete stage IV: Final rinsing**

1. Fully rinse the product in **3 rinses** of **2 minutes** each. Use fresh water for each rinse.
2. Thoroughly rinse all surfaces of products with hidden crevices, lumens with working channels, or complex geometries at least **5 times**.

Use fresh water for each rinse.

3. Allow residual water to drip off for a sufficient length of time.

**Complete stage V: Drying**

1. Dry the product using a soft, lint-free wipe.

Use medical-quality filtered compressed air (pmax = 0.5 bar) to dry the areas you cannot reach with the wipe.

2. After drying, visually inspect the products in a well-lit area (see [section 3.3.1](#)) and make sure they are undamaged and free of visible residue.

Use a magnifying glass. Repeat the cleaning process if necessary.

Immediately set aside any damaged products.

**4.7 Automated Cleaning and Thermal Disinfection****4.7.1 Overview**

| Stage | Work step            | Temperature (°C/°F) | Time (min) | Water quality  | Cleaning solution / notes   |
|-------|----------------------|---------------------|------------|----------------|---|
| I     | Pre-rinsing          | < 25/77             | 2          | Drinking water | ---   |
| II    | Cleaning             | 55/131              | 10         | Drinking water | Alkaline detergent, e.g. 0.5 % neodisher® MediClean forte (5 ml/L), ph > 10 |
| III   | Rinsing I            | > 10/50             | 1          | Drinking water | ---   |
| IV    | Rinsing II           | > 10/50             | 1          | DI water*      | ---   |
| V     | Thermal disinfection | > 90/194            | 5          | DI water*      | A <sub>0</sub> value > 3000   |
| VI    | Drying               | ---                 | ---        | ---            | ---   |

Table 4-2: Overview of automated cleaning and thermal disinfection.

\* DI water = Deionized water (demineralized, low-germ, max. 10 germs/ml, as well as low in endotoxins, max. 0.25 endotoxin units/ml)

**Note:** Automated reprocessing can cause color-anodized or plastic components (e.g. serial rings and the eyepiece cap) to fade.

### 4.7.2 Performing Automated Cleaning and Thermal Disinfection

#### AUTOCLAVE

*Thermal disinfection should be performed only for products that are labeled as autoclavable.*

*Use a legally marketed washer/disinfector that has been validated as effective. The device must meet the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).*

*The manufacturer recommends dry removal of the products prior to automated cleaning.*

*In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.*

*Use only legally marketed reprocessing trays that are approved for automated cleaning and thermal disinfection. Observe the instructions for use of the tray.*

1. Place all parts in the reprocessing tray.
2. Seal the tray.
3. Place the tray into the washer/disinfector.

Apply a validated loading plan when doing so.

Follow the manufacturer's instructions and instructions for use for the device.

When loading, take care to avoid creating rinsing blind spots.

4. Start the cleaning/disinfection cycle in accordance with the manufacturer's instructions and instructions for use for the device.

**CAUTION!** Risk of scalding when unloading the device. Be sure to wear gloves.

5. Remove the reprocessed product from the washer/disinfector.
6. After drying, visually inspect the product in a well-lit area. Products should be dry, undamaged and free of visible residue.

Use a magnifying glass.

Repeat the cleaning and disinfection process if necessary.

Immediately set aside any damaged products.

## 4.8 Sterilization

### 4.8.1 Performing Steam Sterilization

**AUTOCLAVE**

*Steam sterilization should be performed only for products that are labeled as autoclavable.*

*Use fully desalinated drinking water that meets the requirements of European Standard EN 285.*

*Use only legally marketed reprocessing trays that are approved for steam sterilization in an autoclave. Observe the instructions for use of the tray.*

**CAUTION!** Residue from cleaning agents and organic material can negatively influence the sterilization result. Perform the sterilization only on a properly cleaned and dried product.

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place all parts in the reprocessing tray.
3. Seal the tray.
4. Double-wrap the tray with two single layers of legally marketed sterilization wrap.

**CAUTION!** Risk of scalding when loading the device. Wear suitable gloves.

5. Load the wrapped tray in the sterilizer as per a validated loading plan. Observe the instructions of the sterilizer manufacturer when doing so.
6. Start the steam sterilization with an autoclave in accordance with the device manufacturer's instructions, using the following parameters:

|              |  |
|--------------|--|
| Cycle type   | Pre-vacuum process                               |
| Pulses       | 4  |
| Temperature  | 134 °C (273 °F)                                  |
| Holding time | at least 3 min<br>(effective sterilization time) |

**Note:** It is the responsibility of the operator to establish that the sterilizer used meets the parameters above.

**CAUTION!** Risk of scalding when unloading the device. Wear suitable gloves.

**NOTICE!** Products can be shock sensitive when hot. Avoid knocking or shaking.

**NOTICE!** Damage resulting from a sudden change in temperature. Allow products to cool to room temperature; do not use additional cooling measures.

7. Remove the reprocessed product from the sterilizer.

Ensure that products remain sterile after reprocessing.

#### 4.8.2 STERRAD® Sterilization Process 100S, NX, 100NX

**NOTICE!** STERRAD® sterilization may result in cosmetic device changes that do not usually affect the device function.

**NOTICE!** The product must be completely dry before it can be sterilized in the STERRAD® sterilizer. Loads containing moisture may cause a cycle cancellation.

*The STERRAD® Sterilization Systems, manufactured by Advanced Sterilization Products (ASP), use low-temperature, hydrogen peroxide gas plasma technology for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.*

*Refer to the STERRAD® Sterilization Systems User Guide for detailed instructions for use of any STERRAD® unit, ASP's STERRAD® Sterility Guide (SSG) at [www.sterradsterilityguide.com](http://www.sterradsterilityguide.com) or contact ASP customer service.*

*Use only legally marketed reprocessing trays that are approved for STERRAD® sterilization. Observe the instructions for use of the tray.*

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place all parts in the reprocessing tray.
3. Place a STERRAD® indicator strip in the tray.
4. Seal the tray.
5. Double-wrap the tray with 2 single layers of legally marketed sterilization wrap.
6. Load the wrapped tray in the sterilizer.  
  
Position the tray so that the sterilization medium can act thoroughly on all areas. Do not allow any item to touch the wall of the sterilizer.
7. Start the sterilization cycle in accordance with the manufacturer's instructions and instructions for use for the sterilizer.

STERRAD® sterilization was validated for the following cycles:

- STERRAD® 100S Short Cycle
  - STERRAD® NX Standard Cycle
  - STERRAD® 100NX Standard Cycle
8. Remove the reprocessed product from the sterilizer.

Ensure that products remain sterile after reprocessing.

## 5 Maintenance and Repairs

### 5.1 Troubleshooting

| Issue                                   | Possible causes   | Fixes  |
|---|---|--|
| Image cloudy                            | Glass surfaces soiled   | Manual cleaning as per <a href="#">section 4.6</a> , then reprocess, check the water quality   |
|   | Leaky, defective lens system  | Send the endoscope in for repair   |
| Image too dark, too little illumination | Glass surfaces soiled   | Manual cleaning as per <a href="#">section 4.6</a> , then reprocess, check the water quality   |
|   | Unsuitable light guide  | Use a suitable light guide   |
|   | Light guide incorrectly attached to endoscope   | Check that the light guide is seated properly  |
|   | Defective fiber optics  | Check fiber optics as per <a href="#">section 3.3</a>  |
|   | Defective light guide or light source   | Check the light guide and light source   |
| Image is yellowed                       | Soiled fiber optics   | Manual cleaning as per <a href="#">section 4.6</a> , then reprocess, check the water quality   |
|   | Light guide soiled or defective   | Check the light guide (e.g. by illuminating a white surface)   |
| Corrosion, staining, discoloration      | Inadequate cleaning (e.g. remaining protein residue)  | Manual cleaning as per <a href="#">section 4.6</a> , if necessary, rub thoroughly, then reprocess  |
|   | Insufficient rinsing between the reprocessing stages, especially prior to sterilization   | Rinse sufficiently between the reprocessing stages   |
|   | Excessive chloride concentration  | Check water quality  |
|   | Excessive concentration of minerals (e.g. limescale) or organic substances  | Check the water quality, only use demineralized water if necessary   |
|   | Heavy metal ions and/or silicates, elevated concentration of iron, copper, and manganese in the water or sterilizer steam                   | Check the water quality, only use demineralized water if necessary   |
|   | Contaminated, too frequently reused cleaning and disinfecting solution  | The cleaning and disinfecting solution should be replaced regularly  |
|   | Extraneous rust, e.g. resulting from rust contained in steam or damaged or non-corrosion-resistant instruments reprocessed at the same time | Check the supply systems, if multiple items are being reprocessed together, pay attention to material compatibility and signs of prior damage; prevent the different items from touching one another |
|   | Contact corrosion   | Avoid contact with other devices   |

Table 5-1: Troubleshooting table.

## 5.2 Repair

Should you need to arrange a repair for the device, please contact our Service Department.

When you return the device, please also enclose a detailed **description of the fault**.

### **WARNING**

#### **Contaminated device**

Risk of infection

- > Reprocess the product prior to shipping  
(*chapter 4*)

Only ship devices that have been properly reprocessed.

The device should be returned in its original packaging, if possible.

Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

## 6 Product Data

### 6.1 Technical Data

|                                     |                 |
|-------------------------------------|-----------------|
| <b>Item no.</b>                     | <b>21.0003a</b> |
| Direction of view                   | 0°              |
| Field of view                       | Wide angle      |
| Working length                      | 175 mm          |
| Outer diameter of insertion portion | 4 mm            |

|                                     |                 |
|-------------------------------------|-----------------|
| <b>Item no.</b>                     | <b>21.0008a</b> |
| Direction of view                   | 30°             |
| Field of view                       | Wide angle      |
| Working length                      | 175 mm          |
| Outer diameter of insertion portion | 4 mm            |

|                                     |                 |
|-------------------------------------|-----------------|
| <b>Item no.</b>                     | <b>21.0011a</b> |
| Direction of view                   | 70°             |
| Field of view                       | Wide angle      |
| Working length                      | 175 mm          |
| Outer diameter of insertion portion | 4 mm            |

|                                     |                 |
|-------------------------------------|-----------------|
| <b>Item no.</b>                     | <b>21.0101a</b> |
| Direction of view                   | 45°             |
| Field of view                       | Wide angle      |
| Working length                      | 175 mm          |
| Outer diameter of insertion portion | 4 mm            |

|                                     |                 |
|-------------------------------------|-----------------|
| <b>Item no.</b>                     | <b>21.0617a</b> |
| Direction of view                   | 0°              |
| Field of view                       | Standard        |
| Working length                      | 187 mm          |
| Outer diameter of insertion portion | 2.7 mm          |

|                                     |                 |
|-------------------------------------|-----------------|
| <b>Item no.</b>                     | <b>21.0619a</b> |
| Direction of view                   | 30°             |
| Field of view                       | Standard        |
| Working length                      | 187 mm          |
| Outer diameter of insertion portion | 2.7 mm          |

|                                     |                 |
|-------------------------------------|-----------------|
| <b>Item no.</b>                     | <b>21.0621a</b> |
| Direction of view                   | 70°             |
| Field of view                       | Standard        |
| Working length                      | 187 mm          |
| Outer diameter of insertion portion | 2.7 mm          |

## 6.2 Ambient Conditions

### Transport and storage conditions

|                       |                   |
|-----------------------|-------------------|
| Temperature           | -20 °C to +70 °C  |
| Relative air humidity | 5 % to 95 %       |
| Atmospheric pressure  | 70 kPa to 106 kPa |

Reprocessed devices must be protected against recontamination, in a dark, dry, well-ventilated, dust-free and light-controlled location with a consistent temperature.

Direct sunlight, high temperatures, high humidity or radiation can damage the device or present a risk of infection.

When placing into storage, make sure that the device cannot be damaged by other instruments. It is therefore best to store the device individually or use containers in which it can be secured.

### Operating conditions

|                       |                   |
|-----------------------|-------------------|
| Temperature           | +15 °C to +35 °C  |
| Relative air humidity | 5 % to 95 %       |
| Atmospheric pressure  | 70 kPa to 106 kPa |

### CAUTION

#### Non-observance of the ambient conditions

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Observe operating conditions as well as transport and storage conditions

## 6.3 Spare Parts and Accessories

Use original spare parts and accessories only.

| Image   | Designation                                  | Item number |
|---|--|-------------|
|  | Endoscope light guide adapter, Storz systems | 05.0114z    |
|  | Endoscope light guide adapter, Wolf systems  | 05.0116b    |

Table 6-1: Spare parts and accessories.

---

## 7 Disposal

### **WARNING**

#### **Contaminated device**

Risk of infection

- > Reprocess the product prior to disposal  
(*chapter 4*)

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

SCHÖLLY FIBEROPTIC GMBH  
Robert-Bosch-Str. 1-3  
79211 Denzlingen  
Germany