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Instructions For Use
Camera Head Full HD autoclavable
Camera Head for Medical Endoscopy
(Not valid for USA)

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Table of Contents		Page
0	Important Information about this Document.....	34
0.1	Scope of Validity, Identification, Purpose.....	34
0.2	Target Group.....	34
0.3	Using and Storing this Document.....	34
0.4	Supplementary Documents.....	34
1	General Information about the Device.....	35
1.1	Scope of Delivery.....	35
1.2	Product Description.....	35
1.2.1	Performance Characteristics and Function.....	35
1.2.2	Visual Overview.....	36
1.3	Usage.....	37
1.3.1	Intended Use.....	37
1.3.2	Indications.....	37
1.3.3	Contraindications.....	37
1.4	Conformity.....	37
1.4.1	Standards and Directives.....	37
1.4.2	Medical Device Classification.....	37
1.5	Marking.....	38
1.5.1	Pictograms and Information on the Device and Packaging.....	38
1.5.2	Pictograms in this Document.....	39
1.6	Service Department Contact Details.....	39
1.7	Obligation to Report Serious Incidents.....	39
2	General Safety Information.....	40
2.1	Safety Messages in this Document.....	40
2.1.1	Safety Messages at the Start of a Chapter.....	40
2.1.2	Safety Messages in the Body of the Text.....	41
2.2	Product Safety.....	41
2.2.1	Basic Safety Notices.....	41
2.2.2	Staff Qualifications.....	42

3 Use.....	43
3.1 Safety Notices.....	43
3.2 Camera Head Buttons.....	44
3.2.1 General Information and Default Button Settings.....	44
3.2.2 Changing the Layout of Camera Head Buttons.....	44
3.2.3 Navigating in the Menu using the Camera Head Buttons.....	45
3.3 Performing a Visual Inspection.....	45
3.4 Connecting up the Camera Head.....	45
3.5 Connecting an Endocoupler.....	46
3.6 Stopping Usage and completing Pre-Cleaning.....	46
3.7 Storage and Transport.....	46
3.8 Troubleshooting.....	47
3.9 Repair.....	48
4 Reprocessing.....	49
4.1 Safety Notices.....	49
4.2 Staff Qualifications.....	50
4.3 Validated Procedures.....	51
4.4 Cleaning and Disinfecting Agents.....	51
4.5 Manual Cleaning and Disinfection.....	52
4.5.1 Overview.....	52
4.5.2 Performing Manual Cleaning and Disinfection.....	52
4.6 Automated Cleaning and Thermal Disinfection.....	55
4.6.1 Overview.....	55
4.6.2 Performing Automated Cleaning and Thermal Disinfection.....	56
4.7 Sterilization.....	57
4.7.1 Performing Steam Sterilization.....	57
4.7.2 STERRAD® Sterilization Systems 100S, NX, 100NX.....	58
5 Product Data.....	59
5.1 Technical Data.....	59
5.2 Ambient Conditions.....	59
6 Disposal.....	60

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0 Important Information about this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following product:

Item name: **Camera Head Full HD autoclavable**

Item number: **95-3904**

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, medical engineers and employees of central sterile supply who are entrusted with the installation, operation, maintenance, and reprocessing of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

0.4 Supplementary Documents

Instructions for use of the controller

Be sure to observe the instructions for use of the controller (95-3960, 95-3970, 95-3980, 95-3981, 95-3985, 95-6930, 95-6980, ELEVCCU, ELEVCCUVR, ELEV4KCCU) in regards to the use and configuration options of the camera head buttons.

Instructions for use of the endocoupler

To safely use this device, also observe the instructions for use of the endocoupler being used.

1 General Information about the Device

1.1 Scope of Delivery

The scope of delivery for the device includes:

- 1x camera head
- 1x instructions for use

Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Product Description

1.2.1 Performance Characteristics and Function

Compatible controllers	<p>The present camera head can be connected to the controllers with the following item numbers and is compatible with these only: 95-3960, 95-3970, 95-3980, 95-3981, 95-3985, 95-6930, 95-6980, ELEVCCU, ELEVCCUVR, ELEV4KCCU.</p> <p>The combination of the controller and camera head result in a camera based on CMOS technology.</p>
Full HD – 1080p	<p>When used with a compatible monitor, the camera delivers a native full HD image resolution using progressive scanning (1080p).</p>
Camera head buttons	<p>The camera head features three fully programmable camera head buttons that can be used to control image display, record images or video, or to navigate and change the settings in the controller's configuration menu.</p>
C-mount connection	<p>The present camera head has a threaded port that allows for the mounting of connection type C endocouplers.</p> <p>In the following, the term “endocoupler” is used as a synonym for “objective”.</p>
Autoclavable	<p>For use under sterile conditions, the camera may be autoclaved.</p>

1.2.2 Visual Overview


Figure 1-1: The “Full HD autoclavable“ camera head model.

- ❶ Programmable camera head buttons
- ❷ C-mount thread for mounting connection type C endocouplers

1.3 Usage

1.3.1 Intended Use

The application of a camera control unit in combination with the camera head and an endocoupler serves the visualization of the intracorporal surgical area during endoscopic diagnostic and surgical procedures. The camera control unit in combination with a camera head is only used for visualization and not for diagnosis. Refer also to the intended use of the relevant additional components.

1.3.2 Indications

Minimally invasive surgical application.

Surgery, general endoscopic procedure.

1.3.3 Contraindications

The use of the camera components and its accessories is contraindicated if endoscopic procedures are contraindicated for any reason. As is appropriate with any surgical procedure, consideration must be given to patient size and workspace volume when using the camera components. Depending on the patient's disease, there may exist contraindications that rely on the general condition of the patient or the specific disease pattern. The decision to perform an endoscopic procedure rests with the responsible surgeon and should be made on the basis of an individual risk-benefit analysis.

1.4 Conformity

1.4.1 Standards and Directives

The device consisting of controller and camera head meets the requirements of the following applicable standards:

- **IEC 60601-1** Medical electrical equipment - Part 1:
General requirements for basic safety and essential performance
- **IEC 60601-1-2** Medical electrical equipment - Part 1-2:
Electromagnetic compatibility
- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:
Particular requirements for basic safety and essential performance of endoscopic equipment

1.4.2 Medical Device Classification

The product is a Class I medical device within the meaning of applicable **EU guidelines**.

1.5 Marking

1.5.1 Pictograms and Information on the Device and Packaging

	Adhere to the instructions for use
	The device is suitable for autoclave sterilization
	CE marking
	Item number
	Serial number
	Medical Device
	Manufacturer
	Date of manufacture
	Caution (IEC 60601-1 3rd edition) / Take note of accompanying documents (IEC 60601-1 2nd edition)
	MR unsafe
	Type CF defibrillation-proof applied part according to IEC 60601-1
	Release
	Lock
	Permissible storage and transport temperature
	Permissible relative air humidity during storage and transport
	Permissible atmospheric pressure during storage and transport
	Caution: Federal (US) law restricts this device to sale by or on order of a physician
	Separate collection for WEEE (waste of electrical and electronic equipment)

1.5.2 Pictograms in this Document

General warning sign

Electricity hazard warning sign

Biohazard warning sign, risk of infection

1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

1.7 Obligation to Report Serious Incidents

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 General Safety Information
2.1 Safety Messages in this Document
2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

⚠ WARNING

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

⚠ CAUTION

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

NOTICE

Nature and origin of the risk of material damage

Potential consequences of non-observance

> Preventive action

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

Keep the original packaging. Transport and store the product in its original packaging and use it to return goods if service support is required.

If you detect any faults or malfunctions, inform us immediately.

-  **WARNING!** Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.
-  **WARNING!** Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.
-  **WARNING!** Risks from the arrangement, setup, combination, or properties of connected or surrounding devices or equipment. Follow the instructions for use of the respective devices. Perform a risk assessment.
-  **WARNING!** Risk of suffocation. Keep packaging material out of reach of children.
-  **WARNING!** Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.
-  **WARNING!** Using the device during operation of a defibrillator. Risk to persons. Before discharging a defibrillator, remove device from the surgical field.
-  **CAUTION!** Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. Do not kink, crush or strongly bend camera cable. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.
-  **WARNING!** Improper handling, maintenance, and use carries a risk to the patient, user and other persons, or can lead to premature wearing of the device.

2.2.2 Staff Qualifications

Specific qualifications are required for installation, operation, and maintenance of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

3 Use

3.1 Safety Notices

WARNING

Use of non-sterile parts

Infection risk for the patient

- > Use only properly reprocessed endocouplers, endoscopes and endoscopic accessories
- > Always carry out a visual inspection prior to use

WARNING

Displaying a recorded image instead of the live image

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera at all times

WARNING

Interactions between devices in simultaneous use (e.g. lasers, electro-surgery)

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the required type BF or type CF classification requirements according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active HF electrodes
- > Do not activate HF electrodes in the vicinity of flammable gases or liquids
- > Before using HF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

3.2 Camera Head Buttons


3.2.1 General Information and Default Button Settings

Each camera head button can be programmed to trigger different functions by tapping briefly or by holding down.

The buttons can be fully programmed in the menu, subject to the following restrictions:

- It is not possible to change the way in which the menu itself is opened;
- White balance, Light source on/off, Selective color enhancement (SCE), Smoke reduction, Noise reduction and Grid removal can only be assigned to buttons 1 and 2, activated by holding down the button.

The **Default settings** for the camera head buttons are as follows:

	Button 1	Button 2	Button 3
Tap briefly	Brightness +	Brightness -	Take photo
Press and hold down	Light source on/off	White balance	Open menu

Table 3-1: Default settings for the camera head buttons.


3.2.2 Changing the Layout of Camera Head Buttons

Navigate in the menu with the arrow buttons and confirm your selection each time using the Menu button.

The following example describes how to program button 3 for edge enhancement.

1. Press the **Menu** button on the camera head for approximately 2 seconds.
 - > The main menu appears on the screen.
2. Navigate to the **Options** menu item and select.
3. Navigate to **Button layout** and select.
4. Navigate to **M short** and select.
5. Navigate to the **Edge enhancement** option and select.
 - > Your selection is displayed in white type.
 - > You can now use button 3 to control edge enhancement.

3.2.3 Navigating in the Menu using the Camera Head Buttons

Entry into menu:	Press and hold button 3.
Upwards in the menu:	Press button 1 briefly.
Downwards in the menu:	Press button 2 briefly.
Confirm selection in menu:	Press button 3 briefly.
Leave the menu:	Press and hold button 3.

3.3 Performing a Visual Inspection

A visual inspection must be performed prior to each medical procedure.

1. Check that the glass surfaces of the camera head are completely clean.

To do this, direct the camera head onto a clean, white surface. If the glass is dirty, you will see spots or shadows on the screen that cannot be seen on the surface.

NOTICE! Scratching of the glass. Use a soft, lint-free wipe and not a brush.

Clean soiled glass surfaces with 70 % ethanol.

2. Check that the camera cable is undamaged.

Do not use the camera head if the camera cable is kinked, broken, or twisted.

3. Check that the contacts on the camera cable connection plug are free from moisture and dirt.

Do not use the camera head if there is any moisture in the plug or contacts are soiled.

4. Inspect the housing for external damage.

Do not use the device if there is any external damage to the housing.

3.4 Connecting up the Camera Head

1. Slide the connection plug fully into the connection socket on the controller until it clicks into place.

3.5 Connecting an Endocoupler

C-mount standard endocouplers with a flange focal distance of 17.5 mm (0.69 in) can be mounted on the camera head.

To achieve a screen-filling circular reproduction of images as supplied by most rigid endoscopes with diameters from 4 mm to 10 mm, zoom endocouplers with a focal length of 13 - 29 mm can be used, such as manufactured among others by W.O.M. (Autoclavable HD-Zoom Coupler, type EZ241, focal length (f) = 14 - 29 mm) or Richard Wolf (RIWO zoom lens, REF 85261504, focal length (f) = 13 - 29 mm).

Comparable endocouplers can also be used.

1. Remove any protective caps from the endocoupler, if present.
2. Screw the endocoupler onto the camera head turning it clockwise.

Note: The user is responsible at all times for maintaining the sterile barrier.

3.6 Stopping Usage and completing Pre-Cleaning

*Complete the pre-cleaning process **immediately after use.***

1. Remove and discard the sterile drape, if present.
2. Unscrew the endocoupler manually by turning it counterclockwise and proceed to reprocessing as per manufacturer's instructions.
3. Uncouple the endoscope and the light guide.
4. Disconnect the camera head from the controller.
5. Switch off the electrical equipment.
6. Pre-clean the camera head including cable and plug at the conclusion of the surgery by wiping it with a lint-free and wetted, but not dripping wipe until no more residue can be seen.
7. Arrange for reprocessing.

Ensure that the product is reprocessed within **6 hours**.

3.7 Storage and Transport

Store the device in a dry, well-ventilated, temperature-controlled room, where it is protected from dust.

When placing into storage or storing temporarily, make sure that the device cannot be damaged by its surroundings.

To avoid damage, keep out of direct sunlight and away from radioactivity and strong electromagnetic radiation.

The ambient conditions for transport, storage, and operation are described in section 5.2.

3.8 Troubleshooting

Issue	Possible causes	Fixes
Image too dark, too little illumination	Glass surfaces on camera head are soiled	Clean glass surfaces (<i>chapter 4</i>)
	Stubborn residue on the glass surfaces	Remove residue (<i>chapter 4</i>)
	Light guide defective	Connect a new light guide. Send in for repair
	Camera brightness set too low	Increase camera brightness
	Endoscope optical system defective	If the image is too dark even without the camera, use a different endoscope and send in the endoscope for repair
	Light output set too dark	Increase the light output at the light source
Image blurred, cannot be focused	Endocoupler does not correspond to C-mount standard	Connect a suitable endocoupler
	Endocoupler unsuitable for application	Connect a suitable endocoupler
	Endocoupler soiled	Follow the manufacturer's instructions for use of the device
No image on the monitor	Connection cable not connected or defective	Connect the camera head to the controller. Check the plug for moisture. If the cable is defective, send it in for repair

Table 3-2: Troubleshooting table.

3.9 Repair

Should you need to arrange a repair for the device, please contact one of our subsidiaries. Contact details can be found on the back of these instructions for use.

When you send in equipment, please enclose as accurate a **fault description** as possible, and record the item number and serial number of the product on the delivery note. You will find these details on the nameplate.

WARNING

Contaminated device

Risk of infection

- > Reprocess the product prior to shipping (*chapter 4*)
- > Clearly identify contaminated products as such

Send only thoroughly cleaned goods in for repair.

For shipping, please use the original packaging whenever possible. Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

4 Reprocessing

4.1 Safety Notices

WARNING

If the product is potentially contaminated with Creutzfeldt-Jakob pathogens, it is not suitable for reprocessing

Risk of cross-contamination during use and reprocessing

- > Do not reprocess products that are potentially contaminated
- > Products that are potentially contaminated must be disposed of

CAUTION

Improper cleaning and disinfection

Risk to patient as a result of inadequate cleaning and disinfection, damage to the device

- > Use an automatic washer/disinfector that meets the performance requirements of ISO 15883-1
- > Load the device with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Ensure that the washer/disinfector is properly maintained
- > Only use cleaning and disinfecting agents that are approved for the product
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing cleaning and disinfecting agents (active ingredients: aldehyde, alcohol)

NOTICE

Improper sterilization

Risk of damage to the device

- > Use only approved sterilization processes

NOTICE

Unsuitable cleaning agents and process chemicals

Corrosion damage, premature aging, and visible material changes

- > Only use cleaning and disinfecting agents that are approved for the product
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials
- > Avoid contact with chloride-based solutions
- > If the products should come into contact with chloride-based solutions, rinse thoroughly with deionized water and dry completely

NOTICE

Use of ultrasonic bath

Damage to the device

- > The product must not be exposed to the ultrasonic bath

4.2 Staff Qualifications

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must always fall under the responsibility of qualified personnel who have the necessary knowledge and expertise.

This knowledge and expertise can be acquired by completing further training in this field or as a result of completing a dedicated qualification with practical experience, supplemented by appropriate further training measures, where required.

4.3 Validated Procedures

If not using a sterile drape, you must employ the reprocessing procedure.

Reprocessing comprises the following steps

- Pre-cleaning immediately after use
- Cleaning and disinfection (manual or automated) and
- Sterilization

The procedures specified in this document were validated for efficacy.

It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

4.4 Cleaning and Disinfecting Agents

The **manual cleaning** process was validated using:

- Cidezyme®/Enzol® (Johnson & Johnson)

The **manual disinfection** process was validated using:

- Cidex® OPA (Johnson & Johnson)

The **automated cleaning** process was validated using:

- neodisher® MediClean forte 0.5 %
(Chemische Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning and disinfecting agents only.

Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

4.5 Manual Cleaning and Disinfection

4.5.1 Overview

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Cleaning/disinfectant solution
I	Cleaning	as per manufacturer's instructions	2 - 5	Drinking water	Enzymatic cleaning agent: Cidezime®/Enzol®, prepared as per manufacturer's instructions
II	Intermediate rinsing	20-30 / 68-86	3 x ≥ 1	Drinking water	---
III	Disinfection	as per manufacturer's instructions	12	---	0.55 % ortho-phthalaldehyde solution: Cidex® OPA
IV	Final rinsing	20-30 / 68-86	3 x ≥ 1	DI water*	---
V	Drying	---	---	---	---

Table 4-1: Overview of manual cleaning and disinfection.

* DI water = Deionized water (demineralized, low-germ, max. 10 germs/ml, as well as low endotoxin contaminated, max. 0.25 endotoxin units/ml)

4.5.2 Performing Manual Cleaning and Disinfection

WARNING

The product is not sterile after manual cleaning and disinfection.

NOTICE

Scratch-sensitive surface

Risk of corrosion

Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads soaked with alcohol or neutral cleaner only.

Note: Use fresh water for each rinse.

Complete stage I: Cleaning

1. Pre-clean the camera head including cable and plug at the conclusion of the surgery by wiping it with a lint-free and wetted, but not dripping wipe until no more residue can be seen.

2. Unscrew the endocoupler from the camera head manually by turning it counterclockwise. Follow the manufacturer's instructions for use in regards to the reprocessing of the endocoupler.
3. Completely immerse the product for **2 to 5 minutes** in the cleaning solution prepared as per manufacturer's instructions.
4. Make sure that the cleaning solution touches all accessible surfaces and that the product remains immersed throughout the entire cleaning time. When placing the product, ensure that all air escapes from the hidden crevices.
5. While the product is immersed in the cleaning solution, clean it using a soft, clean, lint-free wipe or with a soft bristled brush. Wipe or brush for at least **1 minute** or until no more residue can be seen.
6. During the cleaning, maneuver mobile components of the completely immersed product in all directions as far as they will go **3 times**.
7. Thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least **5 times** with the cleaning solution. Use a disposable 50 ml syringe.

Complete stage II: Intermediate rinsing

1. Completely immerse the product in drinking water (20-30 °C / 68-86 °F). Fully rinse all accessible surfaces in **3 rinses** of at least **1 minute** each.
2. During the rinsing, maneuver mobile components in all directions as far as they will go **3 times**.
3. Thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least **6 times** with drinking water. Use a disposable 50 ml syringe.
4. Allow residual water to drip off for a sufficient length of time.

Complete stage III: Disinfection

1. Completely immerse the product in the disinfectant solution for at least **12 minutes**.
2. Make sure that all accessible surfaces remain immersed in the disinfectant solution throughout the entire disinfecting time. Remove all adherent air bubbles from the component surfaces.
3. Wipe or brush the product, while it is completely immersed in the cleaning solution, using a soft, clean, lint-free wipe or with a soft bristled brush for at least **1 minute**.
4. During the cleaning, maneuver mobile components of the completely immersed in the disinfectant solution product in all directions as far as they will go **3 times**.
5. Thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least **5 times** with disinfectant solution. Use a disposable 50 ml syringe.

Complete stage IV: Final rinsing

1. Completely immerse the product in deionized water (20-30 °C / 68-86 °F) and rinse all accessible surfaces in **3 rinses** of at least **1 minute** each.
2. During the rinsing, maneuver mobile components in all directions as far as they will go **3 times**.
3. Thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least **3 times** with deionized water. Use a disposable 50 ml syringe.
4. Allow residual water to drip off for a sufficient length of time.

Complete stage V: Drying

1. Dry the product using a soft, lint-free wipe. Use medical-quality filtered compressed air ($p_{max} = 0.5 \text{ bar}$) to dry the areas you cannot reach with the wipe.
2. Make sure the product is clean and undamaged. Repeat the cleaning and disinfection process if necessary.

Visual inspection

If any damage is detected, set aside the product immediately.

4.6 Automated Cleaning and Thermal Disinfection

4.6.1 Overview

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Cleaning solution
I	Pre-rinsing	Cold	2	Drinking water	---
II	Cleaning	55/131	10	---	Alkaline cleaner: neodisher® MediClean forte, prepared as per manufacturer's instructions
III	Intermediate rinsing 1	per the equipment manufacturer's standard cycle	1	Drinking water	---
IV	Intermediate rinsing 2	per the equipment manufacturer's standard cycle	1	Drinking water	---
V	Thermal disinfection	90/194	5	Deionized water (A ₀ value: > 3000)	---
VI	Drying	High (98.8/210)	30	---	---

Table 4-2: Overview of automated cleaning and thermal disinfection.

Note: Color anodized parts or plastic components may bleach out during automated cleaning.

4.6.2 Performing Automated Cleaning and Thermal Disinfection

AUTOCLAVE

Thermal disinfection should be performed only for products that are labeled as autoclavable.

Use a washer/disinfector that has been validated as effective. The device must meet the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).

The manufacturer recommends dry removal of the products prior to automated cleaning.

In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.

Use only legally marketed reprocessing trays¹ that are approved for automated cleaning and thermal disinfection. Observe the instructions for use of the tray.

1. Place the camera head in the reprocessing tray.
2. Seal the tray.
3. Place the tray into the washer/disinfector.

Apply a validated loading plan when doing so.

Follow the manufacturer's instructions and instructions for use for the device.

When loading, take care to avoid creating rinsing blind spots.

4. Start the cleaning/disinfection cycle in accordance with the manufacturer's instructions and instructions for use for the device.

CAUTION! Risk of scalding when unloading the device. Be sure to wear gloves.

5. Remove the reprocessed product from the washer/disinfector.
6. After drying, visually inspect the product in a well-lit area. Products should be dry, undamaged and free of visible residue.

Use a magnifying glass.

Repeat the cleaning and disinfection process if necessary.

Immediately set aside any damaged products.

¹ for automated cleaning validation the camera head was placed in a standard reprocessing tray (OM-1002-SY)

4.7 Sterilization

4.7.1 Performing Steam Sterilization

Steam sterilization should be performed only for products that are labeled as autoclavable.

Use fully desalinated drinking water that meets the requirements of European Standard EN 285.

Use only legally marketed reprocessing trays¹ that are approved for steam sterilization in an autoclave. Observe the instructions for use of the tray.

CAUTION! Residue from cleaning agents and organic material can negatively influence the sterilization result. Perform the sterilization only on a properly cleaned and dried product.

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place the camera head in the reprocessing tray.
3. Seal the tray.
4. Double-wrap the tray with two single layers of legally marketed sterilization wrap.

CAUTION! Risk of scalding when loading the device. Wear suitable gloves.

5. Load the wrapped tray in the sterilizer as per a validated loading plan. Observe the instructions of the sterilizer manufacturer when doing so.
6. Start the steam sterilization with an autoclave in accordance with the device manufacturer's instructions, using the following parameters:

Cycle type	Pre-vacuum process
Pulses	3
Temperature	134 °C (273 °F)
Holding time	3 - 5 min (effective sterilization time)

Note: It is the responsibility of the operator to establish that the sterilizer used meets the parameters above.

CAUTION! Risk of scalding when unloading the device. Wear suitable gloves.

NOTICE! Products can be shock sensitive when hot. Avoid knocking or shaking.

NOTICE! Damage resulting from a sudden change in temperature. Allow products to cool to room temperature; do not use additional cooling measures.

7. Remove the reprocessed product from the sterilizer.
Ensure that products remain sterile after reprocessing.

¹ for steam sterilization validation the camera head was placed in a standard reprocessing tray (OM-1002-SY) that was double-wrapped in polypropylene wrap (Kimguard® KC200 Sterilization Wrap, Kimberly Clark, PC 68030)

4.7.2 STERRAD® Sterilization Systems 100S, NX, 100NX

Performing sterilization

NOTICE! STERRAD® sterilization may result in cosmetic device changes that do not usually affect the device function.

NOTICE! The product must be completely dry before it can be sterilized in the STERRAD® sterilizer. Loads containing moisture may cause a cycle cancellation.

The STERRAD® Sterilization Systems, manufactured by Advanced Sterilization Products (ASP), use low-temperature, hydrogen peroxide gas plasma technology for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.

Refer to the STERRAD® Sterilization Systems User Guide for detailed instructions for use of any STERRAD® unit, ASP's STERRAD® Sterility Guide (SSG) at www.sterradsterilityguide.com or contact ASP customer service.

Use only legally marketed reprocessing trays¹ that are approved for STERRAD® sterilization. Observe the instructions for use of the tray.

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place the camera head in the reprocessing tray.
3. Place a STERRAD® indicator strip in the tray.
4. Seal the tray.
5. Double-wrap the tray with two single layers of legally marketed sterilization wrap.
6. STERRAD® 100S und STERRAD® 100NX: Load the wrapped tray on the top shelf or the bottom shelf of the sterilization chamber.

STERRAD® NX: Load the wrapped tray on the bottom shelf of the sterilization chamber. The top shelf must be removed.

Position the tray so that the sterilization medium can act thoroughly on all areas. Do not allow any item to touch the wall of the sterilizer.

7. Start the sterilization cycle in accordance with the manufacturer's instructions and instructions for use for the sterilizer.

STERRAD® sterilization was validated for the following cycles:

- STERRAD® 100S Short Cycle (top shelf or bottom shelf)
- STERRAD® NX Standard Cycle (bottom shelf only)
- STERRAD® 100NX Standard Cycle (top shelf or bottom shelf)

8. Remove the reprocessed product from the sterilizer.

Ensure that products remain sterile after reprocessing.

¹ for sterilization validation the camera head was placed in a standard reprocessing tray (OM-1002-SY) that was double-wrapped in polypropylene wrap (H400 Sterilization Wrap, HALYARD* Health, PC 68236)

5 Product Data

5.1 Technical Data

Dimensions without endocoupler (W x H x D)	40 x 45 x 150 mm
Weight without endocoupler	230 g without cable
Image sensor	1 x 1/3" CMOS, progressive scan
Resolution	1920 x 1080 pixels
Endocoupler connection type	C-Mount
Camera cable length	3.5 m
Degree of protection	IP X7
Applied part as per IEC 60601-1	CF defibrillation-proof

5.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	50 kPa to 106 kPa

Operating conditions

Temperature	+10 °C to +34 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

CAUTION

Non-observance of the ambient conditions

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

> Observe operating conditions as well as transport and storage conditions

6 Disposal

WARNING

Contaminated device

Risk of infection

> Reprocess the product prior to disposal (*chapter 4*)

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

The product is packed in a polypropylene (PP) plastic case (0.5 kg) with a polyurethane (PU) soft foam inlay (0.12 kg), polyethylene (PE) bags (0.02 kg) as well as a corrugated cardboard box (0.54 kg).

Any product carrying this symbol must be disposed of separately through dedicated electrical and electronic devices recycling. Within the EU, such disposal is taken care of by the manufacturer free of charge.

SCHÖLLY FIBEROPTIC GMBH
Robert-Bosch-Str. 1-3
79211 Denzlingen
Germany