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
**NIR FI Kamerakopf Full HD Zoom**

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# **Instructions For Use**

**NIR FI Camera Head Full HD Zoom**

Camera Head for NIR Fluorescence Imaging

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---

	<b>Content</b>	<b>Page</b>
<b>0</b>	<b>Important Information about this Document.....</b>	<b>42</b>
0.1	Scope of Validity, Identification, Purpose.....	42
0.2	Target Group.....	42
0.3	Using and Storing this Document.....	42
0.4	Supplementary Documents.....	42
<b>1</b>	<b>General Information about the Device.....</b>	<b>43</b>
1.1	Scope of Delivery.....	43
1.2	Product Description.....	43
1.2.1	Performance Characteristics and Function.....	43
1.2.2	Visual Overview.....	44
1.3	Usage.....	45
1.3.1	Intended Use.....	45
1.3.2	Indications.....	45
1.3.3	Contraindications.....	45
1.3.4	User Profile.....	45
1.3.5	Compatibility with HF Surgical Devices.....	45
1.3.6	Components Required for Operation.....	45
1.4	Conformity.....	46
1.4.1	Standards and Directives.....	46
1.4.2	Protection against Electric Shock.....	46
1.4.3	Medical Device Classification.....	46
1.5	Marking.....	47
1.5.1	Pictograms and Information on the Device and Packaging.....	47
1.5.2	Pictograms in this Document.....	48
1.6	Service Department Contact Details.....	48
1.7	Obligation to Report Serious Incidents.....	48
<b>2</b>	<b>General Safety Information.....</b>	<b>49</b>
2.1	Safety Messages in this Document.....	49
2.1.1	Safety Messages at the Start of a Chapter.....	49
2.1.2	Safety Messages in the Body of the Text.....	50
2.2	Product Safety.....	50
2.2.1	Basic Safety Notices.....	50
2.2.2	Staff Qualifications.....	51

---

<b>3 Use.....</b>	<b>52</b>
3.1 Safety Notices.....	52
3.2 Camera Head Buttons.....	52
3.2.1 General Information and Default Button Settings.....	52
3.2.2 Changing the Button Function.....	53
3.2.3 Navigation in the Menu of the camera control unit using the camera head buttons.....	54
3.3 Visual Inspection.....	54
3.4 Connection to the Controller.....	55
3.5 Connection of an Endoscope.....	55
3.6 Adjustment of the Image Focus and Zooming.....	56
3.7 Post-Use and Pre-Cleaning Activities.....	56
3.8 Storage and Transport.....	56
3.9 Troubleshooting Table.....	57
3.10 Repair.....	57
<b>4 Reprocessing.....</b>	<b>58</b>
4.1 Safety Notices.....	58
4.2 Staff Qualifications.....	59
4.3 Validated Procedures.....	60
4.4 Cleaning and Disinfecting Agents.....	60
4.5 Manual Cleaning and Disinfection.....	61
4.5.1 Overview.....	61
4.5.2 Performance.....	61
4.6 Automated Cleaning and Thermal Disinfection.....	64
4.6.1 Manual Pre-Cleaning.....	64
4.6.2 Overview.....	65
4.6.3 Performance.....	65
4.7 Sterilization.....	66
4.7.1 STERRAD® Sterilization Systems 100S, NX, 100NX.....	66
4.7.2 STERIS Sterilization Systems V-PRO®.....	68
4.7.3 STERIS Sterilization SYSTEM 1® / SYSTEM 1E®.....	69
4.7.4 Ethylene Oxide Gas Sterilization.....	69
<b>5 Product Data.....</b>	<b>71</b>
5.1 Technical Data.....	71
5.2 Ambient Conditions.....	71
<b>6 Disposal.....</b>	<b>72</b>

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## **0 Important Information about this Document**

### **0.1 Scope of Validity, Identification, Purpose**

These instructions for use apply to the following product:

Item name: **NIR FI Camera Head Full HD Zoom**

Item number: **95-3908**

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

### **0.2 Target Group**

These instructions for use are intended for physicians, medical assistants, medical engineers and employees of central sterile supply who are entrusted with the installation, operation, maintenance, and reprocessing of the device.

### **0.3 Using and Storing this Document**

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

### **0.4 Supplementary Documents**

Instructions for use of the camera control unit

Be sure to observe the instructions for use of the camera control unit in regards to the use and configuration options of the camera head buttons (95-3981, 95-3985).

# 1 General Information about the Device

## 1.1 Scope of Delivery

The scope of delivery for the device includes:

- one camera head
- one set of instructions for use

Check scope of delivery!

Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact your distributor. You will find contact details on the back of this document.

## 1.2 Product Description

### 1.2.1 Performance Characteristics and Function

Compatible controllers	The present camera head can be connected to the controllers with the following item numbers and is compatible with these only: 95-3981, 95-3985.
CMOS camera	The combination of camera head and camera control unit result in a camera based on CMOS technology that is used to provide endoscopic real-time video for display on a monitor.
NIR fluorescence	The camera head is sensitive in the visible and near-infrared spectrum to provide visualization of near-infrared light for NIR fluorescence imaging.
Full HD – 1080p	When used with a compatible monitor, the camera delivers a native full HD image resolution using progressive scanning (1080p).
Camera head buttons	The camera head features three fully programmable camera head buttons that can be used to control image display, record images or video, or to navigate and change the settings in the camera control unit's configuration menu.
Parfocal zoom	The camera head features a manually adjustable focal length of 14.25 – 28 mm, which corresponds to a 2x zoom.  The integrated objective is designed parfocally, so that the image remains sharp at any time.

1.2.2 Visual Overview


Figure 1-1: The “NIR FI Full HD Zoom” camera head model.

- ① Zoom ring
- ② Focusing ring
- ③ Eyepiece locking mechanism
- ④ Protection against inadvertent opening of the locking mechanism
- ⑤ Camera head buttons

---

## 1.3 Usage

### 1.3.1 Intended Use

The camera head is connected to a camera control unit and endoscope to provide real-time video during diagnostic or therapeutic endoscopic procedures.

The camera head is also used for intraoperative fluorescence imaging in the near-infrared range with a suitable fluorescence dye (principally Indocyanine Green) that promotes the fluorescing appearance of the camera image.

### 1.3.2 Indications

The camera head is indicated for use in surgical endoscopic procedures that physicians deem appropriate for the patient wherever a camera is required for intraoperative imaging.

The decision to perform an endoscopic procedure depends on the condition of the patient and remains with the responsible physician.

In order to assure a safe application, standards/guidelines released by official organizations such as scientific societies of endoscopic surgery have to be considered.

### 1.3.3 Contraindications

The use is contraindicated in all cases where endoscopic procedures are contraindicated for any reason.

There are no known contraindications linked directly to the devices.

Medical contraindications that are listed for suitable fluorescence dyes that may be used apply.

For any interventions the responsible physician must decide applying professional standards and based on the general condition of the patient and an individual benefit-risk analysis whether the intended application of the device is admissible.

### 1.3.4 User Profile

Intended to be used by trained healthcare professionals in a professional healthcare environment (hospital).

### 1.3.5 Compatibility with HF Surgical Devices

The device is approved for use in conjunction with HF surgical devices. Applying HF surgery does not alter the conditions of the intended use.

### 1.3.6 Components Required for Operation

The camera head can be connected to the camera control units with the following item numbers and is be compatible with these only:

- NIR FI CCU Full HD (95-3981)
- NIR FI CCU 4K (95-3985)

For the complete system requirements for achieving the intended results when using the camera for NIR fluorescence imaging, observe the instructions for use supplied with the light source (05--0761nir).

## 1.4 Conformity

### 1.4.1 Standards and Directives

The device consisting of controller and camera head meets the requirements of the following applicable standards:

- **IEC 60601-1** Medical electrical equipment - Part 1:  
General requirements for basic safety and essential performance
- **IEC 60601-1-2** Medical electrical equipment - Part 1-2:  
Electromagnetic compatibility
- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:  
Particular requirements for basic safety and essential performance of endoscopic equipment

### 1.4.2 Protection against Electric Shock

In combination with the camera control unit the product meets the - defibrillation-proof type CF classification requirements according to IEC 60601-1 for protection against electric shock and is approved for use in combination with applied parts intended for direct conductive contact with the heart.

### 1.4.3 Medical Device Classification

The product is a Class I medical device within the meaning of applicable **EU guidelines**.

## 1.5 Marking

### 1.5.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the product or packaging.

Adhere to the instructions for use

CE marking

Item number

Serial number

Medical Device

Manufacturer

Date of manufacture

Caution (IEC 60601-1 3rd edition) /

Take note of accompanying documents (IEC 60601-1 2nd edition)

MR unsafe

Type CF defibrillation-proof applied part according to IEC 60601-1

Release

Lock

Permissible storage and transport temperature

Permissible relative air humidity during storage and transport

Permissible atmospheric pressure during storage and transport

Caution: Federal (US) law restricts this device to sale by or on order of a physician

Separate collection for WEEE (waste of electrical and electronic equipment)

### 1.5.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Electricity hazard warning sign

Biohazard warning sign, risk of infection

### 1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact your distributor.

You will find contact details on the back of this document.

### 1.7 Obligation to Report Serious Incidents

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**2 General Safety Information**  
**2.1 Safety Messages in this Document**  
**2.1.1 Safety Messages at the Start of a Chapter**

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

<p><b>⚠ WARNING</b></p> <p><b>Nature and origin of the injury risk</b></p> <p>Potential consequences of non-observance</p> <p>&gt; Preventive action</p>
--

The following kind of message refers to a risk that could lead to minor or moderate injury:

<p><b>⚠ CAUTION</b></p> <p><b>Nature and origin of the injury risk</b></p> <p>Potential consequences of non-observance</p> <p>&gt; Preventive action</p>
--

The following kind of message refers to a risk that could lead to material damage:

<p><b>NOTICE</b></p> <p><b>Nature and origin of the risk of material damage</b></p> <p>Potential consequences of non-observance</p> <p>&gt; Preventive action</p>
---

### 2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

**WARNING!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

**CAUTION!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

**NOTICE!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

## 2.2 Product Safety

### 2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Risks despite highest of quality standards

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

Keep the original packaging. Transport and store the product in its original packaging and use it to return goods if service support is required.

If you detect any faults or malfunctions, inform us immediately.

-  **WARNING!** Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.
-  **WARNING!** Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.
-  **WARNING!** Risks from the arrangement, setup, combination, or properties of connected or surrounding devices or equipment. Follow the instructions for use of the respective devices. Perform a risk assessment.
-  **WARNING!** Risk of suffocation. Keep packaging material out of reach of children.
-  **WARNING!** Interference to the live image caused by electromagnetic disturbances (white light and NIR fluorescence imaging mode). Risk to the patient due to altered image representation, system locked in FI mode, live image failure (e.g. delay, frozen image), image interference (e.g. flickering, stripes), altered image orientation. Remove all sources of interference. Observe the recommended minimum distances according to the Information on Electromagnetic Compatibility. Ensure that the monitor is correctly displaying the live image. Ensure that the live image is displayed in the correct orientation (not mirrored).
-  **WARNING!** Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.
-  **WARNING!** Using the device during operation of a defibrillator. Risk to persons. Before discharging a defibrillator, remove device from the surgical field.
-  **CAUTION!** Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. Do not kink, crush or strongly bend camera cable. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.
-  **WARNING!** Improper handling, maintenance, and use carries a risk to the patient, user and other persons, or can lead to premature wearing of the device.

### 2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for installation, operation, and maintenance of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

## 3 Use

### 3.1 Safety Notices

#### **WARNING**

##### **Use of non-sterile parts**

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Always carry out a visual inspection prior to use

#### **WARNING**

##### **Displaying a recorded image instead of the live image**

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera at all times

#### **WARNING**

##### **Interactions between devices in simultaneous use (e.g. lasers, electro-surgery)**

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the required type BF, CF or CF defibrillation-proof classification requirements according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active HF electrodes
- > Do not activate HF electrodes in the vicinity of flammable gases or liquids
- > Before using HF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

## 3.2 Camera Head Buttons


### 3.2.1 General Information and Default Button Settings

Each camera head button can be programmed to trigger different functions by tapping briefly or by holding down.

The buttons can be fully programmed in the menu, subject to the following restrictions:

- It is not possible to change the way in which the menu itself is opened;
- Certain functions (e.g. White balance, Light source on/off, Advanced Views and Noise reduction) can only be assigned to buttons 1 and 2, activated by holding down the button.

The **Default settings** for the camera head buttons are as follows:

	Button 1	Button 2	Button 3
Tap briefly	<p><b>Light source is on:</b> Toggle imaging modes</p> <p><b>Light source is off:</b> Activate light source</p>	Video capture start / stop	Take photo
Press and hold down	<p><b>NIR FI mode is on:</b> Directly return to white light imaging</p> <p><b>White light imaging is on:</b> Deactivate light source</p>	White balance	Open menu

Table 3-1: Default settings for the camera head buttons.

### 3.2.2 Changing the Button Function

#### Changing the function of camera head buttons


Navigate in the menu with the arrow buttons and confirm your selection each time using the Menu button.

The following example describes how to program button 3 for edge enhancement.

Proceed as follows:

1. Press the **Menu** button on the camera head for approx. 2 seconds.
    - > The main menu appears on the screen.
  2. Navigate to the **Options** menu item and select.
  3. Navigate to **Button function** and select.
  4. Navigate to **M short** and select.
  5. Navigate to the **Edge enhancement** option and select.
    - > Your selection is displayed in white type.
    - > You can now use button 3 to control edge enhancement.
- The button function has been changed.

### 3.2.3 Navigation in the Menu of the camera control unit using the camera head buttons

Entry into menu:	Press and hold button 3.
Upwards in the menu:	Press button 1 briefly.
Downwards in the menu:	Press button 2 briefly.
Confirm selection in menu:	Press button 3 briefly.
Leave the menu:	Press and hold button 3.

## 3.3 Visual Inspection

### Performing a visual inspection

*A visual inspection must be performed prior to each medical procedure.*

*Proceed as follows:*


- |                                      |   |
|--------------------------------------|---|
| Are the glass surfaces clean?        | <p>1. Check that the glass surfaces of the camera head are completely clean.</p> <p>To do this, direct the camera head onto a clean, white surface. If the glass is dirty, you will see spots or shadows on the screen that cannot be seen on the surface.</p> <p><b>NOTICE!</b> Scratching of the glass. Use a soft, lint-free wipe and not a brush.</p> <p>Clean soiled glass surfaces with 70 % ethanol.</p> |
| Are the adjusting rings operational? | <p>2. Check that the adjusting ring(s) is/are operational.</p> <p>It should be possible to rotate each ring without it catching.</p> <p>Every ring should respond to turning with slight, uniform resistance.</p> <p>Do not use the camera head if the adjusting ring catches when turned or is too loose.</p>  |
| Is the cable undamaged?              | <p>3. Check that the camera cable is undamaged.</p> <p>Do not use the camera head if the camera cable is kinked, broken, or twisted.</p>  |
| Is the plug dry and undamaged?       | <p>4. Check that the contacts on the camera cable connection plug are free from moisture and dirt.</p> <p>Do not use the camera head if there is any moisture in the plug or contacts are soiled.</p>   |
| Is the housing undamaged?            | <p>5. Inspect the housing for external damage.</p> <p>Do not use the device if there is any external damage to the housing.</p>   |
- This concludes the visual inspection.

### 3.4 Connection to the Controller

#### Connecting up the camera head to the camera control unit

**NOTICE!** Damage to the camera cable by kinking, crushing or sharp bending. Do not kink, crush or strongly bend camera cable.

*Proceed as follows:*

1. Slide the connection plug fully into the connection socket on the camera control unit until it clicks into place.
    - > *The name of the connected camera head model is displayed on the screen.*
- The camera head is connected.

### 3.5 Connection of an Endoscope

#### Connecting an endoscope

*All endoscopes with an ACMI-compatible eyepiece can be connected to the camera head.*

*Proceed as follows:*

1. Attach a sterile drape to the endoscope if necessary.
  - If applicable, observe the instructions for use supplied with the sterile drape.

Note: The user is responsible at all times for maintaining the sterile barrier.

2. Press together the two locking mechanisms on the endocoupler to spread apart the retaining clips.
  3. Insert the endoscope including the eyepiece cup and release the locking mechanisms.
- 
- The diagram shows two U-shaped locking mechanisms, one on the left and one on the right. A double-headed horizontal arrow is positioned between them, indicating that they are to be pushed together to engage the locking mechanism.
4. Secure the endoscope against inadvertent opening. Use the locking device on the endocoupler to do so.
  5. If applicable, draw the sterile drape over the camera head and camera cable toward the camera control until the drape is completely unfolded.
- The endoscope is now connected.

### 3.6 Adjustment of the Image Focus and Zooming

#### Adjusting the image focus and zooming

*Proceed as follows:*

1. Adjust the desired image size using the rear zoom ring.
2. Direct the endoscope at a structure with a sharp outline at the working distance.
3. Adjust the front focusing ring on the camera head until the image becomes sharp.

Parfocal        Once the image focus is sharp, the image also remains sharp when zooming.

- The image focus has been set.

### 3.7 Post-Use and Pre-Cleaning Activities

#### Stopping usage

Note: The endocoupler is permanently attached to the camera head and cannot be detached.

*Proceed as follows:*

1. Remove and discard the sterile drape, if present.
  2. Uncouple the endoscope and the light guide.
  3. Disconnect the camera head from the camera control unit.
  4. Remove any visible surgical residue as completely as possible using a moist but not dripping lint-free wipe.
  5. Switch off the electrical equipment.
- This concludes use and pre-cleaning.

### 3.8 Storage and Transport

Store the device in a dry, well-ventilated, temperature-controlled room, where it is protected from dust.

When placing into storage or storing temporarily, make sure that the device cannot be damaged by its surroundings.

To avoid damage, keep out of direct sunlight and away from radioactivity and strong electromagnetic radiation.

The ambient conditions for transport, storage, and operation are described in section 5.2.

### 3.9 Troubleshooting Table

Issue	Possible causes	Fixes
Image too dark, too little illumination	Glass surfaces on camera head are soiled	Clean glass surfaces ( <i>chapter 4</i> )
	Stubborn residue on the glass surfaces	Remove residue ( <i>chapter 4</i> )
	Light guide defective	Connect a new light guide. Send in for repair
	Camera brightness set too low	Increase camera brightness
	Endoscope optical system defective	If the image is too dark even without the camera, use a different endoscope and send in the endoscope for repair
	Light output set too dark	Increase the light output at the light source
No image on the monitor	Connection cable not connected or defective	Connect the camera head to the controller. Check the plug for moisture. If the cable is defective, send it in for repair

Table 3-2: Troubleshooting table.

### 3.10 Repair

Should you need to arrange a repair for the device, please contact one of our subsidiaries. Contact details can be found on the back of these instructions for use.

When you send in equipment, please enclose as accurate a **fault description** as possible, and record the item number and serial number of the product on the delivery note. You will find these details on the nameplate.

#### **WARNING**

##### **Contaminated device**

Risk of infection

- > Reprocess the device properly prior to shipping (*chapter 4*)
- > Clearly identify contaminated products as such

Send only thoroughly cleaned goods in for repair.

For shipping, please use the original packaging whenever possible. Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

## 4 Reprocessing

### 4.1 Safety Notices

#### **WARNING**

**If the product is potentially contaminated with Creutzfeldt-Jakob pathogens, it is not suitable for reprocessing**

Risk of cross-contamination during use and reprocessing

- > Do not reprocess products that are potentially contaminated
- > Products that are potentially contaminated must be disposed of

#### **CAUTION**

##### **Improper cleaning and disinfection**

Risk to patient as a result of inadequate cleaning and disinfection, damage to the device

- > Use an automatic washer/disinfector that meets the performance requirements of ISO 15883-1
- > Load the washer with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Ensure that the washer/disinfector is properly maintained
- > Only use cleaning and disinfecting agents that are approved for the product
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing cleaning and disinfecting agents (active ingredients: aldehyde, alcohol)

#### **NOTICE**

##### **Improper sterilization**

Risk of damage to the device

- > Use only approved sterilization processes
- > Do not sterilize control unit
- > Do not sterilize camera head with steam or in an autoclave

**NOTICE**

**Unsuitable cleaning agents and process chemicals**

Corrosion damage, premature aging, and visible material changes

- > Only use cleaning and disinfecting agents that are approved for the product
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials
- > Avoid contact with chloride-based solutions
- > If the products should come into contact with chloride-based solutions, rinse thoroughly with DI water and dry completely

**NOTICE**

**Use of ultrasonic bath**

Damage to the device

- > The product must not be exposed to the ultrasonic bath

**4.2 Staff Qualifications**

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must always fall under the responsibility of qualified personnel who have the necessary knowledge and expertise.

This knowledge and expertise can be acquired by completing further training in this field or as a result of completing a dedicated qualification with practical experience, supplemented by appropriate further training measures, where required.

### 4.3 Validated Procedures

If not using a sterile drape, you must employ the reprocessing procedure.

Reprocessing comprises the following steps

- Pre-cleaning immediately after use
- Cleaning and disinfection (manual or automated) and
- Sterilization

The procedures specified in this document were validated for efficacy.

Operator responsibility

It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

### 4.4 Cleaning and Disinfecting Agents

The **manual cleaning** process was validated using:

- Cidezyme®/Enzol® (Johnson & Johnson)

The **manual disinfection** process was validated using:

- Cidex® OPA (Johnson & Johnson)

The **automated cleaning** process was validated using:

- neodisher® MediClean forte 0.5 %  
(Chemische Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning and disinfecting agents only.

Follow the  
manufacturer's  
instructions!

Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

## 4.5 Manual Cleaning and Disinfection

### 4.5.1 Overview

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Cleaning/disinfectant solution
I	Cleaning	as per manufacturer's instructions	2 - 5	Tap water	Enzymatic cleaning agent: Cidezyme®/Enzol®, prepared as per manufacturer's instructions
II	Intermediate rinsing 2x	20-30 / 68-86	2 x ≥ 1	Tap water	---
III	Disinfection	as per manufacturer's instructions	12	---	0.55 % ortho-phthalaldehyde solution: Cidex® OPA, prepared as per manufacturer's instructions
IV	Final rinsing	20-30 / 68-86	3 x ≥ 1	DI water*	---
V	Drying	---	---	---	---

Table 4-1: Overview of manual cleaning and disinfection.

\*DI water = Deionized water (demineralized, low-germ, max. 10 germs/ml, as well as low endotoxin contaminated, max. 0.25 endotoxin units/ml)

### 4.5.2 Performance

#### Complete stage I: Cleaning

Proceed as follows:

**NOTICE!** Damage resulting from rough handling. Handle the product with care. Do not knock, crush or bend. Set down carefully.

**NOTICE!** Scratch-sensitive surface. Risk of corrosion. Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads soaked with alcohol or neutral cleaner only.

Note: The endocoupler is permanently attached to the camera head and cannot be removed.

1. Pre-clean the camera head including cable and plug at the conclusion of the surgery by wiping it with a lint-free and wetted, but not dripping wipe. Before wetting the wipe with the enzymatic cleaning solution (Cidezyme®/Enzol®), it should be prepared per the manufacturer's instructions. Wipe the product until it is free of visible residue, but for at least **one** minute.

Immerse for 2-5 min	2. Completely immerse the product for <b>two</b> to <b>five</b> minutes in the cleaning solution prepared as per manufacturer's instructions.
Allow air bubbles to escape	Make sure that the cleaning solution touches all accessible surfaces and that the product remains immersed throughout the entire cleaning time. When placing the camera head into the bath, ensure that all air escapes completely from the hidden crevices.
Wipe or brush external surfaces	3. While the product is soaking in the cleaning solution, clean it using a soft, clean, lint-free wipe or with a soft bristled brush.
Wipe or brush for at least 1 min	Wipe or brush for at least <b>one</b> minute or until no more residue can be seen.
Move mobile components	During the cleaning, maneuver mobile components of the completely immersed product in all directions as far as they will go <b>three</b> times.  <b>NOTICE!</b> Scratch-sensitive surface. Risk of corrosion. Do not use metal brushes, metal objects, or abrasive cleaners.
Rinse	4. Then thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least <b>five</b> times with the cleaning solution. Use a disposable 50 ml syringe.

### Complete stage II: Intermediate rinsing

*Proceed as follows:*

2 rinses of at least 1 min each	1. Completely immerse the product in tap water (20-30 °C / 68-86 °F). Fully rinse all accessible surfaces <b>two</b> times for at least <b>one</b> minute each.
Move mobile components	During the rinsing, maneuver mobile components in all directions as far as they will go <b>three</b> times.  Use fresh water for each rinse.
Rinse surfaces	2. Thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least <b>three</b> times with tap water. Use a disposable 50 ml syringe.  Use fresh water for each rinse.  Allow residual water to drip off for a sufficient length of time.

### Complete stage III: Disinfection

*Proceed as follows:*

Immerse for 12 min	1. Completely immerse the product in the disinfectant solution for at least <b>twelve</b> minutes.  Make sure that all accessible surfaces remain immersed in the disinfectant solution throughout the entire disinfecting time.
Remove air bubbles	Remove all adherent air bubbles from the component surfaces.

Wipe or brush product	<p>2. Then wipe or brush the product, while it is completely immersed in the cleaning solution, using a soft, clean, lint-free wipe or with a soft bristled brush.</p> <p>Wipe or brush for at least <b>one</b> minute.</p>
Move mobile components	<p>Make sure that the disinfectant solution reaches all accessible surfaces by maneuvering the mobile components in all directions as far as they will go <b>three</b> times.</p> <p><b>NOTICE!</b> Scratch-sensitive surface. Risk of corrosion. Do not use metal brushes, metal objects, or abrasive cleaners.</p> <p>3. Thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least <b>five</b> times. Use a disposable 50 ml syringe.</p>
<p><b>Complete stage IV: Final rinsing</b></p> <p><i>Proceed as follows:</i></p>	
3 rinses of at least 1 min each	<p>1. Completely immerse the product in DI water (20-30 °C / 68-86 °F) and thoroughly rinse all accessible surfaces in <b>three</b> rinses of at least <b>one</b> minute each.</p>
Move mobile components	<p>During the rinsing, maneuver mobile components in all directions as far as they will go <b>three</b> times.</p>
Rinse surfaces	<p>2. Thoroughly rinse all surfaces of products with hidden crevices or complex geometries at least <b>three</b> times with DI water. Use a disposable 50 ml syringe.</p> <p>Use fresh water for each rinse.</p> <p>Allow residual water to drip off for a sufficient length of time.</p>
<p><b>Complete stage V: Drying</b></p> <p><i>Proceed as follows:</i></p>	
	<p>1. Dry the product using a soft, lint-free wipe.</p> <p>Use medical-quality filtered compressed air (pmax = 0.5 bar) to dry the areas you cannot reach with the wipe.</p>
	<p>2. Make sure the product is clean and undamaged. Repeat the cleaning and disinfection process if necessary.</p>
Visual inspection	<p>If any damage is detected, set aside the product immediately.</p> <p><input checked="" type="checkbox"/> This concludes the manual cleaning and disinfection process.</p>

**WARNING!** The product is not sterile after manual cleaning and disinfection.

## 4.6 Automated Cleaning and Thermal Disinfection

### 4.6.1 Manual Pre-Cleaning

#### Performing manual pre-cleaning

*Manual pre-cleaning must be performed prior to automated cleaning.*

*Proceed as follows:*

**NOTICE!** Damage resulting from rough handling. Handle the product with care. Do not knock, crush or bend. Set down carefully.

**NOTICE!** Scratch-sensitive surface. Risk of corrosion. Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads soaked with alcohol or a neutral cleaner only.

Note: The endocoupler is permanently attached to the camera head and cannot be detached.

1. Pre-clean the camera head including cable and plug at the conclusion of the surgery by wiping it with a lint-free and wetted, but not dripping wipe.  
Before wetting the wipe with the enzymatic cleaning solution (Cidezyme®/Enzol®), it should be prepared per the manufacturer's instructions.  
Wipe the product until it is free of visible residue, but for at least **one** minute.
  - Immerse for 10-30 min 2. Completely immerse the fully disassembled product for **ten** to **thirty** minutes in an alkaline cleaning solution (neodisher® MediClean forte) prepared as per manufacturer's instructions.  
Allow air bubbles to escape Make sure that the cleaning solution touches all accessible surfaces and that the product remains immersed throughout the entire cleaning time. When placing hollow components into the bath, ensure that all air escapes completely from the hidden crevices.
  - Wipe or brush external surfaces 3. While the product is soaking in the cleaning solution, clean all external surfaces using a soft, clean, lint-free wipe or with a soft bristled brush.  
Wipe or brush for at least **one** minute or until no more residue can be seen.  
Maneuver mobile components of the completely immersed in the cleaning solution product in all directions as far as they will go at least **three** times.
  - 2 rinses of at least 1 min each 4. Completely immerse all parts in tap water (20-30 °C / 68-86 °F). Fully rinse all accessible surfaces **two** times for at least **one** minute each.  
During the rinsing, maneuver mobile components in all directions as far as they will go **three** times.  
Use fresh water for each rinse.  
Allow residual water to drip off for a sufficient length of time.
- This concludes the manual pre-cleaning process.

**4.6.2 Overview**

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Cleaning solution
I	Pre-rinsing	Cold	2	Tap water	---
II	Cleaning	55/131	10	---	Alkaline cleaner: neodisher® MediClean forte, prepared as per manufacturer's instructions
III	Intermediate rinsing 1	per the equipment manufacturer's standard cycle	1	Tap water	---
IV	Intermediate rinsing 2	per the equipment manufacturer's standard cycle	1	Tap water	---
V	Thermal disinfection	90/194	5	Deionized water (A <sub>0</sub> value: > 3000)	---
VI	Drying	High (98.8/210)	30	---	---

Table 4-2: Overview of automated cleaning and thermal disinfection.

**Note:** Color anodized parts or plastic components may bleach out during automated cleaning.

**4.6.3 Performance****Performing automated cleaning and thermal disinfection**

*Use a washer/disinfector that has been validated as effective. The device must meet the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).*

*Dry removal of the products is recommended prior to automated cleaning.*

*In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.*

*Use only legally marketed reprocessing trays for automated cleaning and thermal disinfection. Observe the instructions for use of the manufacturer.*

*Proceed as follows:*

- |   |   |
|---|---|
| Follow the manufacturer's instructions! | <ol style="list-style-type: none"> <li>1. Place the reprocessing tray* into the washer/disinfector.</li> </ol> <p>Follow the manufacturer's instructions and instructions for use for the device. When loading, take care to avoid creating rinsing blind spots.</p>  |
| Visual inspection                       | <ol style="list-style-type: none"> <li>2. Start the cleaning/disinfection cycle in accordance with the manufacturer's instructions and instructions for use for the device.</li> <li>3. After drying, visually inspect the products in a well-lit area and make sure they are dry, undamaged and free of visible residue.</li> </ol> <p>Use a magnifying glass. Repeat the cleaning and disinfection process if necessary.</p> <p>Immediately set aside any damaged products.</p> |
-  **CAUTION!** Risk of scalding when unloading the device. Be sure to wear gloves.
4. Remove the products from the device.
- This concludes the automated cleaning and thermal disinfection process.
- \* for automated cleaning and thermal disinfection validation the camera head was placed in a standard reprocessing tray (OM-1002-SY).

## 4.7 Sterilization

### 4.7.1 STERRAD® Sterilization Systems 100S, NX, 100NX

#### Performing sterilization

**NOTICE!** STERRAD® sterilization may result in cosmetic device changes that do not necessarily affect the device function.

**NOTICE!** The product must be completely dry before it can be sterilized in a STERRAD® sterilizer. Loads containing moisture may cause a cycle cancellation.

**NOTICE!** Damage resulting from rough handling. Handle the product with care. Do not knock, crush or bend. Set down carefully.

*The STERRAD® Sterilization Systems, manufactured by Advanced Sterilization Products (ASP), use low-temperature, hydrogen peroxide gas plasma technology for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.*

*Refer to the STERRAD® Sterilization Systems User Guide for detailed instructions for use of any STERRAD® unit, ASP's STERRAD® Sterility Guide (SSG) at [www.sterradsterilityguide.com](http://www.sterradsterilityguide.com) or contact ASP customer service.*

*Use only reprocessing trays that are approved for STERRAD® sterilization. Observe the instructions for use of the manufacturer.*

*Proceed as follows:*

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place a STERRAD® indicator strip in the reprocessing tray.

Double-wrap the tray\* with two single layers of legally marketed sterilization wrap.

3. STERRAD® 100S und STERRAD® 100NX: Load the wrapped tray on the top shelf of the sterilization chamber. The bottom shelf must be empty.

STERRAD® NX: Load the wrapped tray on the bottom shelf of the sterilization chamber. The top shelf must be removed.

Position the tray so that the sterilization medium can act thoroughly on all areas. Do not allow any item to touch the wall of the sterilizer.

Follow the  
manufacturer's  
instructions!

4. Start the sterilization cycle in accordance with the manufacturer's instructions and instructions for use for the sterilizer.

STERRAD® sterilization was validated for the following cycles:

- STERRAD® 100S Short Cycle (top shelf only)
- STERRAD® NX Standard Cycle (bottom shelf only)
- STERRAD® 100NX Standard Cycle (top shelf only)

5. Remove the sterilized products from the sterilizer.

Ensure that the products remain sterile after reprocessing.

- This concludes sterilization.

\* for sterilization validation the camera head was placed in a standard reprocessing tray (OM-1002-SY) that was double-wrapped in polypropylene wrap (Kimguard® KC400 Sterilization Wrap, Kimberly Clark, PC 68248).

#### 4.7.2 STERIS Sterilization Systems V-PRO®

##### Performing sterilization

**NOTICE!** STERIS sterilization may result in cosmetic device changes that do not necessarily affect the device function.

**NOTICE!** The product must be completely dry before it can be sterilized in a STERIS sterilizer. Loads containing moisture may cause a cycle cancellation.

**NOTICE!** Only use trays that have been approved by STERIS for the hydrogen peroxide gas sterilization procedure.

*Validated STERIS Low Temperature Sterilization Systems: V-PRO®, V-PRO® 1, V-PRO® 1 Plus, V-PRO® maX, V-PRO® maX 2, V-PRO® 60, V-PRO® s2.*

*The V-PRO® Sterilization Systems, manufactured by STERIS, use low-temperature, hydrogen peroxide gas plasma technology for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.*

*Refer to the instructions for use of your respective STERIS unit or contact STERIS customer service.*

*Use only reprocessing trays that are approved for STERIS Low Temperature Sterilization. Observe the instructions for use of the manufacturer.*

*Proceed as follows:*

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place an indicator strip on the reprocessing tray.

Double-wrap the tray\* with two single layers of legally marketed sterilization wrap.

Follow the  
manufacturer's  
instructions!

3. Load the wrapped tray in the sterilizer as per a validated loading plan and observing the instructions of the sterilizer manufacturer.
4. Start the STERIS Low Temperature Sterilization in accordance with the manufacturer's instructions.
5. Remove the sterilized products from the sterilizer.  
Ensure that the products remain sterile after reprocessing.

This concludes sterilization.

\* for sterilization validation the camera head was placed in a reprocessing tray (OM-1002-SY) that was double-wrapped in polypropylene wrap (HALYARD\* H600 ONE-STEP\*).

### 4.7.3 STERIS Sterilization SYSTEM 1® / SYSTEM 1E®

#### Performing sterilization

**NOTICE!** STERIS sterilization may result in cosmetic device changes that do not necessarily affect the device function.

*Validated STERIS Sterilization Systems: SYSTEM 1E® Liquid Chemical Sterilant Processing System, SYSTEM 1® endo Liquid Chemical Sterilant Processing System, SYSTEM 1® Sterilization System, SYSTEM 1® PLUS Sterilization System.*

*The 1®/1E® Sterilization Systems manufactured by STERIS use peracetic acid for sterilization (S40™ Sterilant Concentrate).*

*Refer to the instructions for use of your respective STERIS unit or contact STERIS customer service.*

*Use only containers that are approved for STERIS System 1/1E sterilization. Observe the instructions for use of the manufacturer.*

*Proceed as follows:*

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place an indicator strip on the reprocessing container\*. Follow the manufacturer's instructions.
3. Load the container in the sterilizer as per a validated loading plan and observing the instructions of the sterilizer manufacturer.
4. Start the STERIS System 1/1E sterilization in accordance with the manufacturer's instructions.
5. Remove the sterilized products from the sterilizer.

Follow the manufacturer's instructions!

Ensure that the products remain sterile after reprocessing.

**NOTICE!** The product must be used directly after sterilization, as there is no way of storing it without losing its sterility.

This concludes sterilization.

\* for sterilization validation the camera head was placed in a dedicated reprocessing container (STERIS part numbers: C1220E, C1140E, C1160E, C1220S1, C1140S1, C1160S1, C1220INT, C1140INT or 1160INT).

### 4.7.4 Ethylene Oxide Gas Sterilization

#### Performing gas sterilization

*The device is validated for Ethylene Oxide Gas Sterilization and can be sterilized and aerated using the parameters given here. Observe the instructions for use of the sterilizer and the parameters appropriate to the institutional, local and national requirements.*

*Use only reprocessing trays that are approved for Ethylene Oxide Gas Sterilization. Observe the instructions for use of the manufacturer.*

*Proceed as follows:*

Follow the  
manufacturer's  
instructions!

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Double-wrap the tray\* with two single layers of legally marketed sterilization wrap.
3. Load the wrapped tray in the sterilizer as per a validated loading plan and observing the instructions of the sterilizer manufacturer.
4. Start the Ethylene Oxide Gas Sterilization in accordance with the device manufacturer's instructions, using the following validated parameters:

**Preconditioning**

Temperature	55 °C (131 °F)
Relative air humidity	70 %
Vacuum set points	0.09 bar (9 kPa)
Conditioning time	30 min

**Sterilization**

Temperature	55 °C (131 °F)
Relative air humidity	70 %
Concentration (100 % EO)	735 mg/l
Conditioning time	3 hours

**Aeration**

Aeration time	12 hours
Temperature	50 - 57 °C (122 - 135 °F)

**WARNING!** Toxic ethylene oxide residue. Risk to persons. Adhere to the full aeration time.

Note: It is the responsibility of the operator to establish that the sterilizer used meets the parameters above.

Note: Ethylene Oxide Gas Sterilization may result in cosmetic device changes that do not necessarily affect the device function.

5. Remove the sterilized products from the sterilizer.  
Ensure that the products remain sterile after reprocessing.

This concludes sterilization.

\* for sterilization validation the camera head was placed in a reprocessing tray (OM-1002-SY) that was double-wrapped in polypropylene wrap (HALYARD\* H300 ONE-STEP\*).

## 5 Product Data

### 5.1 Technical Data

Dimensions incl. endocoupler (W x H x D)	45 x 46.5 x 149 mm
Outer diameter of endocoupler	52.1 mm
Weight incl. endocoupler	330 g without cable
Image sensor	1/3" CMOS, progressive scan
Resolution	1920 x 1080 pixels
Focal length	14.25 - 28 mm
Camera cable length	3.5 m
Degree of protection	IP X7
Applied part as per IEC 60601-1	Defibrillation-proof type CF

### 5.2 Ambient Conditions

#### Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	50 kPa to 106 kPa

#### Operating conditions

Temperature	+10 °C to +35 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

#### CAUTION

##### Non-observance of the ambient conditions

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

> Observe operating conditions as well as transport and storage conditions

## 6 Disposal

### **WARNING**

#### **Contaminated device**

Risk of infection

> The device must be reprocessed prior to disposal (*chapter 4*)

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

The product is packed in a polypropylene (PP) plastic case (0.5 kg) with a polyurethane (PU) soft foam inlay (0.12 kg), polyethylene (PE) bags (0.02 kg) as well as a corrugated cardboard box (0.54 kg).

Any product carrying this symbol must be disposed of separately through dedicated electrical and electronic devices recycling. Within the EU, such disposal is taken care of by the manufacturer free of charge.

SCHÖLLY FIBEROPTIC GMBH  
Robert-Bosch-Str. 1-3  
79211 Denzlingen  
Germany