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In order to assure a safe application, standards/guidelines released by official organizations such as scientific societies of endoscopic surgery have to be considered.

2.1.3 Contraindications

Use is contraindicated if endoscopic procedures are contraindicated for any reason.

For any interventions the responsible physician must decide whether the foreseen application of the device is admissible according to professional standards based on the general condition of the patient and an individual benefit-risk analysis.

2.2 Scope of Delivery

The scope of delivery for the device includes:

- 1x light guide adapter, light-source-side
- or** 1x light guide adapter, endoscope side
- 1x instructions for use

2.3 Marking

Observe the symbols used on product and packaging.

Find a list with all symbols used in product, packaging and documents online: https://ifu.schoelly.de

2.4 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact our Service Department. You will find contact details at the end this document.

3 Operation

3.1 Staff Qualifications

This device may be operated only by medical specialists with the requisite training, advanced training, knowledge, and practical experience in the endoscopic discipline concerned as defined by the provisions in place at the site of operation.

3.2 Inspection Prior to Use

Performing a visual inspection

- Make sure that the product and all the components to be used are externally undamaged.
- Check that all the surfaces are clean and smooth.
- Do not use the product if it has sharp corners or edges, bulges, or rough surfaces.

3.3 Connecting the Light Guide Adapter

Take care when combining light guide and endoscope that the fiber optic cross-section matches.

► Screw in the matching adapter for the endoscope used on the distal end of the light guide.

► Screw in the matching adapter for the light source used onto the proximal end of the light guide.

► Connect the light guide to the light source and ensure that the light guide is firmly connected.

► Connect the endoscope to the light guide and ensure that the light guide is firmly connected.

NOTE: If the existing connector at the light guide is compatible with the used light sources or endoscopes an adapter is not required.

3.4 Stopping Usage and Pre-cleaning the Product

Complete the pre-cleaning process **immediately after use**.

- Detach the adapters used if applicable on both sides.
- Remove visible surgical residue as completely as possible using a moist, lint-free wipe.
- Dry the product using a soft, lint-free wipe.
- Place the product into a container for dry removal and seal.
- Arrange for reprocessing, ensuring that the product is reprocessed within 6 hours.

4 Reprocessing

4.1 Information on Reprocessing

Light Guide Adapters can be reprocessed with the same methods as the corresponding Light Guides.

Find an extensive description of the reprocessing in the instructions for use of the Light Guide (Document no. TPA037 and others).

Observe especially the warning and safety notices listed there.

4.2 Preparing Automated Cleaning

In addition to the processes described, light guide adapters must be prepared for automated cleaning as follows.

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Manual pre-cleaning must be performed prior to automated cleaning.

- Completely immerse the product for **10 to 30 minutes** in a cleaning solution (neodisher® MediClean forte) prepared as per manufacturer’s instructions.
- Make sure that the cleaning solution touches all accessible surfaces and that the product remains immersed throughout the entire cleaning time. Place products with cavities in the solution in a way that allows all air to escape from the cavities.
- After the appropriate exposure time: While the product is immersed in the cleaning solution, clean it using a soft, clean, lint-free wipe or with a soft bristled brush. Wipe or brush for at least **1 minute** and until no more residue can be seen.
- Move moving parts by at least **360°** or **3 times** in all directions to the end stop. Doing so the product must be completely immersed in the cleaning solution.
- Thoroughly rinse all surfaces of the product with hidden crevices or complex geometries at least **5 times** with the cleaning solution. Use a disposable 50 ml syringe.
- Then perform the automated cleaning as described in the instructions for use of the corresponding light guide.

5 Disposal

--	--

Observe the national regulations governing waste disposal and recycling when recycling and disposing of the product and its components.

--	--

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EN Instructions For Use - Light Guide Adapter

0 Scope of Validity, Identification

This document is valid for the following light-source-side and endo-scope-side light guide adapters:

Ref.	Designation
05.0100z	Light guide adapter, light-source-side, Storz systems
05.0102b	Light guide adapter, light-source-side, Wolf systems
05.0106c	Light guide adapter, light-source-side, ACMI systems (Circon)
05.0101o	Light guide adapter, light-source-side, Olympus systems
05.0108z	Light guide adapter, endoscope-side, Storz systems

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https://ifu.schoelly.de

Ref.	Designation
05.0110b	Light guide adapter, endoscope-side, Wolf systems
05.0112c	Light guide adapter, endoscope-side, ACMI systems

Also observe the respective instructions for use of light guide, endoscope and light source.

1 General Safety Information

1.1 Basic Safety Notices

This product meets the state of the art and was designed and produced according to highest quality standards. Nonetheless in all life phases dangers may occur.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

Operate the device according to the intended use only. Doing so, observe the instructions for use.

Before each use, make sure the product and the accessories used are undamaged and functional.

Keep the original packaging. Transport and store the product in its original packaging and use it to return goods if service support is required.

Observe the instructions for use of all devices and instruments used during the application of the product.

Inform us immediately in case you should notice any faults or malfunctions.

1.2 Hazard Sources

Nature and origin of the risk	Explanation	Preventive action
Risk of injury due to unauthorized modifications of the product	Unauthorized modifications influence product safety	Do not make any unauthorized modifications to the product.
Risk to the patient due to component failure	Components may fail during a surgical procedure	Keep spares on hand for emergency replacement.
Risk to the patient resulting from damaged device	The product may be damaged by improper handling	Handle the product with care. Do not continue to use the product after heavy mechanical stress or falls. Send product in for inspection after any heavy mechanical stress or fall.
Risk to the patient due to improper use	Risk to the patient,user and other persons, or premature wear of the device	Only use product according to intended use. Always follow instructions for use.
Risk of suffocation by packaging material	Children can suffocate on packaging materials	Keep packaging material out of reach of children.
Risk of injury by electromagnetic interactions	During magnetic resonance imaging (MRI) electromagnetic interactions occur. This may cause heating of metal components.	Do not operate the device in the vicinity of MRI scanners.
Risk of infection by unsterile shipping	The product is delivered in an unsterile condition	Carry out reprocessing of the product and its accessories prior to first use.
Risk of burns from hot surfaces	The product may heat up during operation	Do not touch hot surfaces. Do not bring hot surfaces into contact with patients.

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