

HD Optics



- High resolution HD optics specially designed for the FlexiVision 100 with camera head
- Rod lens system for excellent images
- Different viewing angles and wide angles, as well as different working lengths

Item no.	Working Ø	Working length	Direction of view	Field of view
ME.27120.0085	2.7 mm	110 mm	0°	95°
ME.27120.3085	2.7 mm	110 mm	30°	85°
ME.27120.7085	2.7 mm	110 mm	70°	80°
ME.27210.0085	2.7 mm	187 mm	0°	95°
ME.27210.3085	2.7 mm	187 mm	30°	85°
ME.27210.7085	2.7 mm	187 mm	70°	80°
ME.40175.00100	4.0 mm	175 mm	0°	100°
ME.40175.30100	4.0 mm	175 mm	30°	100°
ME.40175.70100	4.0 mm	175 mm	70°	100°
ME.500312.0070	5.0 mm	312 mm	0°	70°
ME.500312.3070	5.0 mm	312 mm	30°	70°
ME.100344.0070	10.0 mm	344 mm	0°	70°
ME.100344.3070	10.0 mm	344 mm	30°	70°